

Product datasheet for **SC304697**

BARHL1 (NM_020064) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	BARHL1 (NM_020064) Human Untagged Clone
Tag:	Tag Free
Symbol:	BARHL1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC304697 representing NM_020064. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCC**CGATCGCC**
ATGGAAGGCTCCAATGGCTTTGGGATCGACTCCATTCTCTCCACCGCGCGGGCAGCCCCGCCCTTCCC
AAGGGGACCCCTTGCTCGGGGACTGCCGTTGCCCTGGAGCTGAGTCCACGCTCAGAGAGCAGCAGC
GACTGCTCTTCGCCAGCCTCTCCAGGAAGGACTGTTTGAGACGGGGACCCACGGCCTGGCGGGCA
TCCGGCCAGGTTTGACTCCCACCTGCAGCCCGGGCAGCTCTCAGCCCCGGCCAGTCGCGCACCGTC
ACCTCCTCCTTTCTGATCAGGGACATCCTTGCCGACTGCAAACCACTCGCGGCTGTGCACCCTACTCT
AGCAGCGGGCAGCCGGCAGCCCCTGAGCCTGGGGCCGCTTGCGGCCAAGGCCGCGAGGACTTTAGA
GACAAGCTGGACAAAAGTGGCAGCAACGCCTCATCGGACTCTGAGTATAAAGTGAAGGAGGAGGGGAC
CGCGAGATCTCCAGCTCCAGGGACAGTCCCCCGGTGCGCCTGAAAAAGCCACGCAAGGCGCGCACGGCC
TTCACCGACCATCAGCTGGCGCAGCTGGAGCGCAGCTTCGAGCGGCAGAAGTACCTGAGCGTGCAGGAC
CGCATGGAGCTCGCCGCTCGCTCAACCTACCGACACGCAGGTCAAGACCTGGTACCAGAACCAGGAG
ACTAAATGGAAGCGACAGACGGCCGTGGGTGGAGCTGCTGGCGGAGGCAGGCAATTACTCAGCGCTC
CAGCGGATGTTCCCGTCGCTTATTTCTACCCGACAGTCTGGTTTCAACCTGGACCCCGCGCGGCG
CTCTACCTGTACCGGGCCCCAGCGCGCCGCGCTGCTCTCAGAGACCTCTGGTGCCCGCATCTCTC
ATCCACGACTCCAGGGCGCCAGCGAGCGCCCCCGCGCTGCCCCCTGGCCGCGCTCTCCACGC
GCCGCGCAGCCTCGG**TGA**
ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Restriction Sites:	SgfI-MluI
ACCN:	NM_020064
Insert Size:	984 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_020064.3</u>
RefSeq Size:	1907 bp
RefSeq ORF:	984 bp
Locus ID:	56751
UniProt ID:	<u>Q9BZE3</u>
Cytogenetics:	9q34.13
Protein Families:	Transcription Factors
MW:	35.1 kDa