

Product datasheet for **SC304710**

PAPOLB (NM_020144) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PAPOLB (NM_020144) Human Untagged Clone
Tag:	Tag Free
Symbol:	PAPOLB
Synonyms:	PAPT; TPAP
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC304710 representing NM_020144.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGATGCCGTTTCCGGTGACAACCCAGGGACCACCGCAGCCGGCCGCCCGCAATCGCTACGGCGTC
TCCTCGCCTATCAGTCTAGCGGTCCCCAAGGAGACGGACTGCCTCCTCACCCAGAGGCTAATAGAAACC
CTCAGGCCCTTCGGGTCTTTCGAAGAGGAAGAGGAACTGCAGCGCAGGATTTTAGTTTTGGAAAAATTA
AATAATCTGGTAAAGGAATGGATACGCGAAATCAGTGAAAGCAAGAGTCTTCCCCAGTCTGTAATTGAA
AACGTTGGAGGAAAGATTTTACGTTTGGCTCTTACAGATTAGGAGTACATACGAAAGGCGCAGATATT
GACGCCTTGTGCGTTGCACCAAGTCATGTGGATCGAAGCGACTTTTTACCTCATTCTATGCTAAACTG
AAACTACAGGAGGAAAGTAAAGATTTAAGGGCTGTCGAGGAGGCATTTGTCCAGTTATCAAAGTGTGT
TTTGATGGGATAGAGATTGATATTTTATTGCAAGATTAGCACTACAGACTATCCAGAAGATTTGGAC
TTAAGAGATGACAGTCTACTTAAAAATTTAGACATTAGATGCATAAGAAGCCTTAATGGTTGCCGGGTA
ACCGATGAAATTTTACATCTAGTGCCAAACATTGACAACCTCAGGCTGACTCTGAGAGCCATCAAAGTGT
TGGGCCAAGTGCCACAATATCTATTCCAATATATTAGTTTTCTCGGAGGTGTTTCTGGGCCATGCTA
GTAGCAAGAATTGTGACGCTTTATCCAATGCAGTAGCGTCAACTCTTGTACGGAATTTCTTGGTA
TTTTCAGAATGGGAATGGCCAAACCCAGTGTACTGAAGGAGCCTGAAGAACGGAATCTTAATTTGCCT
GTATGGGACCCAAGAGTAAATCCAGTGATAGGTACCATCTTATGCCTATCATCACACCAGCATACCCA
CAGCAGAACTCCACATACAACGTGTCTATTTCAACCAGGATGGTCATGATTGAGGAGTTTAAACAGGGG
CTTGCTATCACACAGGATTTTGCTAAGTAAGGACAGTGGTCCAAACTCTTTGAAGCTCCAAGCTTC
TTTCAAAGTACAAGCATTATTTGTAATCTGCAAGTGCATCAACAGAAAAACAACATTTAGAATGG
GTGGGCTTGGTGAATCAAAGATCCGAATCCTGGTTGGGAGCTTGGAGAAGAATGAATTTATTACTG
GCACATGTGAATCCACAGTCATTTCCAGCACCCAAAGAAAATCCTGATATGGAAGAATTTGTACAATG
TGGGTGATTGGTTAGGGCTAAAAAGCCAGATAATTCTGAAATCTCAGCATTGATCTCACCTATGAT
ATCCAGTCTTTCACAGATACTGTTTATAGGCAAGCAGTGAATAGTAAGATGTTTGGATGGGTATGAAA
ATTACTGCAATGCATTTAAGAAGAAAGAACTTACCAGCTGCTGCCTCATCATGTGCTTCAGGACAAG
AAAGCACACTCAACAGAAGGTAGAAGATTGACAGATTTGAACGACAGCAGCTTTGACTTGTCTGCAGGC
TGTGAAAACAGCATGTCTGTGCCTTCATCTACTAGCACTATGAAGACAGGCCATTGATTAGCAGTTCT
CAGGGTAGAAAACAGTCTGCCTTGGCTGTAATGACAGCATCTGTGGCTAACATACAGGCTACTGAATTT
TCCTTGCAACAGGTGAATACCAATGAAAGTTTCAAGGGTTGCATTAACGAAAAGTATTCCTCAGCTGTC
TCTCAGCCTGCCATTTCTCCATACCAAAGGCCATGGTCCAGAGTTGTTTCTCAACATGTCTCATA
AGCCATCCAGACCTTCAGGAAACTCAGCAACAAACATATCTAATCCTATAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_020144

Insert Size: 1914 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_020144.4](#)

RefSeq Size: 4271 bp

RefSeq ORF: 1914 bp

Locus ID: 56903

UniProt ID: [Q9NRJ5](#)

Cytogenetics: 7p22.1
Protein Pathways: RNA degradation
MW: 71.8 kDa