

## Product datasheet for **SC304721**

### Cezanne (OTUD7B) (NM\_020205) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Cezanne (OTUD7B) (NM_020205) Human Untagged Clone
Tag:	Tag Free
Symbol:	Cezanne
Synonyms:	CEZANNE; ZA20D1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC304721 representing NM\_020205.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGACCCTGGACATGGATGCTGTTCTGTCCAGATTTTGTCCGTCCACAGGAGCAGAGCCAGGGCTAGCC
CGAGATCTCCTAGAAGGAAAGAATTGGGATGTGAATGCCGCCCTCAGTGATTTTGAACAGCTACGTCAA
GTCCATGCTGAAACCTACCCCATCCTTTAGTGAGGGGAGTGGTGGCTCCAGGACCCCTGAAAAAGGG
TTTTCTGACAGAGAGCCTACTCGCCCTCCCGACCCATCCTCCAGCGGCAGGATGACATCGTTCAAGAA
AAACGCCTGTCTAGGGGCATCTCCACGCCAGCTCCAGCATTGTTTCCCTGGCCCGGTCCCATGTCTCC
TCCAATGGTGGGGTGGGGGAGCAATGAGCACCCCTGGAAATGCCATCTGTGCCTTCCAGCTTCCA
GATCTCACTGTATACAATGAAGACTTCCGCAGCTTCATAGAGAGAGACCTCATTGAGCAGTCCATGCTG
GTTGCCTTGGAACAGGCAGGGCCTTTGAACTGGTGGTGAGTGTGGATCCACCTCTCAGAGGCTGCTT
CCTTTGGCAACTACTGGAGATGGAACTGCCTCCTGCATGCAGCCTCCCTTGGAAATGTGGGTTTCCAT
GATCGGGACTTGATGCTCGGAAAGCTTTGTATGCACTGATGGAGAAGGGAGTTGAGAAGGAAGCGTTG
AAAAGGCCTGGAGGTGGCAGCAGACACAGCAATAAAGAGTCAGGGCTGGTATACACAGAAGATGAA
TGGCAGAAGGAGTGGAATGAAGTATCAAGCTTGCCTCAAGTGAACCCCGAATGCATCTAGTACCAAT
GGAGCCAACTGTGGTGGGGTGGAGAGTTCTGAGGAGCCTGTATATGAGAGCCTTGAAGAGTTTACAGTC
TTTGTCTTGCTCATGTGCTTAGGAGGCCATAGTCGTCGTGGCAGACACCATGCTGAGGGACTCCGGA
GGGAAGCATTGGCCCTATCCCTTTGGAGGAATCTATCTGCCTTTGGAGTCCCAGCCAGCCAGTGT
CACCGCTCCCCTCTGGTCTCGCTATGATCAGGCCACTTTTCTGCACTCGTGTCCATGGAGCAGAAG
GAGAATACCAAGGAACAAGCTGTGATCCCACTTACAGATTACAGATATAAGCTGCTGCCCTGCCTTT
GCTGTGGACCTGGAAAGGGCTGGGAGTGGGCAAGATGATAGTACAATGTCCGATTGGCCAGTTA
ATTCTGTCCCTAGAGGTCAAATTGCATCTGCTGCATAGCTACATGAATGTGAAGTGGATCCCACTGTCC
TCTGATGCACAGGCTCCTTGGCCAGCCTGAGTCCCCACCCTCAGCTGGAGATGAGCCCCGGTCC
ACTCCTGAGTCTGGAGACTCAGACAAGGAGTCAAGTGGCAGCAGTCCACCAGCAACGAGGGCGGCCGG
CGGAAGGAGAAGTCAAAGCGAGATCGGGAGAAGGACAAGAAGAGAGCAGATTCTGTGGTAAACAACTG
GGCAGCTTTGGCAAAACCTTGGCAGCAAGCTCAAGAAGAATGGGGGGCTGATGCACAGCAAGGGT
TCAAAGCCTGGAGGGTGGGACAGGGTGGGAGGAAGCAGCGGCACTGAGACACTGGAGAAGAAGAAG
AAAACTCACTGAAGAGCTGGAAGGGTGGCAAGGAGGAGGAGCAGTGGGGTGGGCTGTGTCTGAGAAG
CCCCAGCTGAGTCTGTTGGTAACGGAGGGAGCAAGTATAGCCAGGAGGTGATGCAGAGCCTGAGCATT
CTGAGGACTGCCATGCAAGGGGAGGGAAAGTTATTTTGTGGAAACCCTGAAGATGGGTACCCGTAC
CAGTATCAGGAGGAAATGATCCAGCGCTACCTTTCTGATGCTGAGGAGAGATTCTTGGCAGAACAGAAG
CAGAAGGAGGCAGAGAGGAAGATCATGAATGGAGGAATAGGGGGTGGCCCTCCTCCAGCCAAAAGCCA
GAGCCAGATGCTAGGGAAGAGCAGCCGACCCGTCCCCAGCAGAGTCCAGGGCAATGGCATTTCCTACT
GGCTACCCTGGGACTTTACTATCCCTCGGCCGTCTGGGGCGGAGTCCACTGCCAGGAACCCCGGAGG
CAGTTGGCAGGGGTCCATGTGTGGGGGCTACCACCATATGCCACCTCCCCAGACAGTGCCTCCT
GGGCGACCCTACCCACCAGGACAGCATCCCTTCTCTGGAGCCAGGCAGCCACTTAAGGATGGACTT
CACAGGGTGCCTTGTACCACCCCTACCAGTGGTGTATCCTATAGCAATGGCTACAGAGAGCCC
CCTGAGCCAGATGGATGGGCTGGAGGTCTCCGGGGCCTTCCCCAACTCAGACCAATGCAAAACCG
AACTGCAGCTTCTATGGACACCCTGAGACAAACAATTCTGTTCTGTTGTACAGGGAAGAAGTGAAG
AGGAGGGAGCGGGAACCGGATGGGGAGCTCCTGGTGCACAGTTCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_020205  
**Insert Size:** 2532 bp

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_020205.3</a>
<b>RefSeq Size:</b>	6415 bp
<b>RefSeq ORF:</b>	2532 bp
<b>Locus ID:</b>	56957
<b>UniProt ID:</b>	<a href="#">Q6GQQ9</a>
<b>Cytogenetics:</b>	1q21.2
<b>MW:</b>	92.5 kDa
<b>Gene Summary:</b>	Negative regulator of the non-canonical NF-kappa-B pathway that acts by mediating deubiquitination of TRAF3, an inhibitor of the NF-kappa-B pathway, thereby acting as a negative regulator of B-cell responses. In response to non-canonical NF-kappa-B stimuli, deubiquitinates 'Lys-48'-linked polyubiquitin chains of TRAF3, preventing TRAF3 proteolysis and over-activation of non-canonical NF-kappa-B. Negatively regulates mucosal immunity against infections (By similarity). Deubiquitinates ZAP70, and thereby regulates T cell receptor (TCR) signaling that leads to the activation of NF-kappa-B (PubMed:26903241). Plays a role in T cell homeostasis and is required for normal T cell responses, including production of IFNG and IL2 (By similarity). Mediates deubiquitination of EGFR (PubMed:22179831). Has deubiquitinating activity toward 'Lys-11', 'Lys-48' and 'Lys-63'-linked polyubiquitin chains (PubMed:27732584). Has a much higher catalytic rate with 'Lys-11'-linked polyubiquitin chains (in vitro); however the physiological significance of these data are unsure (PubMed:27732584). Hydrolyzes both linear and branched forms of polyubiquitin.[UniProtKB/Swiss-Prot Function]