

Product datasheet for **SC304748**

DAZ2 (NM_020363) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DAZ2 (NM_020363) Human Untagged Clone
Tag:	Tag Free
Symbol:	DAZ2
Synonyms:	pDP1678
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC304748 representing NM_020363.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTCTGCTGCAAATCCTGAGACTCCAACTCAACCATCTCCAGAGAGGCCAGCACCCAGTCTTCATCA
GCTGCAGCTAGCCAAGGCTGGGTGTTACCAGAAGGCAAAATCGTGCCAAACACTGTTTTTGTGGTGGA
ATTGATGCTAGGATGGATGAAACTGAGATTGGAAGCTGCTTTGGTAGATACGGTTCAAGTAAAGAAGTG
AAGATAATCACGAATCGAACTGGTGTGTCCAAAGGCTATGGATTGTTTCGTTTGTAAATGACGTGGAT
GTCCAGAAGATAGTAGGATCACAGATACATTTCCATGGTAAAAAGCTGAAGCTGGGCCCTGCAATCAGG
AAACAAAAGTTATGTGCTCGTCATGTGCAGCCACGTCCTTTGGTAGTTAATCCTCCTCCTCCACCACAG
TTTCAGAACGTCGGCGGAATCCAAACTGAAACCTACCTGCAGCCCCAAATCACGCCGAATCCTGTA
ACTCAGCACGTTCAAGGCTATTCTGCTTATCCACATTACCAGGTCAGGTCATCACTGGATGTCAGTTG
CTTGATATATAATTATCAGGAATATCCTACTTATCCCGATTCAAGTATTCAGGTCACCACTGGATATCAG
TTGCTGTATATAATTATCAGCCATTTCTGCTTATCCAAGATCACCATTTCAGGTCAGTCTGGATAT
CAGTTGCCTGTATATAATTATCAGGCATTTCTGCTTATCCAAATTCACCATTTCAGTCCGCACTGGA
TATCAGTTCCTGTATACAATTATCAGCCATTTCTGCTTATCCAAGTTCACCATTTCAGGTCAGTCTGCT
GGATATCAGTTGCCTGTATATAATTATCAGGCATTTCTGCTTATCCAAATTCACCATTTCAGTCCGCA
ACTGGATATCAGTTCCCTGTATACAATTATCAGGCATTTCTGCTTATCCAAATTCACCAGTTCAGGTC
ACCACTGGATATCAGTTGCCTGTATACAATTATCAGGCATTTCTGCTTATCCAAATTCACCAGTTCAG
GTCACCACTGGATATCAGTTGCCTGTATACAATTATCAGGCATTTCTGCTTATCCAAATTCACCATTTC
CAGGTCACCACTGGATATCAGTTGCCTGTATATAATTATCAGGCATTTCTGCTTATCCAAATTCACCA
TTTCAGGTCACCACTGGATATCAGTTGCCTGTATATAATTATCAGGCATTTCTGCTTATCCAAATTCAC
CCATTTTCAGGTCACCACTGGATATCAGTTGCCTGTATATAATTATCAGGCATTTCTGCTTATCCAAAT
TCACCAATTTCAAGTTCACCACTGGATATCAGTTGCCTGTATATAATTATCAGGCATTTCTGCTTATCCA
AGTTCACCAATTTCAAGTTCACCACTGGATATCAGTTGCCTGTATATAATTATCAGGCATTTCTGCTTAT
CCAAGTTCACCAATTTCAAGTTCACCACTGGATATCAGTTGCCTGTATATAATTATCAGGCATTTCTGCT
TATCCAAATTCAGCAGTTCAGGTCACCACTGGATATCAGTTCCATGTATACAATTACCAGATGCCACCG
CAGTGCCCTGTTGGGGAGCAAAGGAGAAATCTGTGGACCGAAGCATACAATGGTGGTATCTGTCTGT
TTAATCCAGAGAAGAGACTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_020363

Insert Size: 1677 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_020363.3</u>
RefSeq Size:	3853 bp
RefSeq ORF:	1677 bp
Locus ID:	57055
UniProt ID:	<u>Q13117</u>
Cytogenetics:	Yq11.223
MW:	63.1 kDa
Gene Summary:	This gene is a member of the DAZ gene family and is a candidate for the human Y-chromosomal azoospermia factor (AZF). Its expression is restricted to premeiotic germ cells, particularly in spermatogonia. It encodes an RNA-binding protein that is important for spermatogenesis. Four copies of this gene are found on chromosome Y within palindromic duplications; one pair of genes is part of the P2 palindrome and the second pair is part of the P1 palindrome. Each gene contains a 2.4 kb repeat including a 72-bp exon, called the DAZ repeat; the number of DAZ repeats is variable and there are several variations in the sequence of the DAZ repeat. Each copy of the gene also contains a 10.8 kb region that may be amplified; this region includes five exons that encode an RNA recognition motif (RRM) domain. This gene contains one copy of the 10.8 kb repeat. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008] Transcript Variant: This variant (1) represents the longest transcript. Variant 1 encodes isoform 1 which includes 15 DAZ repeats.