

Product datasheet for **SC304765**

TMEM63C (NM_020431) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TMEM63C (NM_020431) Human Untagged Clone
Tag:	Tag Free
Symbol:	TMEM63C
Synonyms:	C14orf171; CSC1; hsCSC1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC304765 representing NM_020431.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGTCTGCCTCACCAGACGACCTGAGTACAGGGGAAGGTTACAGAACATGACAGTGGATGAATGCTTC
CAGTCTCGGAACACCGTCTCCAGGGGACGCCCTTTGGGGGTGTCCACCGTGTGCTCAACATC
GCCCTGTGGGTGCTCGTCCTTGTGGTTTACTCCTTCTCCGAAAAGCTGGTGGGACTATGGGCGCCTG
GCTCTGCTGATACACAATGACAGCCTGACCTCGTGATCTATGGGGAGCAGAGCGAGAAGACATCTCCC
TCGGAGACTTCTTGGAGATGGAACGCAGAGACAAGGGATTCTGTTCTGGTTCTTCAACAGCATAACA
ATGAAGGACGAGGATCTGATTAACAAGTGTGGGACGACGCGCATCTACATCGTGTCCAGTACCAC
CTCATCATCTTTGTGCTCATCTGTATCCCCTCCCTGGGCATCATTTTGCCCATCAACTATACTGGA
TCTGTTCTGGACTGGAGCAGTCACTTTGCTCGGACCACCATTGTCAATGTCTCCAGAGAGCAAGCTC
CTGTGGCTGCATAGCCTGCTGCCTTCTTCTACTTCATCACCACCTTCATGTTTCATGGCTCATCACTGC
CTGGGGTTTGACCCAGGAATAGCCAAAAGGTCACAAGGACACTAATGATCACCTATGTGCCAAGGAC
ATTGAAGACCCAGAACTCATATTAAGCATTTTACAGAGGCCATCCAGGCAGTGTGTCGACAAGAGTC
CACTTCTGCTACGACGTAGGAACCTGATCGACTTGGACGATCAGAGGCCCATGCCATGCGGGGCCGG
CTTTTCTATACAGCCAAGGCCAAGAAGACTGGGAAGGTGATGATCAGGATCCACCCTGTGCCCGCTG
TGCTTCTGCAAGTGTGGACCTGCTTCAAGGAGGTGGATGCAGAGCAGTATTACAGCGAGCTAGAGGAG
CAGCTAACGGACGAGTTCAACGCCGAGCTCAACCGCGTCCGCTCAAGCGGCTGGACCTGATCTTTGTC
ACCTTCCAGGACTCCAGGATGGCCAAGCGTGTCCGTAAGGATTACAAGTATGTCCAGTGTGGTGTGCAA
CCCCAGCAGTCTCAGTGACCACCATCGTCAAATCATATTACTGGAGGGTCACTATGGCCCCACACCCC
AAAGACATATTTGGAACACCTGTCTGTCGGCGCTTCTTTTGGTGGGCCGCTTTATCGCAATCAAC
ACCTTCTCTTCTTCTTCTTCTTCTTCTCACCACGCTGCCATCATGAACACTATCGACATGTAC
AACGTCACCCGCCCATCGAGAAGCTGCAGAACCAATTGTGACCCAGTCTTCCCTCTGTGATGCTC
TGGGGCTTACAGTGATACTGCCTCTGATTGTCTACTTCTCCGCCTTCTCGAGGCCACTGGACCAGA
TCAAGTCAGAATCTGGTCATGGTGCACAAGTGTACATCTTTCTGGTGTTCATGGTAGTATTCTGCC
TCTATGGGACTGACCAGTTGGATGTCTTCTCCGCTGGCTTTTGACATCTACTATCTAGAGCAAGCA
TCCATCAGTTCCAGTGTGTCTTCTGCCAGACAACGGCGCTTCTTTGCAACTACGTGATCACGGCA
GCTTTACTTGGCACAGGCATGGAGCTGCTGCGTCTGGGTCACTTCTGCTACAGCACCCGCTCTTC
TTCTCTAGATCAGAGCCAGAGAGTCAACATCAGAAAGAACCAGGCCATAGACTTCCAGTTTGGGCGT
GAGTATGCGTGGATGATGAACGTGTTACGCGTGGTGTGGGTACAGCATCACTTCCCCATCATTGTG
CCTTTTGGGTTGCTTACCTGTGCATGAAGCACTTGACGGATCGTATAACATGTACTACTCCTTTGCA
CCCACAAACTGAACGAGCAGATCCACATGGCTGCCGCTCCCAGGCCATCTTGGCCACTCTTGGGT
CTGTTCTGGATGCTGTTCTTCTCCATCTGCGGTTGGGTTCTCTCCACGCCATCACCATCTTTCCCTG
TCCACCCTCTCATTGCCATGGTATTGCCTTGTGGCATTCTTGGGAAGCTTCGGATGGTTGCC
GACTACGAGCCCAGGAGGAGGAGATCCAGACAGTGTGACATGGAGCCAAGCAGCACCTCCACAG
CCCACCTCCCTCTGTATGTGGCCACCGTGTGCAAGAACCGGAGTTGAATCTGACCCCCGCTCTCC
CCAGCCAGGCACACCTATGGCACCATGAACAACCAGCCGGAAGAGGGAGAAGAAGAGAGTGGTCTGAGG
GGCTTTGCGAGGGAGCTAGACTCGGCCAGTTCAGGAAGGGCTGGAACGGAGGCCAGAACAGTAC
CACTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI
ACCN: NM_020431
Insert Size: 2421 bp

OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	<p>This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.</p>
Components:	<p>The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).</p>
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_020431.3
RefSeq Size:	5434 bp
RefSeq ORF:	2421 bp
Locus ID:	57156
UniProt ID:	Q9P1W3
Cytogenetics:	14q24.3
Protein Families:	Transmembrane
MW:	93.3 kDa
Gene Summary:	Acts as an osmosensitive calcium-permeable cation channel.[UniProtKB/Swiss-Prot Function]