

Product datasheet for **SC304795**

GALNTL1 (GALNT16) (NM_020692) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GALNTL1 (GALNT16) (NM_020692) Human Untagged Clone
Tag:	Tag Free
Symbol:	GALNT16
Synonyms:	GalNAc-T16; GALNACT16; GALNTL1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC304795 representing NM_020692.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAGGAAGATCCGCGCCAATGCCATCGCCATCCTGACCGTAGCCTGGATCCTGGGCACTTCTACTAC
TTATGGCAGGACAACCGAGCCACGCAGCATCCTCCGGCGCCGGGGCGCGCAGAGGGCAGGCAGGAGG
TCGGAGCAGCTCCGCGAGGACCGCACCATCCCGCTCATTGTGACAGGAACTCCCTCGAAAGGCTTTGAT
GAGAAGGCCTACCTGTCCGCCAAGCAGCTGAAGGCTGGAGAGGACCCCTACAGACAGCACGCCTTCAAC
CAGCTGGAGAGTGACAAGCTGAGCCAGACCGCCATCCGGACACCCGCCATTACAGCTGCCATCT
GTGTCTACTCCTCGGACCTGCCAGCCACCAGCGTCATCATCACCTTCCACAATGAGGCCGTTCCACC
CTGCTGCGCACAGTGAAGAGTGTCTGAACCGAACTCCTGCCAATTGATCCAGGAGATCATTTTAGTG
GATGACTTCAGCTCAGATCCGGAAGACTGTCTACTCCTGACCAGGATCCCCAAGGTCAAGTGCCTGCGC
AATGATCGGCGGAAGGGCTGATCCGGTCCGAGTGCCTGGGGCGGACGTGGCTGCAGCTACCGTTCCTC
ACTTTTCTGGATAGCCACTGCGAAGTGAACACCGAGTGGCTGCCGCCATGCTGCAGCGGTGAAGGAG
GACCACACCCCGCTGGTGAAGTCCCATCATTGATGTCATCAGTCTGGATAATTTTGCCTACCTTGCAGCA
TCTGCTGACCTTCGTGGAGGGTTCGACTGGAGCCTGCATTTCAAGTGGGAGCAGATCCCTCTTGAGCAG
AAGATGACCCGGACAGACCCACCAGGCCATAAAGGACGCCTGTCATAGCTGGAGGAATCTTCGTGATC
GACAAGTCTGGTTAAACCACTTGGGAAAGTATGATGCCAGATGGACATCTGGGGGGGAGAGAATTTT
GAGCTCTCCTCAGGGTGTGGATGTGTGGTGGCAGTCTGGAGATCGTCCCCTGCAGCCGGGTGGGCCAT
GTCTTCAGGAAACGGCACCCCTACAACCTCCCTGAGGGTAATGCCCTCACCTACATCAGGAATACTAAG
CGCACTGCAGAAGTGTGGATGGATGAATACAAGCAATACTACTATGAGGCCCGGCCCTCGGCCATCGGG
AAGGCCTTCGGCAGTGTGGCTACGCGGATAGAGCAGAGGAAGAAGATGAACTGCAAGTCTTCGCTGG
TACCTGGAGAACGTCTACCCAGAGCTCACGGTCCCGTGAAGGAAGCACTCCCGGCATCATTAAAGCAG
GGGGTGAAGTCTTAGAATCTCAGGGCCAGAACACAGCTGGTGAATTCCTGCTTGAATGGGATCTGC
AGAGGGTCTGCCAAGAACCCGCAGCCCGCCAGGCATGGCTGTTGAGTACCACCTCATCCAGCAGCAG
GGGAAGTGCCTGGCTGCCACCTCCACCTTAATGTCTCCCTGGATCCCCAGTCACTGCAGATGTGC
AACCTAGAGAAGGCAAGCAGAAATGGAGGAGAAAAGGATCTTTCATCCAGCATTGAGTGCAGTGGCTC
TGCCTGGAGACAAAGCCTGCCAGCTGGTACCAGCAAGTGTGAGGCTGACGCCAGGCCAGCAGTGG
CAGCTGTTGCCACACATGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_020692

Insert Size: 1677 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_020692.2](#)

RefSeq Size: 3148 bp

RefSeq ORF: 1677 bp

Locus ID: 57452

UniProt ID: [Q8N428](#)

Cytogenetics: 14q24.1

Protein Families:	Transmembrane
Protein Pathways:	Metabolic pathways, O-Glycan biosynthesis
MW:	63.1 kDa
Gene Summary:	<p>Catalyzes the initial reaction in O-linked oligosaccharide biosynthesis, the transfer of an N-acetyl-D-galactosamine residue to a serine or threonine residue on the protein receptor. [UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (2) differs in the 3' UTR compared to variant 1. Both variants 1 and 2 encode the same protein.</p>