

Product datasheet for **SC304797**

KCNS2 (NM_020697) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	KCNS2 (NM_020697) Human Untagged Clone
Tag:	Tag Free
Symbol:	KCNS2
Synonyms:	KV9.2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC304797 representing NM_020697.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGACCGGCCAGAGCCTGTGGGACGTGTCGGAGGCTAACGTGAGGACGGGGAGATCCGCATCAATGTG
GGCGGCTTCAAGAGGAGGCTGCGCTCGCACACGCTGCTGCGCTTCCCGAGACGCGCCTGGGCCGTTG
CTGCTCTGCCACTCGCGGAGGCCATTCTGGAGCTCTGCGATGACTACGACGACGTCAGCGGGAGTTC
TACTTCGACCGCAACCCCTGAGCTCTTCCCTACGTGCTGCATTTCTATCACACCGCAAGCTTCAGTTC
ATGGCTGAGCTATGTGTCTTCTCCTTCAGCCAGGAGATCGAGTACTGGGCATCAACGAGTTCCTCATT
GACTCCTGCTGCAGTACAGTACCATGGCCGCAAAGTAGAGCCCGAGCAGGAGAAGTGGGACGAGCAG
AGTGACCAGGAGACACCACGTCTTCTTCGATGAGATCCTTGCCTTCTACAACGACGCCTCAAGTTC
GATGGGACGCCCTCGGCAACTCCGACGGCAGCTGTGGTGGCCTGGACAACCCCGGCTACTCAGTG
CTGAGCAGGTCTTCAGCATCCTGTCCATCCTGGTGGTATGGGTCCATCATCACCATGTGCCTCAAT
AGCCTGCCGATTTCAAATCCCTGACAGCCAGGGCAACCCTGGCGAGGACCCTAGTTCGAAATCGTG
GAGCACTTTGGCATTGCCTGGTTCACATTTGAGCTGGTGGCCAGGTTTGCTGTGGCCCTGACTTCCTC
AAGTTCCTCAAGAATGCCCTAAACCTTATTGACCTCATGTCCATCGTCCCTTTTACACTACTCTGGTG
GTGAACCTGGTGGTGGAGAGCACACTACTTTAGCCAACCTGGGCAGGTTGGCCAGGTCCTGAGGCTG
ATGCGGATCTTCCGCATCTTAAAGCTGGCCAGGCACTCCACTGGCCTCCGCTCCCTGGGGCCACTTTG
AAATACAGCTACAAGAAGTAGGGCTGCTTGTCTACTCTCCGTGGGATTTCCATCTTCTCCGTG
GTGGCTACACCATTTGAAAGGAGGAGAACGAGGGCTGGCCACCATCCCTGCCTGCTGGTGGTGGCT
ACCGTCAGTATGACCACAGTGGGTACGGGGATGTGGTCCAGGGACCACGGCAGGAAAGCTGACTGCC
TCTGCCTGCATCTTGGCAGGCATCCTCGTGGTGGTCTGCCATCACCTTGATCTTCAATAAGTTCTCC
CACTTTTACCGGCGCAAAGCAACTTGAGAGTGCCATGCGCAGCTGTGACTTTGGAGATGGAATGAAG
GAGGTCCCTTCGGTCAATTTAAGGGACTATTATGCCCATAAAGTTAAATCCCTTATGCAAGCCTGACG
AACATGAGCAGGAGCTCACCAAGTGAAGTCAAGTAAATGATTCCCTACGTTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites:

SgfI-MluI

Plasmid Map:



ACCN: NM_020697

Insert Size: 1434 bp

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_020697.3</u>
RefSeq Size:	3219 bp
RefSeq ORF:	1434 bp
Locus ID:	3788
UniProt ID:	<u>Q9ULS6</u>
Cytogenetics:	8q22.2
Protein Families:	Druggable Genome, Ion Channels: Potassium, Transmembrane
MW:	54.2 kDa
Gene Summary:	Potassium channel subunit that does not form functional channels by itself. Can form functional heterotetrameric channels with KCNB1 and KCNB2; modulates the delayed rectifier voltage-gated potassium channel activation and deactivation rates of KCNB1 and KCNB2. [UniProtKB/Swiss-Prot Function]