

## Product datasheet for **SC304800**

### HHATL (NM\_020707) Human Untagged Clone

#### Product data:

|                           |                                        |
|---------------------------|----------------------------------------|
| Product Type:             | Expression Plasmids                    |
| Product Name:             | HHATL (NM_020707) Human Untagged Clone |
| Tag:                      | Tag Free                               |
| Symbol:                   | HHATL                                  |
| Synonyms:                 | C3orf3; GUP1; MBOAT3; MSTP002; OACT3   |
| Mammalian Cell Selection: | Neomycin                               |
| Vector:                   | pCMV6-Entry (PS100001)                 |
| E. coli Selection:        | Kanamycin (25 ug/mL)                   |



[View online »](#)

**Fully Sequenced ORF:** >SC304800 representing NM\_020707.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGGGCATCAAGACAGCATTGCCGGCGGCTGAGCTGGGCCTCTACTCTCTGGTGCTGAGTGGGGCCCTG
GCCTATGCTGGCCGGGGCCTCCTTGAGGCTTCACAAGATGGGGCCACAGGAAGGCCTCCGGGAGTCT
GTGCGACCTGGCTGGGAGTACATTGGCCGGAAGATGGATGTGGCTGACTTCGAGTGGGTGATGTGGTTC
ACCTCCTTTTCGCAACGTCATCATCTTTGCCCTCTCCGGACATGTGCTGTTTGCTAAACTCTGCACGATG
GTTGCCCCAAAGCTCCGCTCCTGGATGTATGCTGTGTACGGGGCCTTGGCTGTATGGGCACAATGGGC
CCTTGGTACCTGCTGCTGCTGCTTGGTCACTGTGTGGGCTCTATGTGGCCTCGCTTTTGGGCCAGCCC
TGGCTCTGTCTTGGCCTTGGCTTGGCCAGCCTGGCCTCCTTCAAGATGGACCCCTAATCTCTTGGCAG
AGCGGGTTTGTAAACAGGCACTTTTGATCTTCAAGAGGTGCTGTTTCATGGGGGCAGCAGCTTCACAGTG
CTGCGTTGCACCAGCTTTGCACTGGAGAGCTGTGCCACCCTGACCGCCACTACTCCTTAGCTGACCTG
CTCAAGTACAACCTCTACCTGCCCTTCTTCTTCTCGGGCCCATCATGACCTTTGATCGCTTCCATGCT
CAGGTGAGCCAGGTGGAGCCAGTGAGACGCGAGGGTGAGCTGTGGCACATCCGAGCCAGGCAGGCCTA
AGCGTGGTGGCCATCATGGCCGTCGACATCTTCTTTCACTTCTTCTACATCCTCACTATCCCCAGCGAC
CTCAAGTTCGCCAACCGCTCCAGACAGTGCCCTCGCTGGCCTAGCCTATTCAAACCTGGTGTATGAC
TGGGTGAAGGCGGCCGCTCCTCTTTGGTGTGTCAACACTGTGGCATGCCTCGACCACTGGACCCACCC
CAGCCTCCCAAGTGCATACCGCACTCTACGCTTTGCGAAACGCACCTTGACCGTGGCATCAACGAC
TGGCTTTGCAATATGTGTATAACCACATTGGTGGGAGCATTCCGCTGTGATCCAGAGCTGGCAGCC
ACAGTGGCCACATTTGCCATCACCACACTGTGGCTTGGCCTTGTGACATTGTCTACCTGTGGTCATTCT
TTAACTGCTTTGGCCTCAACTTTGAGCTCTGGATGCAAAACTGGCAGAGTGGGGGCCCTAGCACGA
ATTGAGGCCTCTCTGTGAGTGCAGATGTCCCGTAGGGTCCGGGCCCTGTTTGGAGCCATGAACCTCTGG
GCCATCATCATGTACAACCTTGTGAGCCTGAACAGCCTCAAATTCACAGAGCTGGTTGCCCGGCGCTG
CTACTCACAGGGTTCCCCAGACCAGCTGTCCATCCTGTTTGTACCTACTGTGGCGTCCAGCTGGTA
AAGGAGCGTGAGCGAACCTTGGCACTGGAGGAGGAGCAGAAGCAGGACAAAGAGAAGCCGGAGTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**ACCN:** NM\_020707

**Insert Size:** 1515 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

|                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol> |
| <b>RefSeq:</b>                | <u>NM_020707.3</u>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>RefSeq Size:</b>           | 1875 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>RefSeq ORF:</b>            | 1515 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Locus ID:</b>              | 57467                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>UniProt ID:</b>            | <u>Q9HCP6</u>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Cytogenetics:</b>          | 3p22.1                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Protein Families:</b>      | Transmembrane                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>MW:</b>                    | 56.7 kDa                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Gene Summary:</b>          | <p>Negatively regulates N-terminal palmitoylation of SHH by HHAT/SKN.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) represents the longer transcript and encodes the functional protein.</p>                                                                                                                                                                                                                                                                                               |