

Product datasheet for **SC304826**

KIAA1328 (NM_020776) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	KIAA1328 (NM_020776) Human Untagged Clone
Tag:	Tag Free
Symbol:	KIAA1328
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC304826 representing NM_020776.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCGGACGTGGCGGGCCCCCTCCGCCCAAGTCCCGCGGCTTCTGGAGCCGGGACTTTTCTGATGAA
GAACAATCAGTAGTATACGTTCCAGGAATTTCTGCTGAAGGAAATGTCAGATCAAGACACAAGCTGATG
AGTCCAAAAGCTGATGTTAAACTTAAGACTTCCAGGGTGACTGATGCTCAATCTCCATGGAGTCCCTTA
AAAGGCACAGGAGATTCACTAGATGAACAGAATTCCTGCAGGGGAGAAAATAAGAGTGCATCATTGAAG
GATTTATGTCTTGAAGACAAAAGACGCATTGCAAACTTAATTAAGAAGTGGCCAGAGTAAGTGAGGAA
AAGGAAGTGACAGAGGAAAGACTGAAAGCTGAGCAGGAGTCATTTGAGAAGAAGTCAGGCAAGTTGGA
GAACAGAATGAAGTATCAAAAGAAAGGGAAGCTCTCAGCTACAGTATAGAGAATGCCAAGAAGCTT
CTAAGCCTGTATCAGAAATTTTATCAGAACACAGGAGAAGCTCACCATGTCTCTCTCAGAACTTGGT
GCTGCTAGAATGCAGGAACAGCAGGTATCCAGTAGAAAAAGCACTCTCCAGTGTTCTGTGGAAGTGA
GATGGTTCTACTTGAGCATAGCCAGACCAGACCTACTATCAAAACCAAGCAAGACCTAAGTCTGCA
GTCCAGGATTGAGTTCAGAATCCCTTATAGCATTAGGAATAATTCTTTGAAACAGTAACCCCTTCAT
CATCCCAAAGATGATCTAGATAAGATACCATCAGAGACCAACATGTAATTGTGAATCTCCAGGGAGA
AAACCTGCAGTCCCAACAGAGAAAATGCCACAAGAAGATTGCACATGAAGGAATGTCCACATCTTAAG
CCTACTCCTAGTCAATGCTGTGGTCATAGACTTGTGCAGATCGTGTTCATGACAGCCATCCTACAAAC
ATGACCCCTCAACATCCTAAGACACATCCAGAATCATGCAGTTATTGTGCGCTTTCTTGGGCATCTCTG
GTGCATGGTGGTGGGCACTGCAACCCATTGAACTTTGAAAAGCAGATCTCAGAAGATAGAAAGCAG
CAACTGATGCTTCAGAAAATGGAAGTGGAAATGAAAAGGAGCGCCTTCAGCATCTGCTGGCCAGCAG
GAGACAAAAGCTTCTTCTAAACAGCAGCAGCTTACCAGTCTCGACTGGATTACAATTGTTTATTGAAG
TCAAACGTGTGATGGCTGGCTGCTTGAACATCATCATCTATTAAGAACCAAGACCCCAAGCAAGT
GGAGAGAAATAGGAAGGAGAGGAAGACAGTTGGGTTTCATTGCGATATGAAAGATGATGCCAGTGGTCA
TGTCAAAGAAAGATACATGTAGACCCCAAGAGGGACAGTGACAGGAGTTAGAAAAGATGCGTCTACA
TCTCCTATGCCAACAGGAAGCCTAAAGGATTTTGTACCACAGCCTCACCATCATTACAGCACACCACC
TCCCGGTATGAGACATCTTTGTTGGATTTGGTTCAGTCTCTGAGCCCAAGCTCTGCGCCCAACCTCAG
CGCTATCCCTCCAGAGAAGCTGGGGCTGGAATCATGGTACTTTCCGACTCAGTCTCTAAATCAACC
CGGAAGAAGATGGGGATGCACAGAACCCTGAAGAGTTGGAGGAGAATCAGATTCTGGAAGATATATT
TTCATTGA
ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_020776

Insert Size: 1734 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_020776.1</u>
RefSeq Size:	4875 bp
RefSeq ORF:	1734 bp
Locus ID:	57536
UniProt ID:	<u>Q86T90</u>
Cytogenetics:	18q12.2
MW:	65.4 kDa
Gene Summary:	Competes with SMC1 for binding to SMC3. May affect the availability of SMC3 to engage in the formation of multimeric protein complexes.[UniProtKB/Swiss-Prot Function]