

## Product datasheet for **SC304850**

### **KIAA1467 (FAM234B) (NM\_020853) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	KIAA1467 (FAM234B) (NM_020853) Human Untagged Clone
Tag:	Tag Free
Symbol:	FAM234B
Synonyms:	KIAA1467
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

**Fully Sequenced ORF:** >SC304850 representing NM\_020853.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGGCGACCGTGCTGCCAGGGCGCTCAAGCTGCCGGGAAGAAGAGCCAGACCTAGGGGAGTATGAT
CCACTTACCCAGGCTGACAGTGATGAGAGCGAAGACGATCTGGTGCTTAACCTGCAGAAGAATGGAGGG
GTCAAAAATGGGAAGAGTCCTTTGGGAGAAGCGCCAGAACCCGACTCAGATGCTGAGTTGCAGAGGCT
GCAAAGCCACATCTTTCAGAAGTCACCACGGAGGGCTACCCCTCAGAACCCTTGGGGCCCTGGAACAG
AAGGCGGCTCCTCCCTGGTGTATATGTGCGCACGTCTGTCTTCTGCTGACTTTGGGGATCTCGATG
ATCCTGGTGTCTGTGTCTTCTGATCCCCTGTCTCCAGAGATCTGCACAGCACCTGGAGCCGC
CACTTGGGCTCCCAGGGAGGTGGGGACCTGTCTCCATTGGAATTGGCTGATGTGAATGGAGATGGCTG
CGTGATGTCTTCTCCTTTGTGATGTCAAGGAACGGGAGTGCAGTAGGTGTCTCAAGACCAGCTGCT
AATCTTGATGCCTTTGGGGATGAATGGCAGCACACTGTGGTCTAGTCTTCTCCCTGAGGAGGCTCGA
GATATCATGTTTGGAGCTGATGCCAGGAAGCTTGGCTGAAACCATCTGCCTTGTGACAGGGACACAC
AAGATGCTCAGCGCATTCAATGCAACGTGAGGAAAGCCATTTGGACTTTAAACCCAAACTACTTGTCC
AACGGTACCTTGGCTGCCCCAGTTGTGGTACTGCCAGACTTGGATGAAGACGGTGTTCGAGACCTTGTG
GTTCTGGCCATTGGGGAATTGCAGCCAGATCTGTGCTTCTGCTGGTGTCTGGCCGACCGGAAATCCA
GTGGGTCGACCTGTGAAGTACAACATCGTTGGAGTTGGGAATCTGATTGGTCTCAGGTTACATCACC
ACAAATGGGGCTGTCTACATCCTGTTGGCTTTGGAAATATAACAAGTGTGCGACTGCGGGACATTTTT
GTTCAGGCCAAAATCGAGACAGCTCACACCTTCTCTGCAGATAGAAGAGCCAGAATGGGAAAAGCGA
AGATCCATCAACCTGTCTGAGCTCATTGATGTTTACAGTGATGGTGTGAACTACTCCAGATGGTGAAG
GCACCAGATTTCAACTGCAGCAACCTTCTGATTACAACCAGACAAGCCTTGTGCTGCTTCGGGGGCAA
AATCTGACACCTTACTGGGCATTGAGACTTCAAGGCCTGCGCAGCCAGCCTACTCCTGGATATTTCACT
GATGATCAGACATTAGACTTCCTTCTGCAGATACAGGATGGAGTTGGGATGAAAAAGATGATGGTTGTG
GATGGTACTCTGGCTCCATTGTTGGAGTTACCGTGCTCCGTGTACATGAAAGAAACGCCAGCCACC
TCAGCAGTTACTTCAGACCAGAAGTCTGTCTTCTTCTGGGCCGAAGGGCTGTCAGCTGCATCTCCC
AATTCCGATATCATCCTAGGAAGTACAGCCGCCAGCCTTACCACCTTACCTCCTGCATCCTGCGTTC
CCCTCCATCCTTCTGGATCTGGCCAACACCACCGGCACAGTGACGGCTTCCAGAGTTGGAATTAACGAC
CTCTGAAAGATGCCTTTTATGTTACCAGGACAACAGGGCCAAGCTCCGAAGGCCATCCAGCAGCCCTG
GTGGTCAGCAAGCTTAGTCTACGGTGGGCACTAATGGAGGGCCAGATGGCTCAGCTACAGGAGTCCACC
CCAAAATTGGCCGTGGGAGCTGCGAAGATTTCTCTCTAGGATAAAGTTTGTGAAGCTCCCTACGAG
ATCTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** Sgfl-Mlul

**Plasmid Map:**


**ACCN:** NM\_020853

**Insert Size:** 1869 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_020853.1](#)

**RefSeq Size:** 4725 bp

**RefSeq ORF:** 1869 bp

**Locus ID:** 57613

**UniProt ID:** [A2RU67](#)

**Cytogenetics:** 12p13.1  
**Protein Families:** Druggable Genome, Transmembrane  
**MW:** 67 kDa