

## Product datasheet for **SC304914**

### CDH19 (NM\_021153) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	CDH19 (NM_021153) Human Untagged Clone
Tag:	Tag Free
Symbol:	CDH19
Synonyms:	CDH7; CDH7L2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

**Fully Sequenced ORF:** >SC304914 representing NM\_021153.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGAACGTGTTATTTACTGCTGCGTTTTATGTTGGGAATTCCTCTCCTATGGCCTTGCTTGGAGCAACA
GAAAACCTCAAACAAGAAAGTCAAGCAGCCAGTGCGATCTCATTGAGAGTGAAGCGTGGCTGGGTG
TGGAAACCAATTTTTGTACCAGAGGAAATGAATACGACTAGTCATCACATCGGCCAGCTAAGATCTGAT
TTAGACAATGGAACAATTCTTTCCAGTACAAGCTTTTGGGAGCTGGAGCTGGAAGTACTTTTATCATT
GATGAAAGAACAGGTGACATATATGCCATACAGAAGCTTGATAGAGAGGAGCGATCCCTCTACATCTTA
AGAGCCCAGGTAATAGACATCGTACTGGAAGGGCTGTGGAACCTGAGTCTGAGTTTGCATCAAAGTT
TCGGATATCAATGACAATGAACCAAAATTCCTAGATGAACCTTATGAGGCCATTGTACCAGAGATGTCT
CCAGAAGAACATTAGTTATCCAGGTGACAGCAAGTGATGCTGACGATCCCTCAAGTGGTAATAATGCT
CGTCTCCTCTACAGCTTACTTCAAGGCCAGCCATATTTTTCTGTTGAACCAACAACAGGAGTCATAAGA
ATATCTTCTAAAATGGATAGAGAACTGCAAGATGAGTATTGGGTAAATCATTCAAGCCAAGGACATGATT
GGTCAGCCAGGAGCGTTGTCTGGAACAACAAGTATTAAATTAACCTTTTCAAGTGTAAATGACAATAAG
CCTATATTTAAAGAAAAGTTTATACCGCTTGACTGTCTCTGAATCTGCACCCACTGGGACTTCTATAGGA
ACAATCATGGCATATGATAATGACATAGGAGAGAATGCAGAAATGGATTACAGCATTGAAGAGGATGAT
TCGCAACATTTGACATTACTAATCATGAACTCAAGAAGGAATAGTTATATTAAGAAAAGAAAGTG
GATTTTGAGCACCAGAACCCTACGGTATTAGAGCAAAAGTTAAAAACCATCATGTTCTGAGCAGCTC
ATGAAGTACCACACTGAGGCTCCACCCTTTCATTAAGATCCAGGTGGAAGATGTTGATGAGCCTCCT
CTTTTCTCCTCCATATATGATTTGAAGTTTTTGAAGAAACCCACAGGGATCATTTGTAGGCGTG
GTGCTGCCACAGACCCAGACAATAGGAAATCTCCTATCAGGTATTCTATTACTAGGAGCAAAGTGTTC
AATATCAATGATAATGGTACAATCACTACAAGTAACTCACTGGATCGTGAATCAGTGCTTGGTACAAC
CTAAGTATTACAGCCACAGAAAATACAATATAGAACAGATCTCTTCGATCCCACTGTATGTGCAAGTT
CTTAACATCAATGATCATGCTCCTGAGTCTCTCAATACTATGAGACTTATGTTTGTGAAAATGCAGGC
TCTGGTCAGGTAATTCAGACTATCAGTGCAGTGGATAGAGATGAATCCATAGAAGAGCACCATTTTAC
TTAATCTATCTGTAGAAGACACTAACAAATCAAGTTTTACAATCATAGATAATCAAGATAACACAGCT
GTCATTTTACTAATAGAAGTGGTTTTAACCTTCAAGAAGAACCTGTCTTCTACATCTCCATCTTAATT
GCCGACAATGGAATCCCCTCACTTACAAGTACAAACACCCTTACCATCCATGTCTGTGACTGTGGTGAC
AGTGGGAGCACACAGACCTGCCAGTACCAGGAGCTTGTGCTTTCATGGGATTCAAGACAGAAGTCATC
ATTGCTATTCTCATTGTCATTATGATCATATTTGGGTTTTTTTTTGGCTTTGGGTTTTAAACAACGG
AGAAAACAGATTCTATTTCTGAGAAAAGTGAAGATTTTCAAGAGAAATATATTCCAATATGATGATGAA
GGGGTGGAGAAGAAGATACAGAGGCCTTTGATATAGCAGAGCTGAGGAGTAGTACCATAATGCGGGAA
CGCAAGACTCGGAAAACCAAGCGCTGAGATCAGGAGCCTATACAGGCAGTCTTTGCAAGTTGGCCCC
GACAGTGCATATTCAGGAAATTCATCTGAAAAGCTCGAAGAAGCTAATACTGATCCGTGTGCCCTT
CCTTTTGATTCCCTCCAGACCTACGCTTTTGGGGAACAGGGTCATTAGCTGGATCCCTGAGCTCCTTA
GAATCAGCAGTCTCTGATCAGGATGAAAGCTATGATTACCTTAATGAGTTGGGACCTCGCTTTAAAGA
TTAGCATGCATGTTTGGTCTGCAGTGCAGTCAAATAATTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**Plasmid Map:**


**ACCN:** NM\_021153

**Insert Size:** 2319 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_021153.3](#)

**RefSeq Size:** 6241 bp

**RefSeq ORF:** 2319 bp

**Locus ID:** 28513

**UniProt ID:** [Q9H159](#)

**Cytogenetics:** 18q22.1

**Protein Families:** Druggable Genome, Transmembrane

**MW:** 87 kDa

**Gene Summary:** This gene is one of three related type II cadherin genes situated in a cluster on chromosome 18. The encoded protein is a calcium dependent cell-cell adhesion glycoprotein containing five extracellular cadherin repeats. Loss of cadherins may be associated with cancer formation. Alternative splicing results in multiple transcript variants for this gene. [provided by RefSeq, Aug 2012]

Transcript Variant: This variant (1) represents the longest transcript and encodes the longer isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.