

Product datasheet for **SC304930**

DMRT3 (NM_021240) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DMRT3 (NM_021240) Human Untagged Clone
Tag:	Tag Free
Symbol:	DMRT3
Synonyms:	DMRTA3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC304930 representing NM_021240.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAACGGCTACGGCTCCCCCTACCTGTACATGGGCGGCCCGGTGTCGCAGCCGCCAGGGCGCCCTG
CAGCGCACGCCCAAGTGCAGCGCTGCCGCAACCATGGCGTCTGTCTGGCTCAAGGGCCACAAGCGT
TACTGCCGCTTCAAGGACTGCACCTGCGAGAAGTGCATCCTCATCATCGAGCGGCAGCGGGTCATGGCT
GCGCAGGTGGCGCTGCGCCGGCAGCAGGCCAACGAGAGCTTGGAGAGCCTCATCCCCGACTCGTGCGC
GCTCTGCCAGGGCCCCCGCCGGGGGACGCCGTGCGCGCCCCGAGCCGCCCGCAGCCTCTCAGCCG
TCGCAGCCGCAGCCGCCGCGCCTGCTGCCGAGTTGGCCGCGCCGCCGCGCTGCGTTGGACTGCCGAG
CCGCAGCCCGGGCTCTGCAGGGCAGCTCGCCAAAGCCAGATTTGACTGAAGAACGACTTGGAGACGGC
AAGTCGGCAGACAATACAGAGGTCTTCAGTGACAAAGACACTGACCAGAGGAGTCCCCAGATGTGGCA
AAGAGTAAGGGCTGCTTACCCCTGAGAGCCCTGAGATAGTGTCCGTGGAGGAAGGGGGATACGCTGTC
CAGAAAAACGGAGCAACCCCGAGAGCCGCCCTGACAGCCCAAGTGTACGCGGAGCAGAATCACCTC
CTGATTGAGGGCCCTCGGGGACTGTTTCTGTGCCCTTACGTTGAAAGCCAACAGACCCGCGCTTGAA
GTGTTAAAAAAGATATCCCCAACCAGAAGCCAACGGTGCTTGGAGCTATCCTCAAGGGGTGTGGCGGG
GACCTGGTGAGCGCCGTGGAAGTCTTCTGTCCAGCCGATCCTCAGTCACGGGAGCAGAGCGAATTCC
GCAGAACCTGAGAGTCTAGCGTTGCCCTCCAATGGGCACATCTTTGAACACACCTTGAGCTCCTACCCC
ATCTCGTCTTCCAAATGGTCTGTGGGATCAGCCTTTCGAGTCCCAGACACGTTGAGGTTTTCTGCCGAC
TCTAGCAACGTTGCCCCAGTCCCTTGGCTGGGCTCTGCAGCCCCCTTTCCCCAGCCACCCCGGTAC
CCGCTGATGCTGAGGAATACTTTGGCGAGAAGCCAGTCGAGCCCCCTTTTGGCCAATGATGTCACCCTG
TGAACACCATGACGCTGCAGCAGCAGTATCAGCTGAGGTCCCAGTATGTCAGTCTTTCCCCAGTAAC
TCTACCAGCGTCTTCAAGAGCTCGCCCGTCTTCTGCCGCGCCACGGAAGACCCTCGGATTTCCATC
CCTGATGATGGGTGCCATTTGTGTCAAAGCAGTCCATTTACACCGAGGACGACTATGACGAGAGGTCT
GACTCCTCAGACTCTAGAACACTCAACACATCATCTTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites:

SgfI-MluI

Plasmid Map:



ACCN: NM_021240

Insert Size: 1419 bp

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_021240.3</u>
RefSeq Size:	2200 bp
RefSeq ORF:	1419 bp
Locus ID:	58524
UniProt ID:	<u>Q9NQL9</u>
Cytogenetics:	9p24.3
MW:	51.2 kDa
Gene Summary:	<p>Probable transcription factor that plays a role in configuring the spinal circuits controlling stride in vertebrates. Involved in neuronal specification within specific subdivision of spinal cord neurons and in the development of a coordinated locomotor network controlling limb movements. May regulate transcription during sexual development (By similarity). [UniProtKB/Swiss-Prot Function]</p>