

Product datasheet for SC304935

IFNA17 (NM_021268) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	IFNA17 (NM_021268) Human Untagged Clone
Tag:	Tag Free
Symbol:	IFNA17
Synonyms:	IFN-alpha1; IFNA; INFA; LEIF2C1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC304935 representing NM_021268. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCCCTGTCCTTTTCTTACTGATGGCCGTGCTGGTGCTCAGCTACAAATCCATCTGTTCTCTAGGC
TGTGATCTGCCTCAGACCCACAGCCTGGGTAATAGGAGGGCCTTGATACTCCTGGCACAAATGGGAAGA
ATCTCTCCTTTCTCCTGCCTGAAGGACAGACATGACTTTGGACTTCCCAGGAGGAGTTTGATGGCAAC
CAGTTCAGAAGACTCAAGCCATCTCTGTCTCCATGAGATGATCCAGCAGACCTCAATCTCTTCAGC
ACAGAGGACTCATCTGCTGCTTGGGAACAGAGCCTCCTAGAAAAATTTTCCACTGAACCTTACCAGCAA
CTGAATAACCTGGAAGCATGTGTGATACAGGAGGTTGGGATGGAAGAGACTCCCCTGATGAATGAGGAC
TCCATCCTGGCTGTGAGGAAATACTTCCAAAGAATCACTCTTTATCTAACAGAGAAGAAATACAGCCCT
TGTGCCTGGGAGGTTGTGAGCAGAAAATCATGAGATCTCTCTTTTTTCAAAACTTGCAAAAAATA
TTAAGGAGGAAGGATTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-MluI



[View online »](#)

Plasmid Map:


ACCN:	NM_021268
Insert Size:	570 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq:	NM_021268.2
RefSeq Size:	980 bp
RefSeq ORF:	570 bp
Locus ID:	3451
UniProt ID:	P01571
Cytogenetics:	9p21.3
Protein Families:	Druggable Genome, Secreted Protein
Protein Pathways:	Antigen processing and presentation, Autoimmune thyroid disease, Cytokine-cytokine receptor interaction, Cytosolic DNA-sensing pathway, Jak-STAT signaling pathway, Natural killer cell mediated cytotoxicity, Regulation of autophagy, RIG-I-like receptor signaling pathway, Toll-like receptor signaling pathway
MW:	21.7 kDa
Gene Summary:	Produced by macrophages, IFN-alpha have antiviral activities. Interferon stimulates the production of two enzymes: a protein kinase and an oligoadenylate synthetase. [UniProtKB/Swiss-Prot Function]