

Product datasheet for **SC304937**

ICEBERG (CARD18) (NM_021571) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	ICEBERG (CARD18) (NM_021571) Human Untagged Clone
Tag:	Tag Free
Symbol:	CARD18
Synonyms:	ICEBERG; pseudo-ICE; UNQ5804
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC304937 representing NM_021571. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCTGACCAACTCTTGCGTAAAAAGAGAAGAATTTTATCCATTCAGTGGGTGCAGGCACAATAAAT
GCCTTGCTGGATTGCCTATTAGAGGATGAAGTTATTAGCCAGGAAGACATGAACAAAGTGAGAGATGAA
AATGACACTGTCATGGATAAGGCTCGAGTCTTGATTGACCTTGTTACTGGAAAAGGACCCAAGTCTTGC
TGCAAATTTATCAAGCATCTCTGTGAAGAAGACCCTCACTTGCCTCAAAGATGGGTTTGCACTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-Mlul



[View online »](#)

Plasmid Map:


ACCN: NM_021571

Insert Size: 273 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_021571.3
RefSeq Size:	670 bp
RefSeq ORF:	273 bp
Locus ID:	59082
UniProt ID:	P57730
Cytogenetics:	11q22.3
Protein Families:	Druggable Genome
Protein Pathways:	NOD-like receptor signaling pathway
MW:	10.1 kDa
Gene Summary:	Inhibits generation of IL-1-beta by interacting with caspase-1 and preventing its association with RIP2. Down-regulates the release of IL1B.[UniProtKB/Swiss-Prot Function]