

Product datasheet for **SC304971**

MUPCDH (CDHR5) (NM_021924) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MUPCDH (CDHR5) (NM_021924) Human Untagged Clone
Tag:	Tag Free
Symbol:	CDHR5
Synonyms:	MLPCDH; MU-PCDH; MUCDHL; MUPCDH
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC304971 representing NM_021924.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGGGCTTTGGGCCCTGCTGTGGCCTCCCCTGCTGTTACCCGGGCTGCTCGTCCGACCCCCGGGGACC
ATGGCCAGGCCAGTACTGCTCTGTGAACAAGGACATCTTTGAAGTAGAGGAGAACAACAATGTACCC
GAGCCCGTGGTGGACATCCACGTCCCGGAGGGCCAGGAGGTGACCCTCGGAGCCTTGTCCACCCCTTT
GCATTTCCGGATCCAGGAAACCAGCTGTTTCTCAACGTGACTCCTGATTACGAGGAGAAGTCACTGCTT
GAGGCTCAGTGTGTGTGAGAGCGGAGGCACATTGGTGACCCAGCTAAGGGTGTTCGTGTGAGTGTG
GACGTCAATGACAATGCCCCGAATTCCTTTAAGACCAAGGAGATAAGGGTGGAGGAGACACGAAA
GTGAACCTCACCGTCACTCCCGAGACGCAACTGCAGGCTGAGGACCGGACAAGGACGACATTCTGTTT
TACACCTCCAGGAAATGACAGCAGGTGCCAGTACTTCTCCCTGGTGTGTAACCGTCCCGCC
CTGAGGCTGGACCGCCCTGGACTTCTACGAGCGGCCGAACATGACCTTCTGGTGTGTTGCGGGAC
ACTCCGGGGGAGAATGTGAACCCAGCCACACTGCCACCGCCACACTAGTGTGAACGTGGTCCCGCC
GACCTGCCGGCCCGTGGTTCTGCCCTGCACCTTCTCAGATGGCTACGCTGCATTCAAGCTCAGTAC
CACGGGGCTGTCCCCACGGGGCACATACTGCCATCTCCCCTCGTCTGCTGCCGACCCATCTACGCT
GAGGACGGAGACCGCGGCATCAACCAGCCCATCTACAGCATCTTTAGGGGAAACGTGAATGGTACA
TTCATCATCCACCAGACTCGGGCAACCTCACCGTGGCCAGGAGTGTCCCAGCCCCATGACCTTCTT
CTGCTGGTGAAGGGCCACAGGCCGACCTTGCCCGCTACTCAGTGACCCAGGTACCGTGGAGGCTGTG
GCTGCGGCCGGGAGCCCGCCCGCTTCCCCAGAGACTGTATCGTGGCACCGTGGCGGTGGCGCTGGA
GGGGCGTGTGGTCAAGGATGCAGTGCCTTCTCAGCCTCTGAGGATCCAGGCTCAGGACCCGGAG
TTCTCGGACCTCAACTCGGCCATCATATCGAATTACCAACCACTCACACTTCCGGATGGAGGGAGAG
GTTGTGCTGACCACCACACACTGGCACAGGCGGAGCCTTCTACGAGAGGTTGAGGCCACAACACG
GTGACCTCTGGCACCGCAACCACAGTATTGAGATACAAGTTTCCGAACAGGAGCCCCCTCCACAGAT
GTCCCCCATCCCCAGAGGCTGGAGGAACAAGTGGGCCCTGGACCAGCACCCTCCGAGGTCCCCAGA
CCCCCTGAGCCCTCCAGGGACCTCCACGACCAGCTCTGGGGGAGGCACAGGCCCTCATCCACCCTCT
GGCACAACCTCTGAGGCCACCAACCTCGTCCACACCCGGGGGGCCCCGGGTGCAGAAAACAGCACCTCC
CACCAACCAGCCACTCCCGTGGGGACACAGCACAGACCCCAAAGCCAGGAACCTCTCAGCCGATGCCC
CCCGGTGTGGAAACCAGCACCTCCACCAACCAGCCACACCCAGTGGGGGCACAGCACAGACCCAGAG
CCAGGAACCTCTCAGCCGATGCCCCCAAGTATGGGAACCAGCACCTCCACCAACCAGCCACACCCGGT
GGGGGCACAGCACAGACCCAGAGGCAGGAACCTCTCAGCCGATGCCCCCGGTATGGGAACCAGCAC
TCCACCAACCAACCACACCCGGTGGGGGCACAGCACAGACCCAGAGCCAGGAACCTCTCAGCCGATG
CCCCTCAGCAAGAGCACCCCATCTTCAGGTGGCGGCCCTCGGAGGACAAGCGCTTCTCGGTGGTGGAT
ATGGCGGCCCTGGCGGGGTGCTGGGTGCGCTGCTGCTGGCTCTCCTTGGCCTCGCCGTCTTGTGTC
CACAAGCACTATGGCCCCGGCTCAAGTGTGCTCTGGCAAAGCTCCGGAGCCCCAGCCCCAAGGCTTT
GACAACCAGGCGTTCTCCCTGACCACAAGGCCAACTGGGCGCCCGTCCCAGCCCCACGCACGACCCC
AAGCCCCGCGGAGGCACCGATGCCCGCAGAGCCCGACCCCGGCCCTGCCTCCCAGGCGGTGCCCT
GAGCCCCCGCAGCGGCCGAGCTGGCGGAAGCCACGCGGTTGAGGTCATCTGACCAAGGAGCGG
CGGCCGGAGGGCGGTACAAGGCTGTCTGGTTTGGCGAGGACATCGGGACGGAGGCAGAGCTGGTCTT
CTCAACGCGCCACCCTGGACGTGGATGGCGCCAGTACTCCGGCAGCGCGATGAGGGCAGGGCGCG
GGGAGGGGTGGGGTCCCTACGATGCGCCCGTGGTGTGACTCTACATCTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_021924

Insert Size: 2538 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_021924.4](#)

RefSeq Size: 3489 bp

RefSeq ORF: 2538 bp

Locus ID: 53841

UniProt ID: [Q9HBB8](#)

Cytogenetics: 11p15.5

Protein Families: Transmembrane

MW: 88.2 kDa

Gene Summary: This gene is a novel mucin-like gene that is a member of the cadherin superfamily. While encoding nonpolymorphic tandem repeats rich in proline, serine and threonine similar to mucin proteins, the gene also contains sequence encoding calcium-binding motifs found in all cadherins. The role of the hybrid extracellular region and the specific function of this protein have not yet been determined. Alternatively spliced transcript variants encoding different isoforms have been described. [provided by RefSeq, Jan 2010]
Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (1, also known as MUCDHLFL). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.