

Product datasheet for SC305017

EDDM3B (NM_022360) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	EDDM3B (NM_022360) Human Untagged Clone
Tag:	Tag Free
Symbol:	EDDM3B
Synonyms:	EP3B; FAM12B; HE3-BETA; HE3B; RAM2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC305017 representing NM_022360. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
 ATGGCATCGTCTCTAAAGATCTGGGACACTCTTGGCCCTACTTTGCATCCTATGCACACTGCTTGTG
 CAGAGCAAAGAAGTTTCTTGGAGAGAATTCATGAAACAGCACTACTTAAGTCCAAGTCGAGAATTCAGA
 GAGTACAAATGTGATGTCCTCATGAGAGAAAATGAAGCTCTGAAAGACAAGAGCTCTCACATGTTTATC
 TATATCTCATGGTACAAAATCGAGCATATATGCACTAGTGACAACTGGATGGATCGCTTCCGAAATGCA
 TATGTATGGGTCCAGAATCCTCTCAAAGTACTCAAGTGTACCAGGAGAATTCAAAAATAGCTACACA
 GAGAGCAGGAGCTTCACTACATTGAATTCCATTGTAGCATGGACGGGTATGTTGATAGCATAGAAGAC
 CTAAAGATGGTAGAACCTATCGGCAACTAG
 ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Restriction Sites:	Sgfl-MluI
ACCN:	NM_022360
Insert Size:	444 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).


[View online »](#)

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_022360.4</u>
RefSeq Size:	897 bp
RefSeq ORF:	444 bp
Locus ID:	64184
UniProt ID:	<u>P56851</u>
Cytogenetics:	14q11.2
Protein Families:	Secreted Protein, Transmembrane
MW:	17.6 kDa
Gene Summary:	<p>Testicular sperm are morphologically differentiated but are not progressively motile nor able to fertilize an egg. Post-testicular maturation requires exposure of spermatozoa to the microenvironment of the epididymal lumen. Spermatozoa undergo extensive changes in the epididymis, including enzymatic modifications, loss of pre-existing components and addition of new glycoproteins from epididymal secretions. These modifying proteins and enzymes are synthesized by epithelial cells lining the epididymal duct and secreted apically into the lumen, where they come into contact with, and may be absorbed onto, the sperm membranes. The proteins encoded by the genes in this cluster are synthesized and secreted by epididymal epithelial cells. [provided by RefSeq, Jul 2008]</p>