

## Product datasheet for **SC305057**

### **NIBAN2 (NM\_022833) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	NIBAN2 (NM_022833) Human Untagged Clone
Tag:	Tag Free
Symbol:	NIBAN2
Synonyms:	bA356B19.6; C9orf88; FAM129B; MEG-3; MINERVA; OC58
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC305057 representing NM\_022833.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGGGGACGTGCTGTCCACGCACCTGGACGACGCCCGGCCAGCACATCGCAGAAAAACCGGGAAG
ATCCTGACGGAGTTCCTCCAGTTCTATGAAGACCAGTATGGCGTGGCTCTCTTCAACAGCATGCCCAT
GAGATTGAGGGCACGGGGCTGCCGACGGCCAGCTGCTCTGGCGCAAGGTGCCACTGGACGAGCGCATC
GTCTTCTCGGGAACTCTTCCAGCACAGGAGGACAGCAAGAAGTGGAGAAACCGCTTACGCTCGTG
CCCCAACAACGAGGGCTGGTGTCTACGAAAAAAGCGGCTATGAGCGGCAGTCCCACCACGAGCC
GTCATCAACAGTGCAGGCTACAAAATCCTCACGTCCTGGACCAATACCTGGAGCTCATTGGCAACTCC
TTACCAGGACCACGGCAAAGTCGGGCAGTGCCCCATCCTCAAGTCCCCACACAGTCCCCTCATC
CTCTGGCATCCTATGCGCGTCACTACTTCTGCATGATGACAGAAGCCGAGCAGGACAAGTGGCAG
GCTGTGCTGCAGGACTGCATCCGGCACTGCAACAATGGAATCCCTGAGGACTCCAAGGTAGAGGGCCCT
GGTTTACAGATGCCATCCGCATGTACCGACAGTCCAAGGAGCTGTACGGCACCTGGGAGATGCTGTGT
GGGAACGAGGTGCAGATCCTGAGCAACCTGGTGTGAGGAGCTGGGCCCTGAGCTGAAGGCAGAGCTC
GGCCCGCGGTGAAGGGGAAACCGCAGGAGCGGCAGCGGAGTGGATCCAGATCTCGACGCCGTGTAC
CACATGGTGTACGAGCAGGCCAAGGCGCGCTTCGAGGAGGTGCTGTCCAAGGTGCAGCAGGTGCAGCCG
GCCATGCAGGCCGTATCCGAATGACATGGACCAAAATACACCTCCAAGGAGCACCTTCCAGCAAG
ATCCGAGCCTTATCCTCCCAAGGCAGAGGTGTGCGTGCAGCAACATGTCCAGCCCTACATCCCATCC
ATCCTGGAGGCCCTGATGGTCCCACCAGCCAGGGCTTCACTGAGGTGCGAGATGCTTCTTCAAGGAG
GTCACGGACATGAACCTGAACGTATCAACGAGGGCGCATTGACAAGCTGGGCGAGTACATGGAGAAG
CTGTCCCGGCTGGCGTACCACCCCTGAAGATGCAGAGCTGCTATGAGAAGATGGAGTCCGCTGCCACTG
GACGGGCTGCAGCAGCGATTTGATGTGTCCAGCAGTCCGTGTTCAAGCAGCGAGCCAGATCCACATG
CGGGAGCAAATGGACAATGCCGTGTATACGTTCCGAGACCTCCTGCACCAGGAGCTGGGGAAGGGGCC
ACCAAGGAGGAGCTGTGCAAGTCCATCCAGCGGTCTGGAGCGGTGCTGAAGAAATACGACTACGAC
AGCAGCTCTGTGCGGAAGAGTTCTTCCGGGAGGCGTGTGCAGATCAGCATCCCGTCTCTGCTCAAG
AAGCTGGCCCTACCTGCAAGTCGGAGCTGCCCGGTTCCAGGAGCTGATCTTCCGAGACTTTGCCAGG
TTCATCCTGGTGGAAAACACGTACGAGGAGGTGGTGTGCAGACCGTATGAAGGACATCCTGCAGGCT
GTGAAGGAGGCCCGGTGCAGAGGAAGCACAACCTTACCGGACAGCATGGTATGCACAACAGCGAC
CCCAACCTGCACCTGCTGGCCGAGGGCGCCCCATCGACTGGGGGAGGAGTACAGCAACAGCGGGGG
GGCGGCAGCCCCAGCCCCAGCCCCGGAGTACGCCACCTCTCGGAAAAGCGACGGCGCCCAAGCAG
GTGGTCTCTGTGGTCCAGGATGAGGAGGTGGGGTGCCTTTGAGGCTAGCCCTGAGTACCACCACCT
GCGTCCCCGACGGTGTCACTGAGATCCGAGGCCTGCTGGCCCAAGGTCTGCGGCCTGAGAGCCCCCA
CCAGCCGGCCCCCTGCTCAACGGGGCCCCGCTGGGGAGAGTCCCAGCCTAAGGCCGCCCCGAGGCC
TCTCGCCGCTGCCTCACCCCTCCAGCATCTCTGCCTGGAAAGGCTGTGGACCTTGGGCCCCCAAG
CCCAGCGACCAGGAGACTGGAGAGCAGGTGTCCAGCCCCAGCAGCCACCCCGCTCCACACCACCAC
GAGGACAGTGCAGGGTGCAGACTGAGTTCTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_022833  
**Insert Size:** 2241 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_022833.3</a>
<b>RefSeq Size:</b>	3991 bp
<b>RefSeq ORF:</b>	2241 bp
<b>Locus ID:</b>	64855
<b>UniProt ID:</b>	<a href="#">Q96TA1</a>
<b>Cytogenetics:</b>	9q34.11
<b>MW:</b>	84.1 kDa
<b>Gene Summary:</b>	May play a role in apoptosis suppression. May promote melanoma cell invasion in vitro. [UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (1).