

Product datasheet for SC305096

ERI3 (NM_024066) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ERI3 (NM_024066) Human Untagged Clone
Tag:	Tag Free
Symbol:	ERI3
Synonyms:	PINT1; PRNPIP
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC305096 representing NM_024066. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCC**CGCATCGCC**
 ATGGCGACAGCCTCTCCGCTGCTGACGGGGGGCGGGGGCGGCCCTGGGAAGGAGGGCTGGTCTCCTGG
 CCCCCCGCCCTCCCTTACTCTCCCTGGACTTGGATGGGCCCGAGTTGGGGGCAACCCCCGGGCAT
 TGGGGCTTCCAGCTCTCACAGAACCTTCAGCATCCCCAGCTGCCGGTCTTGGCATCTTCGAAGTAAGG
 AGAGTTTTAGATGCTTCTGGATGTTCAATGCTAGCACCTTACAGACTGGAGCAGCTCGATTTTCTTCG
 TATTTACTTTCAAGAGCAAGAAAAGTGTGGGCTCCCACTTATTTCTCCCTGTGGTGTCCGGAGTTC
 TGCTCCATATCCACCAGAAAGCTGGCGGCCACGGCTTTGGCGCATCCATGGCGGCAATGGTGTCTCTC
 CCTCCCCAGAGGTATCACTACTTTTAGTGCTGGACTTTGAGGCCACGTGCGACAAGCCACAGATTCAT
 CCTCAGGAAATCATCGAGTTCCTCATCCTAAAGCTAAATGGCCGGACCATGGAGATTGAGTCTACCTTT
 CACATGTATGTCCAGCCTGTAGTCCATCCACAGCTTACCCATTCTGTACAGAGCTACCCGGGATTATT
 CAAGCCATGGTGGATGGTCAGCCAAGCCTGCAGCAAGTGTGGAGAGGGTCGATGAATGGATGGCGAAG
 GAAGGCCCTCTAGATCCAAACGTCAAGTCAATTTTGTACCTGTGGAGACTGGGACTTAAAGTCATG
 CTCCCAGGCCAGTGCCAGTACTTGGGCTTGCCAGTGGCGGATTACTTCAAGCAGTGGATTAATCTGAAA
 AAGGCTTACAGCTTCGCCATGGGCTGCTGGCCCAAGAATGGACTTCTAGACATGAACAAGGCCTCAGC
 CTGCAACACATAGGCCGGCCCCACAGCGCATTGACGACTGCAAGAACATTGCCAACATCATGAAGACA
 CTCGCTATCGAGGCTTCATCTTCAAGCAGACATCGAAGCCGTTCTGA
 ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Restriction Sites:	SgfI-MluI
ACCN:	NM_024066
Insert Size:	1014 bp


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_024066.2</u>
RefSeq Size:	1781 bp
RefSeq ORF:	1014 bp
Locus ID:	79033
UniProt ID:	<u>O43414</u>
Cytogenetics:	1p34.1
MW:	37.2 kDa