

## Product datasheet for **SC305138**

### **C10orf68 (CCDC7) (NM\_024688) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	C10orf68 (CCDC7) (NM_024688) Human Untagged Clone
Tag:	Tag Free
Symbol:	CCDC7
Synonyms:	BioT2-A; BioT2-B; BioT2-C; C10orf68
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC305138 representing NM\_024688.  
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGACTGAAGAACAACCTTACCAAGCAGCAGAAAAATCTCAAGCTCATAATGAAGTACCAATGAAAGG
CTTGTAGTTGAGCATCAAGAATCATTGTCAAAAACCAAAATTACAATAAAGAAACAAGAACTTCTACA
GAGCAACCACTCACCCTCTGATAAAGAACCAAAATGAAAATCTTATACTTAGGCATCAAGACTCAATG
TCAAAATCAGAAATGCAAGTGAAGGAACAAGAACTCTCAAAGGGCAAAGAATTATTACTCATGATGAA
GAACCAGGCAAAAATCTGTGCTTGAACATCAAGATTCAGTGTCAAAACTGGAAATGCAAAATGAAAAA
ACCAAAAAACTTCTAGAGAGAAAAGACATAGTACTCATGATGAAGAATCAGGTGAAAATCCTATGCTT
AAACATCAAGATTCAGTGTCAAAAATCCAAGTGAATTAGAGATTCAGAAACTTCTGAAGGAGAAGGA
CGTAGCATACCAGATAAAAATTTCTATGTTTGTTCATCAAGATTCAGTGTCAAAACTCCAATGCAAGAA
AAGAAAAAATAACTCCTGGAAGGGAAAGCGTAATACTCGTATTGTAGTACCAATGAAAATGTGATT
TCTGTTCAAGATTCAAAGTCAAAACTCCAATGCAAGAAAAGAAACAATAAATTCTGGAGTGGAA
AGACACAAGACGTTTCCTTTAGAAATCAAAAAAAGGATATATCACTTGAACATCTGTTGCCTGAGGAG
AAAGTTTTATTATCAAGAAGCGAATCTCAAACGAAGAACTTCAAGCTAAAGTAACCTCAAGAAAAAAG
ACTGATATAGAAAGCTTGAGAGGTGCTTTGGGAAGACGCTTATTAATGATGAATTCAAGACACAGTCA
AAGAGTTTCCCTGGGCTGATATAGAACAATTGACAGATGCATTTGGAAGAGATACTAAAGGATGAA
TTCAAGACACGATCAAAGAGTCTCCCTGAGACTGATGAACGATTCATAGTACAACCTGAGAGAGGTACA
ATAATGATGCAATTAAGACGCAGTTAAGAGAAAAGGTTACCCTGAGACTGTCTAAAACACTTGAAA
GGTGTAAATGGAAAAGATATAATAAAGCATCTAATCAACATACAGTCAAAGAGTCATGGAGAGACTGAT
AAGAACTTCTTTGCCTATGCTACTGGAAGAGGTCTAATGAAGGAGTCAACCACGACACAATAAAAAGT
CACCCAGAGACTGATAAGGAATTCCTGGCAGATGCTATTGGAAGAGGTATAATAATAGGGCAATCACT
ACACAACCTGAAGAGTCACCGAGAGACTGATAAAGAACTCTTAAAAGATGCTATTGGAAGAGATATAATA
AAGGGACCAATCAGTGCACAATTAAGAGTCACCAGGAGACTGATGTAGAACCCTTGACAAATGCCATT
GGTTCAAGTAAAACAATTGGTGAGATAAAGACACAACCTAAGGACTCACTATGATGTGAACCTATTTAAA
AACAAAGATATGTCGGTACAACGTCAAGAAGGTATCTTTACTAGGAGCATTACACCAAGTAAATCCCA
ACAAAAGTGATCAATTTATCACCTTTGAAAATAAAGAGGAACTTATGAGTACTCTTACCCTATGTG
ACTGCACCTTCAAAGCAATTTATAGAACATATAGGGCTGGTCCAAGCTTTCTAAAGACATCCACCTT
CCTTTACTGAATCAACTCCCTTCAGGGCATTGAAAGTTGTTACCTTAAGCCAAAAACTATTGAATTT
ACTCTGCCTACTGTGACCAATACAGTTGAAAACCTACCTATAAAGTGTTACATGCTGCTGCCCGAAAA
TCTGTACCACACCATACTTTTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

**Plasmid Map:**


**ACCN:** NM\_024688

**Insert Size:** 1887 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_024688.2](#)

**RefSeq Size:** 2730 bp

**RefSeq ORF:** 1887 bp

**Locus ID:** 79741

**Cytogenetics:** 10p11.22

**MW:** 71.5 kDa

**Gene Summary:** May play a role in tumorigenesis.[UniProtKB/Swiss-Prot Function]  
Transcript Variant: This variant (5) differs in the 5' UTR, lacks a large portion of the 5' coding region, and initiates translation at an alternate start codon. This results in a protein (isoform b), which is distinct from isoform a.