

Product datasheet for **SC305155**

TUBAL3 (NM_024803) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	TUBAL3 (NM_024803) Human Untagged Clone
Tag:	Tag Free
Symbol:	TUBAL3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC305155 representing NM_024803. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAGGGAGTGCCTTTCCATCCACATCGGTCAAGCTGGCATCCAGATTGGGGACGCCTGCTGGAACTC
TATTGCCTGGAACATGGAATCCAGCCAAATGGCGTTGTTCTTGACACTCAACAGGATCAGCTGGAAAAT
GCAAAAATGGAGCACACAAAATGCATCTTTGCATACCTTCTTGTGAGACAAGAGCTGGGAAGCATGTG
CCTAGAGCACTCTTCGTGGACTTGGAGCCAACCTGTTATAGATGGGATCCGGACGGCCAGCACCGTTCA
CTCTTCCACCCCGAGCAGCTCCTTAGCGGAAAGGAGGATGCTGCTAACAATTACGCGCGAGGCCGTTAC
TCTGTGGGTCGGAGGTCATCGACCTTGTGCTGGAGAGGACCCGGAAGCTGGCAGAACAGTGTGGTGGG
CTTCAGGGATTTTTGATTTTCCGAAGCTTTGGAGGAGGCACTGGTTCAGGGTTTACGTCTCTCTTAATG
GAGAGGCTCACAGGAGAATATAGCAGAAAAGACTAAGCTGGAGTTCTCGGTCTACCCAGCCCCAGGATC
TCCACTGCTGTGGTAGAGCCTTATAACTCTGCTCCTCACCACCCACTCCACCACAGAGCACACGGACTGT
ACCTTCATGGTGGACAACGAGGCCGTCTATGATATATGCCATCGTAAACTCGGTGTTGAATGCCCTCT
CATGCCAGCATCAATAGATTGGTGGTTCAGGTGGTATCTTCCATCACTGCCTCCCTCCGTTTGAAGGG
CCCTTGAATGTAGACCTAATTGAATCCAGACCAACCTGGTACCTATCCGAGAATACATTTCCCATG
ACAGCCTTCGCCCCATCGTCTCTGCTGACAAAGCCTACCATGAGCAGTCTCTGTGTGACACATCACC
ACTGCCTGCTTTGAGTCTCCAACCAGCTGGTCAAGTGTGATCCTCGGCTTGGGAAGTACATGGCCTGC
TGCTACTCTATAGAGGGGATGTGGTCCCAAGGAAGTGAATGCAGCAATCGCAGCCACGAAGTCGAGG
CACTCTGTTCAAGTTTGTAGATTGGTGTCCAACCTGGTTTCAAGGTGGGCATCAACAATCGGCCGCCACG
GTGATGCCGGGTGGGACCTGGCCAAAGTCCACCGGTCCATCTGCATGCTGAGCAACACCACGGCGATT
GTGGAGGCTGGGCCCGCTGGACCACAAGTTTGACCTCATGTACGCCAAGAGAGCATTTCTGCACTGG
TACCTCAGAGAAGGATGGAAGAAGCAGAGTCTTGGAGGCCAGGAAGATCTGGCAGCCCTGGAGAGG
GACTATGAGGAAGTGGCGCAAAGTTTTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_024803

Insert Size: 1341 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_024803.2](#)

RefSeq Size: 1815 bp

RefSeq ORF: 1341 bp

Locus ID: 79861

UniProt ID: [A6NHL2](#)

Cytogenetics: 10p15.1

Protein Families:	Druggable Genome
Protein Pathways:	Gap junction, Pathogenic Escherichia coli infection
MW:	49.9 kDa
Gene Summary:	<p>Tubulin is the major constituent of microtubules. It binds two moles of GTP, one at an exchangeable site on the beta chain and one at a non-exchangeable site on the alpha chain (By similarity).[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>