

Product datasheet for **SC305175**

DOK3 (NM_024872) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DOK3 (NM_024872) Human Untagged Clone
Tag:	Tag Free
Symbol:	DOK3
Synonyms:	DOKL
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC305175 representing NM_024872.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGACTCGGGGAGCCAGACTGCGATCAGACGCGCGTCCCAGCTGAACCAGCTGTCTCTAGACGGAGGG
ACGGGAAGTGGCCAGAAGGGGAAGTGTGAGGAGTCCCGTCCAGCCTGTATCAGTCTCCCAGGTCTT
GAAGCGGCGGCCCTGCTCCTGGCCGTGACCATGGACCCTCTGGAGACCCCTATCAAGGATGGCATCCTC
TACCAGCAGCATGTCAAGTTTGGCAAGAAGTGCTGGCGGAAGGTGTGGGCTCTGTGTATGCAGGAGGC
CCATCAGGCGTGGCAGGCTGGAGAGCTGGGAGGTCCGGGATGGTGGCCTGGGAGCAGCGGGTGACAGG
TCGGCAGGGCCTGGCCGGCAGGGGAGCGACGGGTATCCGCCTGGTGACTGTGTGTCCGTGCTGCCG
GCTGACGGCAGAGCTGCCCCGGGACACCGTGCCTTCTGCTCACCACCAGCAGCGAAGCCATCTA
CTGGCTGCTCAGCACCAGCAGGCTGGATGGGCCCATCTGCCAGCTGGCCTTCCCGGGACAGGGGAG
GCCTCCTCAGGATCCACAGATGCCAGTCTCCCAAGAGGGGCTGGTCCCATGGAGAAAACCTCCATC
TACTCCTCCTGGCAGGAAGTGGCGGAGTTCCCGTGGTGGTGCAGAGGACTGAGGCGCCACCCGCTGC
CAGCTGAAGGGGCGGCCCTGCTGGTGTGGGCCAGACGCCATCCAGCTGAGGGAGGCCAAGGGCACC
CAGGCCCTCTACAGCTGGCCCTACCACTTCTGCGCAAGTTCCGGTCCGACAAGGGCGTGTCTCCCTTT
GAGGCCGGCCGTGCTGCCACTCGGGTGAAGGGCTCTTTGCCTCAGCACCCCTGTGCCCTGACCTG
TGCAGGGCTGTGGCCGGGGCCATCGCCGCCAGCGGGAGCGGCTGCCAGAGCTGACCAGGCCCCAGCCC
TGCCCCCTGCCACGGGCCACCTCTGCCCTCCCTGGACACCCCGGAGAGCTTCGGGAGATGCCACCA
GGACCTGAGCCACCCACGTCCAGAAAAATGCACCTGGCCGAGCCCGACCCAGAGCCTGCCGCTACTG
CTAGGCCCGGAGCCCAACGATCTGGCGTCCGGGCTCTACGCTTCAAGTGTGCAAGCGTCCAGTGGGCC
CCAGGCAATGAGCACCTCTATGAGAACCCTGTGTGTGCTGGAGGCCAGCCACAGCTGCAGGTGGGGAA
CCTGAGCCGACAGGGGCCCGGACGCCAGCCCAACCAAGTCCCATCTACCAACAGGCCAGGAC
TTGAGCTGGCCCGCCGGCCAAACGACAGTACCCTGGAGGCCAGTACCGCGGCTGCTGGAGCTGGAT
CAGGTGGAGGGCACAGGCCGCCCTGACCCTCAGGCAGGTTTCAAGGCCAAGCTGGTACCCTGCTGAGT
CGTGAAGCGGAGGAAGGGCCAGCCCTTGTGACCGGCCTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_024872

Insert Size:	1491 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_024872.2
RefSeq Size:	1747 bp
RefSeq ORF:	1491 bp
Locus ID:	79930
UniProt ID:	Q7L591
Cytogenetics:	5q35.3
Protein Families:	Druggable Genome
MW:	53.3 kDa
Gene Summary:	DOK proteins are enzymatically inert adaptor or scaffolding proteins. They provide a docking platform for the assembly of multimolecular signaling complexes. DOK3 is a negative regulator of JNK signaling in B-cells through interaction with INPP5D/SHIP1. May modulate ABL1 function (By similarity).[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) encodes the longest isoform (1).