

Product datasheet for **SC305183**

LPPR3 (PLPPR3) (NM_024888) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LPPR3 (PLPPR3) (NM_024888) Human Untagged Clone
Tag:	Tag Free
Symbol:	PLPPR3
Synonyms:	LPPR3; LPR3; PRG-2; PRG2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC305183 representing NM_024888.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGATCTCCACCAAGGAGAAGAACAAGATCCCCAAGGACAGCATGACGCTTCTGCCCTGCTTCTACTTC
GTGGAGCTGCCATAGTGGCTTCTCCATCGTATCCTTGTACTTCTGGAGCTGACCGACCTCTTCAAG
CCGGCCAAGGTGGGCTTCCAGTGTATGACCGCACTCTCTCCATGCCCTACGTGGAGACCAACGAGGAG
CTCATCCCCTGCTGATGCTGCTCAGCTTGGCCTTCGCGGCCCTGCCGCCTCGATCATGGTGGCCGAG
GGCATGTTGACTGTCTGCAGTCCCGGTGTGGGGCCGTGCCGGGGGCCCGCGGGGGCGGAGGGCAGC
ATCAACGCCGGCGGTGCAACTTCAACTCTTCTGCGGCGTACGGTGGGTTTGTGGTGTCCACGTG
TTCGCCCTGTGTCCACAGCCCTGGTACGACGCTGATCCAGCTGGCCACGGTTACCACACTCCCTTC
TTCCTCACCGTCTGCAAGCCCACTACACTCTCTGGGCAGCTCTGCGAGGTCAACCCCTACATCAG
CAGGACATCTGCTCCGGCCAGACATCCACGCCATCTGTCTGCACGGAAGACCTTCCCGTCCCAGCAC
GCCACGCTGTAGCCTTCGCCGGTCTATGTGTGGTGTGAGTCCGGCCCTCACTGCCCTTCCCAGGCC
CTTTGCTGACCGTGGGAGGCCCTCCCTGACCCCAACCCCATGCCCCAGATGACTTCAACTCGGTC
ATCTCGGACACCAAGCTGCTGAAGCCATCCTGGTCTTCGCCTTTCGCATCGCCGGGGCGTATGC
GGGCTCACGCAGATCACGCAGTACCGACGCCACCTGTGGACGTGATGCCGGCTTCTCATCGGGCG
GGCATCGTGCCTACCTGGCCTGCCACGCGGTGGGCAACTTCCAGGCCACCTGCAGAGAAGCCCGCG
GCCCGGCCCGCCCAAGGACGCGCTGCGGGCCCTGACGACGCGGGGCCACGACTCGGTTTATCAGCAG
AATAAGTCGGTGAACACGAGCTGGGGCCCCAGGGCGGCTGGAGGGCGGCCCGGCCCGTGGCC
CGCGAGAAGACCTCGCTGGCAGCCTGAAGCGCGCCAGCGTGGACGTGGACCTGCTGGCCCCGCGCAG
CCCATGGCCAAGGAGAACATGGTGACCTTCAGCCACAGCTGCCAGGGCCAGCGCCCTCGCTGGAC
GACCCCGCGCGCCACATGACCATCCAGTGGCGCTGGACGCCCTCGCGCTCCAAGCAGCTCATCAGC
GAGTGGAAAGCAGAAGAGCCTGGAGGGCCGCGCCTGGGGCTGCCGACGACGCCAGCCCCGGCACCTG
CGCGCGCCCGCCAAACCATGGCGGAGGAGGAGGAAGAGGAGGAGGACGAAGAGGAAGAGGAGGAG
GAAGAGGAGGAGGACGAGGGCCCGCCCGCCCTCGCTTACCCACCGTGCAGGCGCGGCCGGGCTG
GGGCTCGGGTATCTCCACCGCGCGGGGCCCGCGCGCTGGTGCACATCCCGGAGGAGGGCGCG
CAGACGGGGCCGGCCTGTCCCCAAAAGCGCGCGGGGTGCGCGCCAAGTGGCTCATGATGGCCGAG
AAGAGCGGGCGGCAGTGGCCAACCTCCGCGGCTGCTGCAGGTATCGCCATGTCCAAGGCTCCGGGC
GCGCCGGGCCCAAGGCGGCCGAGACGGGCTCGTCCAGCGCCAGCTCCGACTCCTCGCAGTACCGG
TCGCCGTGCGACCGGACTCCGCCAGCATCGTGACCATGACGCGCACGCGCCGACCACCCCGTGGTG
CACCTGTGCGCCGGCGCGCCCTGGGAGTGGAAAGCGCGGGCGGGCGGGGCCAAGGCGGAGGCCGAC
GGCGGCTACGAGCTGGGGACCTGGCGCGCGGCTTCCGCGCGGGGCCAAGCCCCGGGCGTGTCCCC
GGCTCGTCGGTACGCGACGTGGACCAGGAGGAGCCGCGGTTTCGGGGCCGTGGCCACCGTCAACCTGGCC
ACGGGCGAGGGGCTGCCCCGCTGGGCGCGCCGATGGGCGCTGGGCCCGGGCAGCCGGGAGTCCACG
CTGCGGCGCCACGCGGGCGCCTGGGGCTGGCGGAGCGGAGGCGGAGGCGGAGGCGGAGGGCTACTTC
CGCAAGATGCAGGCGCGCCGCTTCCCCGACTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_024888

Insert Size: 2241 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_024888.2](#)

RefSeq Size: 2405 bp

RefSeq ORF: 2241 bp

Locus ID: 79948

UniProt ID: [Q6T4P5](#)

Cytogenetics: 19p13.3

Protein Families: Transmembrane

MW: 78.9 kDa

Gene Summary: The proteins in the lipid phosphate phosphatase (LPP) family, including PRG2, are integral membrane proteins that modulate bioactive lipid phosphates including phosphatidate, lysophosphatidate, and sphingosine-1-phosphate in the context of cell migration, neurite retraction, and mitogenesis (Brauer et al., 2003 [PubMed 12730698]).[supplied by OMIM, Mar 2008]

Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (1).