

## Product datasheet for **SC305203**

### GPR157 (NM\_024980) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	GPR157 (NM_024980) Human Untagged Clone
Tag:	Tag Free
Symbol:	GPR157
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>OriGene SC305203 ORF sequence for NM_024980, the custom clone sequence may differ by one or more nucleotides

```

ATGCAGCCGTCCTCCGCGCCACCGAGCTGGTGGCCGTCGGAGCGCGCCGTGGTGCTGCTGCTGCGTGCGCAC
TCTCCGCGCTCGGCTCGGGCCTGCTGGTGGCCACGCACGCCCTGTGGCCGACCTGCGCAGCCGGGCACG
GCGCCTGCTGCTCTTCTGTGCTGGCCGACCTGCTCTCGGCCGCTCTACTTCTACGGAGTGCTGCAG
AACTTCGCGGGCCCGTCGTGGGACTGCGTGCTGCAGGGCGCGCTGTCCACCTTCGCCAACACCAGCTCCT
TCTTCTGGACCGTGGCCATTGCGCTCTACTTGTACCTCAGCATCGTCCGCGCCGCGCGGGCCCTCGCAC
AGATCGCCTGCTTTGGGCCTTCCATGTCGTGCTGACGCTGGGGGTCCCGTTGGTCATCACTGTGGCGCCGTC
GCCCTGAAGAAGATTGGCTATGACGCTCGGACGTGCTGTGGGCTGGTGCTGGATCGACTGGAGGCCA
AGGACCATGTCTGTGGATGCTGCTGACGGGAAGCTGTGGGAGATGCTGGCATATGTGCTGCTGCCTCT
GCTGTACCTCCTGGTCCGGAAGCACATCAACAGAGCGCACACGGCACTCTNTGAGTACCGGCCATCCTC
TCCCAGGAGCACCCTGCTGCGCCACTCCTCCATGGCGGACAAGAAGCTGGTGCTCATCCCGCTCATCT
TCATCGGCCTCAGGGTCTGGAGCACCCTGCGGTTCTGCTGACCCCTCTGTGGTCCCCGGCCGTGCAGAC
GCCGGTGTGGTGGTCTGCATGGTATCGGGAACACGTTTCAGGGAGGTGCCAACTGCATCATGTTCTGTC
CTCTGCACCCGCGCCGTCGAACTCGGCTCTTCTCTCTGTTGCTGCTGCTCTTCTCAGCCTCCCA
CCAAGAGCCCGGCTGGCACTCCCAAGGCTCCCGCGCTTCCAAGCCAGGAGAATCTCAGGAATCCCAAGG
GACCCAGGGAACTTCCAAGCACCTGA

```

Restriction Sites:	Sgfl-MluI
ACCN:	NM_024980
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_024980.3</a></u> , <u><a href="#">NP_079256.3</a></u>
<b>RefSeq Size:</b>	1153 bp
<b>RefSeq ORF:</b>	1008 bp
<b>Locus ID:</b>	80045
<b>UniProt ID:</b>	<u><a href="#">Q5UAW9</a></u>
<b>Cytogenetics:</b>	1p36.22
<b>Protein Families:</b>	Druggable Genome, Transmembrane
<b>Gene Summary:</b>	Orphan receptor that promotes neuronal differentiation of radial glial progenitors (RGPs). The activity of this receptor is mediated by a G(q)-protein that activates a phosphatidylinositol-calcium second messenger.[UniProtKB/Swiss-Prot Function]