

Product datasheet for **SC305238**

ARMC9 (NM_025139) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ARMC9 (NM_025139) Human Untagged Clone
Tag:	Tag Free
Symbol:	ARMC9
Synonyms:	ARM; JBTS30; KU-MEL-1; NS21
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC305238 representing NM_025139.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGGGGGACATTCTGGCTCATGAATCTGAATTACTTGGACTAGTGAAGAGTATTTAGATTTTGCTGAA
TTTGAAGACACCTTGAAAACATTTTCAAAAGAATGCAAAATAAAAGGAAAACCATTTGTGTAAACAGTA
GGCGGATCTTTTCAAGAGACTCCAAATCATTGACAATTCAGAAGGATCTTGTGCTGCAATTTGACAACGGA
GACCAGAAGGTGTTCTTCGATCTGTGGGAGGAGCACATTTCAAGTTCATCCGAGATGGGGACTCCTTT
GCCCAGAAGCTGGAATTCTATCTCCACATCCATTTTGCCATCTATCTTTTGAAGTACTCTGTGGGAGA
CCGGACAAAGAGGAGCTGGATGAAAAGATTTCTACTTCAAAACCTACCTGGAGACCAAGGGGACGCC
TTGAGCCAGACCAGAGTTTCTTCTTTCTATGCCCTTCTTTTGTTCACACCTATGGTGCACCCC
TCATTTAAAGAACTCTTCCAGGATTCCTGGACTCCAGAGTTAAAGTTGAAGTTGAAAAGTTTCTAGCT
TTAATATCTAAAGCCAGCAACACGCCAAAGCTTTTAAACAATATATAAGGAGAATGGACAAAGTAACAAA
GAAATCTTGCAGCAGCTCCACCAGCAGCTGGTTGAAGCTGAACGTAGGTCACTGACATACCTCAAACGG
TACAATAAGATCCAGGCCGACTACCACAATCTCATTGGAGTCACAGCAGAGCTGGTGGATTCTCTAGAG
GCCACAGTCAGCGGAAGATGATCACCCCTGAGTACCTCCAGAGCGTCTGTGTCCGCTGTTTCAGTAAC
CAGATGCGGCAGAGCCTGGCGCATAGTGTGGACTTCACGAGGCCTGGGACGGCATCCACCATGTTACGA
GCCTCCTTGGCACCCGTGAAATTGAAGGATGTCCCATTTACTGCCCTCCTTGGATTATGAGAACTGAAG
AAGGATTTGATTTGGGGAGTGACCGCTTGAAAGCCTTCTTGTTCAGGCTCTGCGCTGGCGCTTGACC
ACATCCCATCTTGAGAGCAGAGGGAGACCGTTCTGCAAGCCTACATCAGCAATGACCTCTTGGACTGT
TATAGCCACAACCAGAGGAGTGCTTTCAGTTGCTGCACTCCACGAGCGAGCTGGTGGGACGTACATG
GCCAGGCTCATCAATGCTTTTGGCTCAGTGGCAGAAGGTGCGCTCTACCTTGGCCAGAACACAAAGGTG
CTGCAGATGCTGGAGGGAAGGCTGAAGGAGGAGGACAAGGATATCATCACCAGGGAAGTGTCTTGGG
GCCCTGCAGAAGTTCAGTCTCAGGCGCCGCTGCAGACAGCGATGATTCAAGACGCGCTCATCTTCTGG
CTGTTGATGTTCTGAAGGACCCTGACTGCCTGTCTGACTACACGCTGGAGTACTCGGTGGCTTTGCTC
ATGAACCTCTGCCTCCGACGACAGGGAAGAACATGTGTGCCAAGGTGGCAGGCTCGTGCTCAAAGTC
CTTTGGATCTTCTTGGCCATGAAAACCATGAGATACAGCCGTATGTGAATGGAGCTCTGTACAGCATC
CTTTCTGTTCCATCCATTCTGTGAGGAAGCAAGAGCAATGGGAATGGAAGACATCTACGCTGCTTCATC
AAAGAAGGCAATGCTGAAATGATCCGCCAGATAGAATTATCATCAAGCAGCTAAATCCGAAGAGCTA
CCAGATGGTGTCTTGAATCTGATGATGATGAAGATGAAGATGATGAAGAGGACCATGACATCATGGAA
GCCGATCTGGACAAAGACGAACTGATCCAGCCCCAGCTCGGAGAACTCTCAGGAGAGAAGCTTCTGACC
ACGGAGTACCTGGGGATCATGACCAACACGGGGAAGACAAGGCGGAAGGGGCTGGCTAATGTGCAGTGG
AGCGGGGATGAGCCCCTGCAAAGGCCCGTCACCCCCGGCGGCCACAGAAACGGGTACCCAGTGTA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: Sgfl-MluI

ACCN: NM_025139

Insert Size: 1998 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_025139.5</u>
RefSeq Size:	2317 bp
RefSeq ORF:	1998 bp
Locus ID:	80210
UniProt ID:	<u>Q7Z3E5</u>
Cytogenetics:	2q37.1
MW:	75.6 kDa
Gene Summary:	Required for ciliogenesis.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (2) differs in the 5' UTR and lacks multiple 3' coding exons. Its 3' terminal exon extends past a splice site used in variant 1, resulting in a different 3' UTR. The encoded isoform (2) has a shorter C-terminus than isoform 1. Variants 2 and