

## Product datasheet for **SC305267**

### **CPEB1 (NM\_030594) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	CPEB1 (NM_030594) Human Untagged Clone
Tag:	Tag Free
Symbol:	CPEB1
Synonyms:	CPE-BP1; CPEB; CPEB-1; h-CPEB; hCPEB-1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

**Fully Sequenced ORF:** >SC305267 representing NM\_030594.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGCGTTCCCCGCTGGAAGAAGAAGCAGGAAGGATAAAAGATTGCTGGGACAACCAGGAAGCACCTGCT
CTCTCCACGTGTAGTAATGCCAATATCTTTCGAAGGATAAATGCCATATTGGATAATTCTCTGGATTC
AGTAGAGTCTGCACTACACCTATAAACCAGGAATTCATGATCATTTGCCAGACTTCCAGGACTCTGAA
GAAACAGTTACAAGCAGGATGCTTTTCCCAACCTCTGCGCAAGAATCTTCCCGTGGCCTCCAGATGCA
AATGACTTGTGCCTTGGCCTGCAGTCCCTCAGTCTGACAGGCTGGGACCGACCCTGGAGCACCCAGGAC
TCAGATTCTCAGCCAGAGCAGCACACTCGGTACTGAGCATGCTCCATAACCCACTGGGAAATGTC
CTAGGAAAACCCCTTGAGCTTCTGCCTCTGGATCCCTTGGGTCTGACTTGGTGGACAAGTTTCCA
GCACCCTCAGTTAGAGGATCACGCCTGGACACCCGGCCATCCTGGACTCTCGATCTAGCAGCCCTCT
GACTCAGACACCAGTGGCTTCTGCTCTGTGACGGGGTGGTCCAGAGACCCCTTAAAGATGGGGTAGGG
TCTCGGATGGACCAAGAGCAAGCTGCTTGTCTGCAAGTCACTCCCTCCCAACCAGTCTCAAAGAGA
TGGCCAGGAGCTTCTGTGTGGCCATCCTGGGACCTCCTCGAAGCTCCCAAGACCCCTTCAGCATAGAG
AGAGAGGCCAGGCTGCACCGACAAGCTGCAGCTGTGAATGAAGCCACCTGTACCTGGAGTGGCCAGCTT
CCTCCCGGAAGTATAAGAACCCATCTACTCTTGAAGGTGTTTCTAGGAGGTGTTCTTGGGATATT
ACAGAAGCTGGATTAGTTAACACCTTCCGTGTTTTGGCTCTTTGAGTGTGGAGTGGCCTGGTAAGGAT
GGCAAGCATCCCGGTGTCTCCAAAGGGTATGTGTATCTGGTCTTCCAAGTATAGAGAAGTCTGTCCGA
TCCTTGCTTCAAGCTTGTCTCATGACCCGCTGAGCCAGATGGCCTGAGTGAATATTATTTCAAGATG
TCCAGCCGAAGGATGCGCTGCAAGGAGGTGCAGGTGATCCCTGGGTATTAGCCGACAGTAACTTTGTC
CGGAGCCCATCTCAGAGGCTTGACCCAGCAGGACGGTGTGTCGGTGTCTGCATGGAATGCTAAAT
GCTGAGGCCCTGGCAGCCATCTTGAACGACCTATTTGGTGGAGTGGTGTATGCCGGGATTGACACAGAT
AAGCACAAGTATCCCATTTGGTCTGGTCTGTGACTTTCAATAACCAACGGAGTTACCTGAAAGCAGTC
AGCGCTGCTTTTGTGGAGATCAAACCACCAAGTTCACAAAGAAGGTTACAGATTGACCCCTACCTAGAA
GATTCTCTGTGCATATCTGCAGTCTCAGCCTGGTCTTTCTTCTGTGAGATCAGGTCTGCTTCAA
TACTTCTGCCGGAGCTGCTGGCACTGGCGGCACAGCATGGAGGGCTGCGCCACCACAGCCCTGATG
CGGAACCAGAAGAACCGAGATTCAGCTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**Plasmid Map:**


**ACCN:** NM\_030594

**Insert Size:** 1686 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_030594.4](#)

**RefSeq Size:** 3241 bp

**RefSeq ORF:** 1686 bp

**Locus ID:** 64506

**UniProt ID:** [Q9BZB8](#)

**Cytogenetics:** 15q25.2

<b>Protein Pathways:</b>	Dorso-ventral axis formation, Oocyte meiosis, Progesterone-mediated oocyte maturation
<b>MW:</b>	62.1 kDa
<b>Gene Summary:</b>	<p>This gene encodes a member of the cytoplasmic polyadenylation element binding protein family. This highly conserved protein binds to a specific RNA sequence, called the cytoplasmic polyadenylation element, found in the 3' untranslated region of some mRNAs. The encoded protein functions in both the cytoplasm and the nucleus. It is involved in the regulation of mRNA translation, as well as processing of the 3' untranslated region, and may play a role in cell proliferation and tumorigenesis. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2014]</p> <p>Transcript Variant: This variant (1) encodes the longest isoform (1).</p>