

Product datasheet for **SC305273**

LRRC27 (NM_030626) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LRRC27 (NM_030626) Human Untagged Clone
Tag:	Tag Free
Symbol:	LRRC27
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC305273 representing NM_030626.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGAGGGAAGCAGCTCCTACGAAGTTCCTCTGTGGCTGCTGATCTGGAGGAGGGTCTGGTTCAG
ACTAGGAGCTTGCCTGCCACCCCTCCAAAGATGTTACAAGGGTGTGGAGGCATCATCTTTCTCTCC
TCACCGATTTTAGACTTGAGTGAAAGTGGTCTGTGCCGTTTGGAGGAGGTCTTTAGAATCCCCAGCCTT
CAACAATTGCATCTGCAAGGAATGCCCTGTGTGTGATTCTCAAGATTTCTTTAGTTGCTTCCGAAC
CTGACTTGGCTGGACCTCCGGTACAATAGAATTAAGCGCTTCTTCTGGGATTGGAGCTCACCAGCAT
TTGAAACTTTGCTTTTGAAGAAATCTATCAAAATGTTACCTGTGGAGCTGGGGAGCGTAACCACG
CTGAAAGCACTGAACCTAAGACACTGCCCTCTGGAATTCCTCTCAGCTCGTTGTGAGAAGGGATTG
GTGGCTATCCAGCGCTTCTGCGGATGTGGCAGTAGAAGCTCTCTCCAGAAATCCAATTCTCAA
GAGGCTCCACCGTTAGAGAGATGACCTCCGTGACCTCCGAGCCAGGAGTGGAGTTGTCTGGAGAC
CACGCGTCTAACCAAGGAGCTGTGAACGCTCAGGACCCAGAGGGGGCTGTGATGAAAGAGAAGGCCAGC
TTTCTCCACCTGTGAAAAGCCAGACCTGAGTGAAGTCTGCGGACTCCTCAGAGAAGTGG
CCCAGCGAGGAGGAGATCAGGCGCTTTTGAAGCTGAGGCAGGAGATTGTTGAGCACGTGAAGGCAGAC
GTTCTGGGAGATCAGCTCTTGACGAGGGAATTACCTCCAAATCTCAAGCGGCCCTTGAACATTGAGAAA
GAAGTACCAAGCAAGACACGTTTTAGAGGAAGACAGCCTCTCCAGGAGCATCTTACCCGACCTC
TTGTACCGTACCAATGGCGATCCGAGCAAAAAGACTGGAAGAGAGCCGAGCGGCGGCTCCGAGAG
CTCCAGGAGAAGCAGGCTCTGATGGAGCAGCAGAGACGAGAGAAAAGGGCACTGCAGGAGTGGAGAGAG
CGAGCCCAGAGGATGAGGAAGAGGAAGGAGAGCTCAGCAAACTCTGCCTCCGCGAGGAGCATGGTG
GCATCAAAAGATTCCCTCTGCCACAGATCTGATAGATAACAGGAAAGTACCACTGAATCCGCTGAAAAA
ATGAAACCAAGCAAGAGAAATCGCCACAAGCAAGTAAAGAAATGAGTGCCTGCAGGAGAGAAATTTA
GAAGAGAAGATAAAACAGCACGTCCTCCAAATGCGTGAGCAAAAGAAATTCCATGGCCAGGCCCCACTG
GAGGAGATGAGGAAGGCTGCCGAGGATCTGAAATTGCCACAGAGCTACAGGATGAAGTATTGAAGCTA
AAATTGGGATTAACCTTGAACAAAGATCGTCGACGGGCGGCCCTCACTGGAAACCTTTGCTTGGCTG
CCGGCAGCACAGCCTCAAATACATTTTTTAACACAAAATATGGAGAATCAGGAAATGTTGCGAGATAC
CAGTGA
AGCGGACCGACGCTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGAT
ATCCTGGATTACAAGGATGACGACGATAAGGTTTAA
  
```

Restriction Sites: SgfI-RsrII

ACCN: NM_030626

Insert Size: 1593 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_030626.2](#)

RefSeq Size: 8052 bp

RefSeq ORF: 1593 bp

Locus ID: 80313

UniProt ID: [Q9C0I9](#)

Cytogenetics: 10q26.3

MW: 60.1 kDa