

Product datasheet for **SC305283**

GJA9 (NM_030772) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GJA9 (NM_030772) Human Untagged Clone
Tag:	Tag Free
Symbol:	GJA9
Synonyms:	CX58; CX59; GJA10
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC305283 representing NM_030772.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGGGGACTGGAATCTCCTTGGAGATACTCTGGAGGAAGTTCACATCCACTCCACCATGATTGGAAG
ATCTGGCTCACCATCCTGTTTCATATTTTGAATGCTTGTCTGGGTGTAGCAGCTGAAGATGTCTGGAAT
GATGAGCAGTCTGGCTTCATCTGCAATACAGAACAACAGGCTGCAGAAATGTATGCTACGACCAGGCC
TTTCTATCTCCCTCATTAGATACTGGGTCTGCAGGTGATATTTGTGTCTTCACCATCCCTGGTCTAC
ATGGGCCATGCATTGTACCGACTGAGAGTCTTGGAGGAGAGAGGCAAGGATGAAAGCTCAGTTAAGA
GTAGAAGTGGAGGAGGTAGAGTTTGAATGCCTAGGGATCGGAGGAGATTGGAGCAAGAGCTTTGTCAG
CTGGAGAAAAGGAACTAAATAAAGCTCCACTCAGAGGAACCTTGTCTTGCACCTATGTATACACATT
TTCACTCGCTCTGTGTTGAAGTTGGATTGATGTTGACAGTACCTTTTATATGATTTCACTTAGAG
CCGCTATTTAAGTGCCATGGCCACCCGTGTCAAATATAATCGACTGTTTGTCTCAAGACCAACAGAA
AAGACAATATTCTATTATTTATGCAATCTATAGCCACTATTTCACTTTTCTAAACATTCTTGAAATT
TTCCACCTAGGTTTTAAAAAGATTAAAAAGAGGGCTTTGGGAAAAATACAAGTTGAAGAAGGAACATAAT
GAATTCATGCAACAAGGCAAAACAAAATGTAGCCAAATACCAGAGCACATCTGCAATTCAGTGAAG
CGACTCCCTTCTGCCCTGATTATAATCTGTTAGTGAAAAGCAACACACACTGCAAGTGTACCCTAGT
TTAAATTCATCTTCTGTATTCCAGCCAAATCCTGACAATCATAGTGTAATGATGAGAAATGCATTTTG
GATGAACAGGAACTGTACTTTCTAATGAGATTCCCACTTAGTACTAGTTGTAGTCATTTTCAACAC
ATCAGTTCAAACAATAACAAAGACACTCATAAAATATTTGAAAAGAACTTAATGGTAACCAAGTTAATG
GAAAAAGAGAACTGAAGGCAAGACAGCAAAAGGAAGTACTACTCTAGAGGTACCCGTTCTATTCCA
GGTGTGCTATAGATGGAGAGAAACATGAGGCAGTCACCCCAACAGTTTCTCCTTGCCAGCTAAC
TGCATTGGAAACCGCGGTGGCTTAGAGCTACATGGGGTTCCTCTACAGAACATGAAAACCGGGGGTCA
CCTCTAAAGGTAACTCAAGGGCCAGTTCAGAAAGGGCACAGTCAGAACCCCTCCTCCTTCACAAGGA
GATTCTCAATCACTTGACATTCCAAACTGCTGATTCTTTGGGAGGGCTGTCTTTGAGCCAGGTTTG
GTCAGAACCTGTAATAATCCTGTTGCTCCTCAAATCACGTAGTGTCCCTAACGAACAATCTCATTGGT
AGGCGGGTCCCACAGATCTTCAGATCTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_030772

Insert Size: 1548 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_030772.4</u>
RefSeq Size:	2340 bp
RefSeq ORF:	1548 bp
Locus ID:	81025
UniProt ID:	<u>P57773</u>
Cytogenetics:	1p34.3
Protein Families:	Transmembrane
MW:	58.8 kDa
Gene Summary:	Connexins, such as GJA9, are involved in the formation of gap junctions, intercellular conduits that directly connect the cytoplasm of contacting cells. Each gap junction channel is formed by docking of 2 hemichannels, each of which contains 6 connexin subunits (Sohl et al., 2003 [PubMed 12881038]).[supplied by OMIM, Mar 2008]