

Product datasheet for SC305292

LRRC3 (NM_030891) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LRRC3 (NM_030891) Human Untagged Clone
Tag:	Tag Free
Symbol:	LRRC3
Synonyms:	C21orf30; C21orf102; LRRC3DN
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC305292 representing NM_030891. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
 ATGGGCACCGTGCGCCACCTCGCCCTCGCTCCTGCTGGTCTCCACCCGGGAGTCTTGTCTCTTCCTC
 CTCTTCTGCCTGCACCTGGGCGCCGCTGCCACAGCCCTGCCGGTGCCCTGACCACGCAGGGGCTGTG
 GCTGTCTTCTGCAGCTTGCGGGGCCTTCAGGAGGTCCCGAGGACATCCCGGCCAACACCGTGCTCCTG
 AAGCTCGATGCCAACAGATCTCCACCTCCCGACGGGGCCTTCAGCACCTGCACCGGCTCAGGGAG
 CTGGATCTGTCTACAACGCCATCGAGGCCATCGGCTCCGCCACCTTCGCGGGCCTGGCCGGGGCCTG
 CGGCTGCTGGACCTGTCTTACAACGCATCCAGAGGATCCCCAAGGACGCCCTGGCAAACCTCAGCGCC
 AAGATACGCCTGTCCACAACCCCTGCACTGCGAGTGCGCCCTGCAGGAGGCCCTGTGGGAGCTGAAG
 CTGGACCCCGACTCTGTGGACGAGATCGCTGCCACACCTCAGTGCAGGAGGAGTTTGTGGGAAGCCT
 CTGGTTCAGGCTCTGGATGCGGGTGCCAGCCTCTGCAGCGTCCCCACAGGACCACAGACGTGGCCATG
 CTGGTCACCATGTTTCGGCTGGTTCGCCATGGTGATCGCCTACGTCGTGTACTATGTGCGCCACAACAG
 GAGGATGCCCGGAGGCACCTGGAGTACCTGAAGTCTTGCCAGCGCCCCGCTCCAAGGACCCCATC
 GGCCCGGGGCCCTAG
 ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Restriction Sites:	SgfI-MluI
ACCN:	NM_030891
Insert Size:	774 bp


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_030891.4
RefSeq Size:	2544 bp
RefSeq ORF:	774 bp
Locus ID:	81543
UniProt ID:	Q9BY71
Cytogenetics:	21q22.3
Protein Families:	Transmembrane
MW:	28.1 kDa