

Product datasheet for **SC305397**

PCDHGB3 (NM_032097) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PCDHGB3 (NM_032097) Human Untagged Clone
Tag:	Tag Free
Symbol:	PCDHGB3
Synonyms:	PCDH-GAMMA-B3
Vector:	<u>pCMV6 series</u>



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_032097, the custom clone sequence may differ by one or more nucleotides

```

ATGGGAAATAGCTCCGGATGGAGGGGCCAGCAGGGCAGAGGCGAATGCTATTTCTCTTC
CTGCTCTCTTTGTAGACCAGGTTCTCTCCGAACCGATCCGCTACGCTATTCCCAGGAG
CTGGACAGGGGCTCGCTGGTAGGAACTCGCCAAGGACCTGGGGTTTGGCGTGGGGGAT
TTACCTACTAGGAACCTGCGGGTTATTGCAGAGAAGAAATCTTTACCGTGAGCCCCGAA
AATGGGAACTTACTTGTGAGCGACCGTATAGACCGAGAGGAGATTTGTGGCAAGAAGTCG
ACGTGTGTTCTGGAATTTGAAATGGTTGCTGAAAAGCCTTTAAACTTTTTTCATGTAAC
GTGCTGATCCAGGATATTAACGACAACCCACCGACCTTTAGCCAAAATATCACTGAGCTG
GAAATCAGCGAACTGGCTCTCACTGGAGCCACATTTGCCCTGGAATCTGCGCAAGATCCT
GATGTAGGTGTC AATTCGCTGCAGCAGTACTACCTCAGCCCTGATCCGCACTTCTCTTTG
ATTGAGAAAGGAGAACCTGGATGGCAGTAGGTACCCAGAGCTAGTACTGAAAGCACCCCTG
GACAGGGAAGAGCAGCCACATCACCACCTGGTCTCACAGCTGTGGATGGGGGCGAGCCC
TCCAGAAGCTGTACCACCAGATCAGGGTAATTGTGCGAGATGCAAATGATAACCCCCCA
GTATTTACTCAGGACATGTACAGGGTCAATGTTGCAGAGAACCTGCCCGTGGCTCCTCC
GTATTTAAAGTGATGGCCATTGACATGGATGAGGGCATCAATGCCGAAATCATCTATGCC
TTCATCAATATTGGCAAGGAAGTGAGACAACCTGTTCAAGCTGGACAGTAAACCGGGGAA
CTCACCCTATTGGGAACTGGACTTTGAAGAGAGAGATAGCTACACAATTGGGGTGAA
GCAAAGGATGGTGGACATCACACTGCATATTGTAAGTACAGATAGATATTTGAGATGAA
AATGACAATGCCCGGAGATAACCTGGCTTCTGAATCCCAACATATAACAAGAAGATGCT
GAGCTGGGGACTGCCGTTGCCCTGATCAAAACACATGATCTAGATTCTGGATTTAATGGA
GAAATCCTATGCCAACTAAAAGGAACTCCCTTTAAATCGTTCAAGATACAAAAAC
ACATACAGTTGGTGACAGATGGAGCCCTGGACCGGGAGCAGATCCCAGAATACAATGTG
ACGATCACAGCTACCGACAAAGGCAATCCACCGCTCTCCTCCAGCAAGACCATCACTCTG
CACATCCTTGATGTAACGACAACGTTCCCGTTTTCCACCAGGCCTCCTACACCGTGCAT
GTAGCTGAGAACAATCCGCTGGAGCCTCCATTGCGCATGTGAGAGCCTCGGATCCCGAC
TTGGGACCTAATGGCCTTGCTCCTACTACATCGTGGCCAGTGACCTGGAGCCGCGGGAG
CTGTGCTCTACGTGCTCGTGAGCGCGGAGCGGGTGGTGTTCGCGCAGCGAGCCTTC
GACCACGAGCAGCTGCGTGCCTTCGAGCTCACTCTGCAGGCCGCGACCAGGGCTCGCT
ACGCTCAGCGCCAACGTGAGCCTGCGCGTGTGGTGGACGACCGCAACGACAATGCACCG
CTGGTGTGTACCCAGCTCTGGGGCCGAAGGCTCTGCGCTCTTCGATATGGTGCCGCGC
TCTGCAGAGCCTGGCTACCTGGTGACCAAGGTGGTGGCGGTGGACGCGAGACTCGGGATAC
AACGCCTGGCTGTCTACCACATTTGTGCAGGCCAGCGAGCCCGGGCTGTTGAGCCTGGGC
CTGCGCACGGGTGAGGTGCGCACGGCGCTACCTTGGGCGACAGGGAGGCCGCCGCCAG
CGCCTGTGGTCACTGTGCGTGTGGAGGACAGCAGCCTCTTTCAGCCACCGTCAATGCTG
CACCTAATCTTCGAGATAGCTTGAAGAGATACAACCTGACCTTAGCGACCGCCCCACT
CCCTCTGACCCTCAGGCGGAGCTACAGTTTACCTAGTAGTGGCGTTGGCCTTGATCTCA
GTGCTCTTCTCCTCGCGGTGATTCTGGCAATCTCCCTGCGCCTGCGATGCTCCTCCAGA
CCCCTACTGAGGGTACTTTGAGCCTGGTGTCTGCTTCAAGACTGTACCTGGAGTTCTC
CCACCTACAGCGAAAGGACTTTGCCTTATTCTACAATCCGTGTGCTGCCTCACATTCC
TCAAACACCGAGTTTAAATTTCTCAATATAAAGGCTGAAAATGCTGCACCACAAGATCTT
CTATGTGATGAAGCCTCTGGTTTGAAGTAATGACAATCCAGAAATGCCTTCTAATTCA
GGCAATTTGAAAAGGTGAGTTTCTTCAAACCTTCTCCATAA
    
```

Restriction Sites: Please inquire

ACCN: NM_032097

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_032097.1 , NP_115268.1
RefSeq Size:	2445 bp
RefSeq ORF:	2445 bp
Locus ID:	56102
UniProt ID:	Q9Y5G1
Cytogenetics:	5q31.3
Protein Families:	Transmembrane
Gene Summary:	<p>This gene is a member of the protocadherin gamma gene cluster, one of three related clusters tandemly linked on chromosome five. These gene clusters have an immunoglobulin-like organization, suggesting that a novel mechanism may be involved in their regulation and expression. The gamma gene cluster includes 22 genes divided into 3 subfamilies. Subfamily A contains 12 genes, subfamily B contains 7 genes and 2 pseudogenes, and the more distantly related subfamily C contains 3 genes. The tandem array of 22 large, variable region exons are followed by a constant region, containing 3 exons shared by all genes in the cluster. Each variable region exon encodes the extracellular region, which includes 6 cadherin ectodomains and a transmembrane region. The constant region exons encode the common cytoplasmic region. These neural cadherin-like cell adhesion proteins most likely play a critical role in the establishment and function of specific cell-cell connections in the brain. Alternative splicing has been described for the gamma cluster genes. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (2) lacks the constant region exons and contains an alternate 3' coding region compared to variant 1. The encoded isoform (2) has a shorter, distinct C-terminus than isoform 1.</p>