

Product datasheet for SC305510

FIBCD1 (NM_032843) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FIBCD1 (NM_032843) Human Untagged Clone
Tag:	Tag Free
Symbol:	FIBCD1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC305510 representing NM_032843. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGTCAACGACCGGTGGAAGACCATGGGGCGGCTGCCCAACTGAGGACCGGCCGCGGACAAGCCG
CAGCGGCCGAGCTGCGGCTACGTGCTGTGCACCGTGTGCTGGCCCTGGCTGTGCTGGCTGTAGCT
GTCACCGGTGCCGTGCTCTTCTGAACACGCCACGCGCCGGGACGGCGCCCCACCTGTCGTGAGC
ACTGGGCTGCCAGCGCAACAGCGCCCTGGTCACTGTGAAAGGGCGGACAGCTCGCACCTCAGCATC
CTCATTGACCCGCGTGGCCCGACCTCACCAGACAGCTTCGCACGCCTGGAGAGCGCCAGGCCTCGGTG
CTGCAGGCGCTGACAGAGCACCCAGGCCACGCGCTGGTGGGGACCCAGGAGCAGGAGCTGCTGGAC
ACGCTGGCCGACCAAGCTGCCCGGCTGCTGGCCGAGCCTCAGAGCTGCAGACGGAGTGCATGGGGCTG
CGAAGGGGATGGCACGCTGGGCCAGGCCTCAGCGCCCTGCAGAGTGAGCAGGGCCGCTCATCCAG
CTTCTCTGAGAGCCAGGGCCACATGGCTCACCTGGTGAACCTCCGTGAGCGACATCCTGGATGCCCTG
CAGAGGGACCGGGGCTGGGCCGCCCCGAACAAGGCCGACCTCAGAGAGCGCCTGCCCGGGAAACC
CGGCCCGGGGCTGTGCCACTGGCTCCCGGCCCGAGACTGTCTGGACGTCTCCTAAGCGGACAGCAG
GACGATGGCGTACTCTGTCTTCCACCCACTACCCGGCCGGCTCCAGGTGACTGTGACATGCCG
ACGGACGGCGGCGGCTGGACGGTGTTCAGCGCCGGGAGGACGGCTCCGTGAACCTCTCCGGGGCTGG
GATGCGTACCGAGACGGCTTTGGCAGGCTCACCGGGAGCACTGGCTAGGGCTCAAGAGGATCCACGCC
CTGACCACACAGGCTGCCTACGAGCTGCAGTGGACCTGGAGGACTTTGAGAATGGCACGGCCTATGCC
CGCTACGGGAGCTTCGGCGTGGGCTTGTCTCCGTGGACCCTGAGGAAGACGGGTACCCGCTCACCGTG
GCTGACTATCCGGCACTGCAGGCGACTCCCTCCTGAAGCACAGCGGCATGAGGTTACCACCAAGGAC
CGTGACAGCGACCATTCAGAGAACAACGTGCCGCTTCTACCGCGGTGCTGGTGTACCGCAACTGC
CACAGTCCAACCTCAATGGGCAGTACCTGCGCGGTGCGCACGCCTCCTATGCCGACGGCGTGGAGTGG
TCCTCTGGACCGGCTGGCAGTACTCACTCAAGTTCTCTGAGATGAAGATCCGGCCGGTCCGGGAGGAC
CGCTAG
ACGGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGGC
```



[View online »](#)

Restriction Sites: Sgfl-MluI

Plasmid Map:



ACCN: NM_032843

Insert Size: 1386 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_032843.4](#)

RefSeq Size: 3257 bp

RefSeq ORF: 1386 bp

Locus ID: 84929

UniProt ID: [Q8N539](#)

Cytogenetics:	9q34.12
Protein Families:	Druggable Genome, Transmembrane
MW:	50.7 kDa
Gene Summary:	<p>FIBCD1 is a conserved type II transmembrane endocytic receptor that binds chitin and is located primarily in the intestinal brush border (Schlosser et al., 2009 [PubMed 19710473]). [supplied by OMIM, Apr 2010]</p> <p>Transcript Variant: This variant (1) represents the longer transcript.</p>