

## Product datasheet for **SC305517**

### FAM83A (NM\_032899) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	FAM83A (NM_032899) Human Untagged Clone
Tag:	Tag Free
Symbol:	FAM83A
Synonyms:	BJ-TSA-9
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC305517 representing NM_032899. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTGTAGTGAACCGTCAGAATTTGTAAATACGACTACTATAGGGCGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGAGCCGGTCAAGGCACCTGGGCAAAATCCGGAAGCGTCTGGAAGATGTCAAGAGCCAGTGGGTCGGG
CCAGCCAGGGCTGACTTTAGTGACAACGAGAGTGCCTGGCTGGCCACGGACGCCCTCTTGGATGGGGGT
TCTGAAGCCTACTGGCGGGTCTCAGCCAGGAAGGCGAGGTGGACTTCTTGTCTCGGTGGAGGCCAG
TACATCCAGGCCAGGCCAGGGAGCCCCGTGTCCCCAGACACCCTGGGAGGGGCGGAAGCAGGCCCT
AAGGGACTGGACTCCAGCTCCCTACAGTCCGGCACCTACTTCCCTGTGGCCTCAGAGGGCAGCGAGCCG
GCCCTACTGCACAGCTGGCCTCAGCTGAGAAGCCCTACCTGAAGGAAAAATCCAGCGCCACTGTGTAC
TTCCAGACCGTCAAGCACAACAACATCAGAGACCTCGTCCGCCGCTGCATCACCCGGACTAGCCAGGTC
CTGGTCATCTGATGGATGTGTTACGGATGTGGAGATCTTCTGTGACATTCTAGAGGCAGCCAACAAG
CGTGGGGTGTTCGTTTGTGTGCTCCTGGACCAGGGAGGTGTGAAGCTCTTCCAGGAGATGTGTGACAAA
GTCCAGATCTCTGACAGTCACTCAAGAACATTTCCATCCGGAGTGTGGAAGGAGAGATACTGTGCC
AAGTCAGGCAGGAAATTCGCTGGCCAAATCCGGGAGAAGTTCATCATCTCGGACTGGAGATTTGTCTGT
TCTGGATCTTACAGCTTACCTGGCTCTGCGGACACGTGCACCGGAACATCCTCTCCAAGTTCACAGGC
CAGGCGGTGGAGCTGTTTGACGAGGAGTCCGCCACCTCTACGCCTCCTCCAAGCCTGTGATGGGCTGTG
AAGTCCCCGCGGCTGGTCCGCCCGTCCCGCCGAGCAGCCCGGCCAATGGCCGCTTAGCAGCAGC
AGTGGCTCCGCCAGTGACCGCAGCTCCTCCAACCCCTTACGCGGCCGCTCGGCAGGCAGCCACCCCGGT
ACCCGAAGTGTGTCGCGTCTTACAGGGCCTGTAGCCCCGCGGCCACACCCGCTCCACCGCCCCGG
TTCCAGCCCCACCAAGGCCCTTGGGGAGCCCCAGTCCCCAGGCCACCTCTCCCCGCGGCCACAGC
GGCCCCGCGCGCTGTCTACAGCAACCTGGGGCCTACAGGCCACGCGGCTGCAGCTGGAGCAGCTG
GGCCTGGTCCGAGGCTGACTCCAACCTGGAGGCCCTTCTGCAGGCCTCCCTCACTTCTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online »](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM\_032899

Insert Size: 1305 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_032899.6](#)

RefSeq Size: 4639 bp

RefSeq ORF: 1305 bp

Locus ID: 84985

**UniProt ID:** [Q86UY5](#)

**Cytogenetics:** 8q24.13

**MW:** 47.5 kDa

**Gene Summary:** Probable proto-oncogene that functions in the epidermal growth factor receptor/EGFR signaling pathway. Activates both RAS/MAPK and PI3K/AKT/TOR signaling cascades downstream of EGFR. Required for the RAS/MAPK signaling cascade activation upon EGFR stimulation, it also activates both signaling cascades independently of EGFR activation. [UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (a). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.