

Product datasheet for **SC305520**

MT4 (NM_032935) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MT4 (NM_032935) Human Untagged Clone
Tag:	Tag Free
Symbol:	MT4
Synonyms:	MT-4; MT-IV; MTIV
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC305520 representing NM_032935. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCC**CGATCGCC**
 ATGGACCCAGGGAATGTGTCTGCATGTCTGGAGGAATCTGCATGTGTGGAGACAAGTCAAATGCACA
 ACCTGCAACTGTAAACATGTCGGAAGAGCTGCTGTCCCTGCTGCCCCCGGGCTGTGCCAAATGTGCC
 CGGGGCTGCATGCAAAGGAGGCTCAGACAAGTGCAGCTGCTGCCCA**TGA**
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Restriction Sites:	SgfI-MluI
ACCN:	NM_032935
Insert Size:	189 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).


[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_032935.2](#)

RefSeq Size: 294 bp

RefSeq ORF: 189 bp

Locus ID: 84560

UniProt ID: [P47944](#)

Cytogenetics: 16q13

MW: 6.4 kDa

Gene Summary: Seems to bind zinc and copper. Could play a special role in regulating zinc metabolism during the differentiation of stratified epithelia.[UniProtKB/Swiss-Prot Function]