

Product datasheet for **SC305563**

MYLK2 (NM_033118) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MYLK2 (NM_033118) Human Untagged Clone
Tag:	Tag Free
Symbol:	MYLK2
Synonyms:	KMLC; MLCK; MLCK2; skMLCK
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC305563 representing NM_033118.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGCGACAGAAAATGGAGCAGTTGAGCTGGGAATTCAGAACCATCAACAGACAAGGCACCTAAAGGT
CCCACAGGTGAAAGACCCTGGCTGCAGGAAAGACCCTGGCCCCCAGACCCAAAGAAAGCTCCGGAT
CCACCCACCCTGAAGAAAGATGCCAAAGCCCTGCCTCAGAGAAAGGGGATGGTACCCTGGCCCAACCC
TCAACTAGCAGCCAAGGCCCAAAGGAGAGGGTGACAGGGCGGGGGGCCCGCGGAGGGCAGTGCTGGG
CCCCCGCAGCCCTGCCCCAGCAGACTGCGACACCTGAGACCAGCGTCAAGAAGCCCAAGGCTGAGCAG
GGAGCCTCAGGCAGCCAGGATCTGAAAGCCAGGGTGGCAAGAAGGCAGCAGAGGGCCAAGCAGCA
GCCAGGAGGGGCTCACCTGCCTTTCTGCATAGCCCAGCTGTCTGCCATCATCTCCAGTTCTGAGAAG
CTGCTGGCCAAGAAGCCCCAAGCGAGGCATCAGAGCTCACCTTTGAAGGGGTGCCCATGACCCACAGC
CCCACGGATCCCAGGCCAGCCAAGGCAGAAGAAGAAAGAACATCTGGCAGAGAGCCAGAAGGAAGTG
GGAGAGAAAACCCAGGCCAGGCTGGCCAGGCTAAGATGCAAGGGGACACCTCGAGGGGGATTGAGTTC
CAGGCTGTTCCCTCAGAGAAATCCGAGGTGGGGCAGGCCCTCTGTCTCACAGCCAGGGAGGAGACTGC
TTCCAGATTTTGGATGATTGCCCGCCACCTCCGGCCCCCTTCCCTCACCCCATGGTGGAGCTGAGGACC
GGGAATGTGAGCAGTGAATTCAGTATGAACTCCAAGGAGGCGCTCGGAGGTGGCAAGTTTGGGGCAGTC
TGTACCTGCATGGAGAAAGCCACAGGCCTCAAGCTGGCAGCCAAGGTCATCAAGAAACAGACTCCCAAA
GACAAGGAAATGGTGTGCTGGAGATTGAGGTCATGAACCAGCTGAACCACCGCAATCTGATCCAGCTG
TATGCAGCCATCGAGACTCCGCATGAGATCGTCTGTTTTCATGGAGTACATCGAGGGCGGAGAGCTTTC
GAGAGGATTGTGGATGAGGACTACCATCTGACCGAGGTGGACACCATGGTGTGTTGTCAGGCAGATCTGT
GACGGGATCCTCTTCATGCACAAGATGAGGGTTTTGCACCTGGACCTCAAGCCAGAGAACATCCTGTGT
GTCAACACCACCGGCATTTGGTGAAGATCATTGACTTTGGCCTGGCACGGAGGTATAACCCCAACGAG
AAGCTGAAGTGAATTTGGGACCCAGAGTTCCTGTACCTGAGGTGGTGAATTATGACAAATCTCC
GATAAGACAGACATGTGGAGTATGGGGGTGATCACCTACATGCTGCTGAGCGGCCCTCCTCCCTTCTG
GGAGATGATGACACAGAGACCCTAAACAACGTTCTATCTGGCAACTGGTACTTTGATGAAGAGACCTTT
GAGGCCGTATCAGACGAGGCCAAAGACTTTGTCTCCAACCTCATCGTCAAGGACCAGAGGGCCCGGATG
AACGCTGCCAGTGTCTCGCCATCCCTGGCTCAACAACCTGGCGGAGAAAGCCAAACGCTGTAACCGA
CGCCTTAAGTCCCAGATCTTGTCTAAGAAATACCTCATGAAGAGGCGCTGGAAGAAAAATTCATTGCT
GTCAGCGCTGCCAACCGTTCAAGAAGATCAGCAGCTCGGGGCACTGATGGCTCTGGGGGTCGTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_033118

Insert Size: 1791 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_033118.3](#)

RefSeq Size: 2793 bp

RefSeq ORF: 1791 bp

Locus ID: 85366

UniProt ID: [Q9H1R3](#)

Cytogenetics: 20q11.21

Protein Families:	Druggable Genome, Protein Kinase
Protein Pathways:	Calcium signaling pathway, Focal adhesion, Regulation of actin cytoskeleton, Vascular smooth muscle contraction
MW:	64.7 kDa
Gene Summary:	This gene encodes a myosin light chain kinase, a calcium/calmodulin dependent enzyme, that is exclusively expressed in adult skeletal muscle. [provided by RefSeq, Jul 2008]