

Product datasheet for **SC305568**

Scratch 2 (SCRT2) (NM_033129) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Scratch 2 (SCRT2) (NM_033129) Human Untagged Clone
Tag:	Tag Free
Symbol:	Scratch 2
Synonyms:	ZNF898B
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC305568 representing NM_033129. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCC**CGCATCGCC**
 ATGCCGCGCTCCTTCCTGGTAAAGAAGATCAAAGGGGACGGCTTCAGTGCAGCGGGGTGCCGGCCCC
 ACCTACCACCCCTTGGAGACAGCCTACGTGCTGCCTGGCGCCCGGGGCCTCCCGGGACAACGGGTAC
 GCGCCGACCGCCTGCCCGGAGCAGCTACGATGCGGACCAGAAGCCGGGCCTGGAGCTGGCCCCGGCC
 GAGCCCGCTACCCGCGCGCGCGCGGGAGGAGTACAGCGACCCCGAAAGCCGCGAGTCGAGCCTGTCG
 GCGCGCTACTTCCGAGGGGAGGCGGCAGTGACCGACAGCTACTCCATGGACGCCTTCTTCATCTCGGAC
 GGGCGCTCGCGGCGCGCGGGCGGGGCGGGGCGGGGACGCGGGGGCTCGGGAGACGCGGGGGCGCC
 GGGGGGCGCGCGGGCGCGCGGGGGCGCAGGCGGGCGCGGGGACCGGCACGCGTGCGCCGAGTGCGGC
 AAGACCTACGCCACGTCTGTCGAACCTGAGCCGCCACAAGCAGACGACCCGAGCCTGGACAGCCAGCTG
 GCGCGCAAATGCCGACGTGCGGCAAGGCCTACGTGTCCATGCCGCGCTCGCCATGCACCTGCTCACG
 CACAACCTGCGCCACAAGTGCGGCGTCTGCGGCAAGGCCTTCTCGCGCCCTGGCTGCTGCAAGGTAC
 ATGCGCTCGCACACCGGCGAAAAGCCGTTTCGGCTGCGCGCACTGCGGCAAGGCCTTCGCCGACCGCTCC
 AACCTGCGCGCGCACATGCAGACGCACTCGGCCTTCAAGCACTACCGCTGCCGCCAGTGCACAAGAGC
 TTCGCGCTCAAGTCCTACCTCCACAAGCACTGCGAGGCGGCCTGCGCCAAGGCGGCCGAGCCACCCCG
 CCGACCCCGCGCGGCCCGGCCAGCTGA
 AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGAT
 ATCCTGGATTACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-RsrII
ACCN:	NM_033129
Insert Size:	924 bp


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_033129.3
RefSeq Size:	3594 bp
RefSeq ORF:	924 bp
Locus ID:	85508
UniProt ID:	Q9NQ03
Cytogenetics:	20p13
MW:	32.6 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]