

Product datasheet for **SC305595**

C20orf114 (BPIFB1) (NM_033197) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C20orf114 (BPIFB1) (NM_033197) Human Untagged Clone
Tag:	Tag Free
Symbol:	BPIFB1
Synonyms:	C20orf114; LPLUNC1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC305595 representing NM_033197.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCCGGCCCGTGGACCTTCACCCTTCTGTGGTTTGTCTGGCAGCCACCTTGATCCAAGCCACCCTC
AGTCCCAGTGCAGTTCTCATCCTCGGCCAAAAGTCATCAAAGAAAAGCTGACACAGGAGCTGAAGGAC
ACAACGCCACCGATCCTGCAGCAGCTGCCGCTGCTCAGTGCCATGCGGGAAAAGCCAGCCGGAGGC
ATCCCTGTGCTGGGCAGCCTGGTGAACACCGTCCTGAAGCACATCATCTGGCTGAAGGTCATCACAGCT
AACATCCTCCAGCTGCAGGTGAAGCCCTCGGCCAATGACCAGGAGCTGCTAGTCAAGATCCCCTGGAC
ATGGTGGCTGGATTCAACACGCCCTGGTCAAGACCATCGTGGAGTTCACATGACGACTGAGGCCCAA
GCCACCATCCGCATGGACACCGTGAAGTGGCCCCACCCGCCTGGTCTCAGTGACTGTGCCACCAGC
CATGGGAGCCTGCGCATCCAAGTCTGCATAAGCTCTCCTTCTGGTGAACGCCTTAGCTAAGCAGGTC
ATGAACCTCCTAGTGCCATCCCTGCCAATCTAGTGAAAAACAGCTGTGTCCCCTGATCGAGGCTTCC
TTCAATGGCATGTATGCAGACCTCTGCAGCTGGTGAAGGTGCCATTTCCCTCAGATTGACCGTCTG
GAGTTTGACCTTCTGTATCCTGCCATCAAGGGTGACACCATTAGCTCTACCTGGGGCCAAGTTGTTG
GACTCACAGGAAAGGTGACCAAGTGGTTCAATAACTCTGCAGCTTCCCTGACAATGCCACCCTGGAC
AACATCCCCTTACGCTCATCGTGAAGTCAAGGAGCTGGTGAAGGCTGAGTGGCTGCTGTCTCTCCA
GAAGAATTCATGGTCTGTTGGACTCTGTGCTTCTGAGAGTGCCATCGGCTGAAGTCAAGCATCGGG
CTGATCAATGAAAAGGCTGCAGATAAGCTGGGATCTACCCAGATCGTGAAGATCCTAAGTCAAGGACT
CCCGAGTTTTTATAGACCAAGGCCATGCCAAGGTGGCCAACTGATCGTGTGGAAGTGTTCCTCC
AGTGAAGCCCTCCGCCCTTTGTTCAACCTGGGCATCGAAGCCAGCTCGGAAGCTCAGTTTTACACCAA
GGTGACCAACTATACTCAACTGAATAACATCAGCTCTGATCGGATCCAGCTGATGAACCTGGGATT
GGCTGGTTCCAACCTGATGTTCTGAAAAACATCATCACTGAGATCATCCACTCCATCCTGCTGCCGAC
CAGAATGGCAAATTAAGATCTGGGGTCCAGTGTATTGGTGAAGGCCTTGGATTTCAGGCAGCTGAG
TCCTCACTGACCAAGGATGCCCTGTGCTTACTCCAGCCTCCTTGTGAAACCCAGCTCTCTGTCTCC
CAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCC
  
```

Restriction Sites:

SgfI-MluI

Plasmid Map:



ACCN:

NM_033197

Insert Size:	1455 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_033197.2
RefSeq Size:	1734 bp
RefSeq ORF:	1455 bp
Locus ID:	92747
UniProt ID:	Q8TDL5
Cytogenetics:	20q11.21
Protein Families:	Druggable Genome, Secreted Protein, Transmembrane
MW:	52.4 kDa
Gene Summary:	The protein encoded by this gene may be involved in the innate immune response to bacterial exposure in the mouth, nasal cavities, and lungs. The encoded protein is secreted and is a member of the BPI/LBP/PLUNC protein superfamily. This gene is found with other members of the superfamily in a cluster on chromosome 20. [provided by RefSeq, Jul 2008]