

Product datasheet for **SC305635**

MVB12B (NM_033446) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	MVB12B (NM_033446) Human Untagged Clone
Tag:	Tag Free
Symbol:	MVB12B
Synonyms:	C9orf28; FAM125B
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC305635 representing NM_033446. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAGAAGCTGCTTCTGCGTGAGACGGAGCCGGGACCCGCCGCCCGCAGCCACCGCCGCCGCCGCC
CAGCGGGGAACAGACCAGTCCACCATGCCTGAAGTCAAAGACCTCTCAGAAGCCTTGCCAGAAACGTCA
ATGGATCCCATCACGGGAGTCGGGGTGGTGGCTTCTCGGAACCGAGCCCCGACAGGCTATGACGTAGTT
GCACAGACAGCAGATGGTGTGGATGCTGACCTCTGGAAAGACGGCTTATTTAAATCCAAGGTTACCAGA
TACCTGTGTTTACAAGATCATTTTCAAAGAAAATAGTCATCTGGGAACGTGTTAGTAGATATGAAG
CTCATTGACATCAAGACACACTGCCTGTGGGCTTCATCCAATTCAGGAGACGGTGGACACACAGGAA
GTGGCTTTTAGGAAGAAGAGGCTGTGCATTAATTTATTCCACGGGATCAACGGAAGCTGCGATTTGT
GACATTCGGATCATGGGCCGGACCAAGCAGGCCCGCCTCAGTACACGTTTATTGGGAACTGAACAGC
ATGGGGATCTGGTATCGAATGGGCAGAGTACCAAGAAATCATGACTCATCTCAACCCACAACGCCTTCC
CAGTCATCAGCTGCCTCCACCCAGCCCCAACCTTCCAGGCACATCTCCCTAACACTTCTGCCACC
TTCGAGGCAGGAACAGCACCCGGAGCGACTACGAGTACCAGCACTCCAATTTGTATGCCATACAGCA
ATGGATGGTGTGCCTTTTATGATTTAGAGAAGTTTCTGTGTCCAGAAAGTATGCAGCCCTTTGAT
CTCCTGGGAATCACCATCAAATCTTAGCAGAAATCGAAAAGAGTACGAGTACAGCTTCCGCACAGAG
CAGAGCGCAGCCCGCAGGCTCCCGCCAGCCCCACCAGGTGTGAGCAGATCCCGCAGCTCTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_033446

Insert Size: 960 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_033446.2](#)

RefSeq Size: 4842 bp

RefSeq ORF: 960 bp

Locus ID: 89853

UniProt ID: [Q9H7P6](#)

Cytogenetics: 9q33.3

Protein Pathways: Endocytosis

MW: 35.6 kDa

Gene Summary: The protein encoded by this gene is a component of the ESCRT-I complex, a heterotetramer, which mediates the sorting of ubiquitinated cargo protein from the plasma membrane to the endosomal vesicle. ESCRT-I complex plays an essential role in HIV budding and endosomal protein sorting. Depletion and overexpression of this and related protein (MVB12A) inhibit HIV-1 infectivity and induce unusual viral assembly defects, indicating a role for MVB12 subunits in regulating ESCRT-mediated virus budding. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2011]
Transcript Variant: This variant (1) represents the predominant transcript, and encodes the longer isoform (1).