

Product datasheet for **SC305691**

LRRC62 (ELFN2) (NM_052906) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LRRC62 (ELFN2) (NM_052906) Human Untagged Clone
Tag:	Tag Free
Symbol:	ELFN2
Synonyms:	LRRC62; PPP1R29
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC305691 representing NM_052906.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCTGCGCCTGGGGCTGTGCGCGGGCGCTGCTGTGCGTGTGCCGGCCGGGTGCCGTGCGTGCCGAC
TGCTGGCTCATTGAGGGCGACAAGGGCTACGTGTGGCTGGCCATCTGCAGCCAGAACCAGCCGCCCTAC
GAGACCATCCCGCAGCACATCAATAGCACCGTGCACGACCTGCGGCTCAACGAGAACAAGCTCAAAGCC
GTGCTCTACTCCTCGCTCAACCGCTTTGGGAACCTCACCGACCTCAACCTCACCAAGAACGAGATCTCC
TACATCGAGGACGGTGCCTTCTGGGCCAGTCGAGCCTGCAGGTCTGCAGCTGGGTACAACAAGCTC
AGCAACCTGACGGAGGGCATGCTGCGAGGCATGAGCCGCCTGCAGTTCCTCTTTGTCCAGACAACCTC
ATCGAGGTGGTGACGCCACCGCCTTCTCCGAGTGCCCGAGCCTCATCAGCATCGACCTGCTCTCAAC
CGCCTCAGCCGCCTGGACGGTGCCACCTTTGCCAGCTCGCCAGCCTGATGGTGTGTGAGCTGGCCGGC
AACCCCTCAACTGTGAGTGCACCTCTCGGCTTCTGGCCTGGCTGGTGGTCTTCAACAACGTACC
AAGAACTACGACCGCCTGCAGTGTGAGTGCAGCGGGAGTTTGGCGGCTACCCGCTGCTGGTGGCCCGG
CCCTACCACAGCCTCAACGCCATCACCGTACTCCAGGCCAAGTGTGGAATGGCTCGCTGCCCGCCCGG
CCCGTGAGCCACCCACGCCCTACTCCACCGACGCCAGAGGGAGCCAGACGAGAAGCTCGGGCTTCAAC
CCCGACGAGATCCTTTCCGGTGGAGCCGCCGCCTCGTCCACCACGGATGCGTCCGGCAGGGCCAGCCATC
AAGCTGCACCACGTACGTTACCTCGGCCACCCCTGGTGGTCAATCCACACCCCTACAGCAAGATG
TACATCCTCGTGCAGTACAACAACAGCTACTTCTCCGACGTATGACCCTCAAGAACAAGAAGGAGATC
GTGACGCTGGACAACTGCGGGCGCACACTGAGTACACCTTCTGCGTGACCTCGTGCACAACAGCCGC
CGCTTCAACCACACCTGCCTGACCTTACCACGCGGGACCCCGTCCCGGAGACTTGGCGCCAGCACC
TCCACCACCCACTACATCATGACCATCTGGGCTGCCTCTTTGGCATGTTATCGTGTGGGAGCC
GTGTAATACTGCCTGCGCAAGCGGCGCATGCAGGAGGAGAAGCAGAAGTCTGTCAACGTCAAGAAGACC
ATCCTGGAGATGCGCTACGGGGTGTGATGTGGATGCCGGTCCATTGTGCACGCGCCAGAGCTGGGC
GAGCCTCCCGTGTGCCGTATCTCGCATGGCCTCCATCCCTCCATGATCGGGGAGAAGCTGCCACC
GCCAAGGGGTTGGAGGCCGGGTGGACACACCAAGGTAGCCACCAAGGCAACTATATCGAGGTGCGC
ACAGGCGCGGGGGACGGTCTGGCTCGGCCGAGGATGACCTCCCGGACCTCGAGAACGGCCAGGGC
TCGGCTGCAGAGATCTCCACATTGCCAAGGAGGTGGACAAGGTCAACCAGATCATTAACTGCATC
GATGCTCTCAAGCTGGACTCGGCCTCTTTCTGGGAGGCGGCAGCAGCAGTGGGACCCCGAGCTGGCC
TTCGAGTGCCAGTCCCTCCCTGCAGCTGCTGCCGCTCCTCAGCCACTGGCCCCGGGGCCCTGGAGCGG
CCAGCTTCTTTGCGCTCCCTACAAGGAGAGCTCCACCCACCCACTACAGCGCCAGCTGAGCGCCGAC
GCGGCCGTGACCCGAAGACCTGCAGCGTGTGCTCCAGTGGTTCATCAAGAGCGCCAAGGTCTTTAGC
CTGGACGTGCCCGACCATCCGGCCGCCACAGGGCTGGTAAGGGCGACTCCAAGTACATCGAGAAGGGC
AGCCCCCTCAACAGCCCGCTGGACCGGCTCCCGCTGGTCCGGCGGGCAGCGCGGGGGCAGCGCGGGG
GGCGGGGGCATCCACCACCTGGAGGTGAAGCCGGCCTACCACTGCAGCGAGCACCAGGACAGCTTTCC
GCCCTGTACTACGAGGAGGGTCCGACAGCCTGAGCCAGCGCGTGTCTTCTCAAGCCGCTGACCCGC
TCCAAGCGTGACTCCACCTACTCGAGCTCTCCCCAGACACTACTACTCAGGGTACTCTCCAGCCCC
GAGTACTCATCCGAGAGCACGACAAGATCTGGGAGCGCTTCCGGCCCTACAAGAAGCACCACCGGGAG
GAGGTGTACATGGCCCGGCTCACGCCCTGCGCAAGAAGTCCAGTTCGCCAAGGACGAGGATCTGCAT
GACATCCTTGATTACTGGAAGGGGTCTCCGCCAGCAGAAGCTGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_052906

Insert Size: 2463 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_052906.4](#)

RefSeq Size: 8379 bp

RefSeq ORF: 2463 bp

Locus ID: 114794

UniProt ID: [Q5R3F8](#)

Cytogenetics: 22q13.1

Protein Families: Transmembrane

MW: 89.7 kDa

Gene Summary: Inhibits phosphatase activity of protein phosphatase 1 (PP1) complexes.[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (1) represents the longer transcript and encodes the protein.
Sequence Note: This RefSeq record was created from transcript and genomic sequence data because no single transcript was available for the full length of the gene. The extent of this transcript is supported by transcript alignments and orthologous data.