

Product datasheet for **SC305789**

COLQ (NM_080538) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	COLQ (NM_080538) Human Untagged Clone
Tag:	Tag Free
Symbol:	COLQ
Synonyms:	CMS5; EAD
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC305789 representing NM_080538. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGACGGGCTCATCATTAGCCTCGCTCATCTGCTCATATTCAGGACTGCTCTGTTACTCGGCAGGC
TGCTTGGCCCTTCCAGCCTGGATCAGAAGAAGCGTGGTGGCCACAAAGCATGCTGCCTGCTGACGCCT
CCTCCACCACCCTGTTCCACCACCATTCTTCAGAGGTGGCCGAAGTCCGCTTCTCTCCCAGACATG
AAGAATCTCATGCTGGAAGTGGAGACCTCGCAGTCCCCGTCATGCAAGGCTCGTAGGCTCCCCTGGG
CCTCCCGGCCCCAGGGTCCACCGGGGCTTCTGGCAAGACAGGACCAAGGGAGAAAAGGGGGAGCTT
GGCCGACCAGGAAGGAAGGGTAGACCTGGCCCCCAGGTGTTCTGGCATGCCTGGGCCATCGGTTGG
CCAGGCCCTGAAGGACCCAGGGGTGAAAAGGTGACCTGGGTATGATGGGCTTGCCAGGGTCAAGAGGA
CCAATGGGCTCCAAGGGTACCCTGGATCCAGAGGGGAAAAGGGATCCAGAGGTGAAAAGGGTGACCTG
GGTCCCAAAGGAGAAAAGGGTTTCCAGGATTTCTGGAATGTTGGGGCAGAAAAGGTGAAATGGGTCCA
AAAGGTGAACCTGGGATAGCAGGACACCCAGGACCCACAGGAAGACCAGGAAAACGAGGCAAGCAGGGA
CAGAAAGGGGATAGTGGAGTTATGGGCCACCAGGCAAGCCTGGGCTTCTGGTCAACCTGGCCGTCG
GGGCCCCAGGCCCCACCTGCAGGACAATTATAATGGGACCCAAAGGGGAAAGAGGATTTCCCGGG
CCTCCAGGAAGATGCTTTTGGACCCACTATGAATGTGAATAACCCTTCTACGGGGAATCTGTGTAT
GGGCCCAGTTCCCGCGAGTTCCTGTGATTTTTGTGGTCAACAACAGGAGGAGCTTGAGAGGCTGAAC
ACCCAAAACGCCATTGCCTTCCGCAGAGACCAGAGATCTCTGTACTTCAAGGACAGCCTTGCTGGCTC
CCCATCCAGCTGACCCCTTTCTACCCTGTGGATTACACTGCAGACCAGCACGGCACCTGTGGGGATGGG
CTCCTGCAGCCTGGGAGGAGTGTGACGACGGTAACAGCGATGTGGGTGACGACTGCATCCGCTGTAC
CGTGCCTACTGTGGAGATGGTACCAGGATGAGGGTGTGGAGGACTGTGACGGCTCTGACTTTGGCTAC
CTGACATGCGAGACCTATCTCCCTGGGTGATATGGAGACCTGCAATGCACCCAGTACTGCTACATCGAC
TCCACGCCCTGCCGCTACTTCACCTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGCGC
```



View online »

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_080538

Insert Size: 1338 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_080538.2](#)

RefSeq Size: 3136 bp

RefSeq ORF: 1338 bp

Locus ID: 8292

UniProt ID: [Q9Y215](#)

Cytogenetics: 3p25.1

Protein Families: Druggable Genome, Transmembrane

MW: 46.5 kDa

Gene Summary: This gene encodes the subunit of a collagen-like molecule associated with acetylcholinesterase in skeletal muscle. Each molecule is composed of three identical subunits. Each subunit contains a proline-rich attachment domain (PRAD) that binds an acetylcholinesterase tetramer to anchor the catalytic subunit of the enzyme to the basal lamina. Mutations in this gene are associated with endplate acetylcholinesterase deficiency. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]
Transcript Variant: This variant (II) utilizes alternate exon 1A that is not present in variant I. The use of this exon results in a unique N-terminus compared to isoform I.