

Product datasheet for **SC305805**

ABHD16B (NM_080622) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ABHD16B (NM_080622) Human Untagged Clone
Tag:	Tag Free
Symbol:	ABHD16B
Synonyms:	C20orf135; dj591C20.1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC305805 representing NM_080622.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGTGCGTCATCTGTTCTGTGAAGGCGCTGGTGCAGCGTGTCAAGATCTACCTGACCGCCAGCTACACC
TACCCATTCCGCGGCTGGCCCGTGGCCTCCCGTGGGATGACGTGCGCGCCGTGGGCCGGAGCAGCAGC
CACCGGGCGCTGACCTGCGCGGACCGCGGGCGGCGTGTGGTTGCTGCGGGACGAGACGCTGGGCGGG
GATGCGCTGGGGCGCCCTCCACGTGGGGCGCGCAGCCAGGCGCAGTGCCTCTTGACGACGCTCCGCGAG
CTGCCCGCCAGCTCGTAGCTACGCGCTGGCCACTCGCTGGCCGCTGGCTCGTGTACCCCGGCTCC
GTGTCCCTGATGACGCGCGCTGCTGCCGCTGCTGCAGCAGGGCCAAGAGCGCCTCGTGGAGCGCTAC
CACGGCCGGCGCCAAAGCTGGTGGCTGTGACGGCAACGAGATCGACACTATGTTTCATGGACCGCCGC
CAGCACC CGGCGCAGCCACGTGCACGGCCGCGCCTCGTCATCTGCTGCGAAGGCAACGCGGGCTTCTAC
GAGATGGGCTGTCTGTGCACCGCTCGAGGCCGGCTACTCCGTGCTGGGCTGGAACCACCCGGCTTC
GGCAGCAGCACTGGCGTGCCTTCCCTCAGCACGACGCAACGCCATGGACGTGGTGGTTCGAGTACGCA
CTGCACCGCTGCACTTCCCGCCCGCGCACCTGGTGGTCTACGGCTGGTCTGTTGGCGGCTTCACGGCC
ACCTGGGCCACCATGACCTACCCGGAGCTGGGTGCACTGGTGGTGGACGCCACCTTCGACGACCTTGTG
CCGCTGGCGCTGAAGGTCATGCCCCACAGTTGGAAGGGGCTGGTGGTGGCGCACCGTGCAGCAGCACTTC
AACCTCAACGTGGCCGAGCAGCTGTGCTGCTACCCGGGGCCGGTGTGCTGCTCCGACGCACGCAGGAT
GACGTGGTCAGCACTTCCGGCCGCTGCCCCCTGTCACTGGTGGTGGAGGGCAACCGGGGCAAT
GAGCTGCTGTGCGCCTGCTGGAGCACCCTACCCGTCGTGATGGCGGAGAGGGCCGCGCCGTCGTC
ACCCGCTGGTGCAGCTGGCAGCTTGGCGCAGGAGGCCGCTTCTATGCACGCTACCGCTGGAGGAGG
GACTGGTGCCTGGCGCTGCTGCGCTCTACCGTGCACGCTGCGAAGAGGAGCTGGAGGGCGAGAGGCC
CTGGGGCCACACGACCCGCTTCCCATGGTGGTGGCCAGGGCCTGAGCTCGCGCGCGCCGGCGC
CTCGACTGTTCTGGCTCGGAAGCACCTCAAGAAGCTGGAGGGGACTCACTTCAGCCCTCTGGAGCCT
GAGGAGTTTCAGTTGCCCTGGCGCTGTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites:

SgfI-MluI

Plasmid Map:



ACCN: NM_080622

Insert Size: 1410 bp

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_080622.3</u>
RefSeq Size:	1776 bp
RefSeq ORF:	1410 bp
Locus ID:	140701
UniProt ID:	<u>Q9H3Z7</u>
Cytogenetics:	20q13.33
MW:	52.6 kDa