

Product datasheet for **SC305820**

Histone H2A Bbd (H2AFB3) (NM_080720) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Histone H2A Bbd (H2AFB3) (NM_080720) Human Untagged Clone
Tag:	Tag Free
Symbol:	Histone H2A Bbd
Synonyms:	H2AB2; H2ABB2; H2AFB; H2AFB3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC305820 representing NM_080720. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTGAACCGTCAGAATTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCCGAGGAGGAGAGACGCCGAGGGTCTCCGGTGTGGCGGCCGGGGCGGACCTGCTCTCGCACC
GTCCGAGCGGAGCTTTCTGTTTTTCAGTGAGCCAGGTGGAGCGCAGTCTACGGGAGGGCCACTACGCTCAG
CGCCTGAGTCGCACGGCGCCGGTCTACCTCGTGCAGTTATTGAGTACCTGACGGCCAAGGTCTGGAG
CTGGCGGGCAACGAGGCCAGAACAGCGGAGAGCGGAACATCACTCCCCTGCTGCTGGACATGGTGGTT
CACAACGACAGGCTACTGAGCACCCTTTTCAACACGACCACCATCTCAAGTGGCCCTGGCGAGGAC
TAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites:	SgfI-MluI
ACCN:	NM_080720
Insert Size:	348 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.



Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_080720.1</u>
RefSeq Size:	539 bp
RefSeq ORF:	348 bp
Locus ID:	83740
UniProt ID:	<u>P0C5Z0</u>
Cytogenetics:	Xq28
Protein Pathways:	Systemic lupus erythematosus
MW:	12.7 kDa
Gene Summary:	Histones are basic nuclear proteins that are responsible for the nucleosome structure of the chromosomal fiber in eukaryotes. Nucleosomes consist of approximately 146 bp of DNA wrapped around a histone octamer composed of pairs of each of the four core histones (H2A, H2B, H3, and H4). The chromatin fiber is further compacted through the interaction of a linker histone, H1, with the DNA between the nucleosomes to form higher order chromatin structures. This gene encodes a replication-independent histone that is a member of the histone H2A family. This gene is part of a region that is repeated three times on chromosome X, once in intron 22 of the F8 gene and twice closer to the Xq telomere. This record represents the most telomeric copy. [provided by RefSeq, Oct 2015]