

Product datasheet for **SC305827**

KRT72 (NM_080747) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	KRT72 (NM_080747) Human Untagged Clone
Tag:	Tag Free
Symbol:	KRT72
Synonyms:	CK-72; K6irs; K6IRS2; K72; KRT6; KRT6IRS2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC305827 representing NM_080747.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAGCCGCCAACTGACCCATTTCCCCCGGGGAGCGCCTGGGTTTCAGCGTTGCTCCGCGGTCTCT
TCTGGCGGGATCGGCAGCAGCTCCGCCTCATTCCGGGCCCGGGTCAAGGGCTCGGCCTCCTTTGGCAG
AAGAGCCTCTCCTGCCTTGGGGGCAGCCGAAGCCTGGCGCTCAGCGCTGCTGCACGGCGGGCGGGC
CGCCTGGGCGGCTTCGTGGGCACCGCCTTCGGCAGCGCCGGGCTGGGGCCCAAGTGTCCCTCCGTGTGC
CCACCCGGGGCATCCCTCAGGTCACCGTCAACAAGAGCCTCTGGCCCCGCTCAACGTGGAGATGGAC
CCCAGATCCAGAGGGTGCAGCGCCAGGAGCGGGAGCAGATCAAGGCGCTAAACAACAAGTTCGCCTCC
TTCATCGACAAGGTGCGTTTCTGGAGCAGCAGAATCAGGTGCTAGAGACCAAGTGAACCTCTACAG
CAGCTGGACTTGAACAACGCAGGAAGAACCTGGAGCCATTTATGAGGGCTACATCAGCAACCTGCAG
AAGCAGCTGGAGATGCTGTCTGGGACGGGGTGGAGCTGGATTCCGAGCTGAGGAACATGCAGGATTTG
GTGGAGGACTACAAGAAGAGGTATGAGGTGGAGATTAACAGACGCACAGCTGCTGAGAATGAGTTTGTG
GTGCTCAAGAAGGACGTGGATGCTGCTTACATGAATAAGTTGAGCTCCAGGCCAAGGTGGACTCCTTG
ACAGATGAGATTAATTTCTCAAGTGCCTTTATGAAGGGGAGATCACTCAGATCCAGTCCCACATCAGC
GACACGTCATCGTCTGTCAATGGACAACAACCGGGATCTGGACCTGGACAGCATCATTGCCGAGGTC
CGTGCCAGTACGAGGAGATTGCCCTAAAGAGCAAGGCCGAGGCTGAGACCCTGTACCAGACCAAGATC
CAGGAGCTGCAGGTCACAGCAGGCCAGCATGGGGATGACCTCAAGCTCACCAAGGCTGAAATCTCTGAG
CTCAACCCCTGATCCAGAGGATCCGCTCAGAGATAGGGAATGTGAAGAAGCAGTGTGCCGATCTGGAG
ACGGCCATCGCCGACGCTGAACAGCGGGGGGACTGCGCCCTGAAAGATGCCCGGGCCAAGCTGGATGAG
CTGGAGGGCGCCCTGCACCAGGCCAAGGAGGAGCTGGCAGGATGCTGCGTGAGTACCAGGAGCTCGTG
AGCCTGAAGCTGGCCCTGGATATGGAGATCGCCACCTACCGCAAGCTGCTGGAGAGCGAGGAGTGCAGG
ATGCTTGCGCAATATCCAATTCTGTGAGCATCTCCGTCATCAGCAGCACAATGCTGGGCAGGAGGG
GCTGGCTTCAGCATGGGCTTTGGCGCTCAAGCAGTTATAGCTACAAAACTGCAGCTGCAGACGTCAAG
ACCAAAGGCAGCTGTGGCAGTGAGCTCAAGGATCCCTTGCCAAAACTCGGGGAGCAGCTGTGCCACC
AAAAGGCCTCCAGATGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites:

Sgfl-MluI

Plasmid Map:



ACCN:

NM_080747

Insert Size:	1536 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_080747.2
RefSeq Size:	2030 bp
RefSeq ORF:	1536 bp
Locus ID:	140807
UniProt ID:	Q14CN4
Cytogenetics:	12q13.13
MW:	55.9 kDa
Gene Summary:	<p>Keratins are intermediate filament proteins responsible for the structural integrity of epithelial cells. The type II keratins consist of basic or neutral proteins which are arranged in pairs of heterotypic keratin chains coexpressed during differentiation of simple and stratified epithelial tissues. This gene encodes a type II keratin that is specifically expressed in the inner root sheath of hair follicles. The type II keratins are clustered in a region of chromosome 12q12-q13. Multiple transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Jun 2009]</p> <p>Transcript Variant: This variant (1) represents the predominant transcript and encodes the longer isoform (1). Variants 1 and 2 encode the same protein.</p>