

Product datasheet for **SC305831**

WFDC10A (NM_080753) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	WFDC10A (NM_080753) Human Untagged Clone
Tag:	Tag Free
Symbol:	WFDC10A
Synonyms:	C20orf146; dj688G8.3; WAP10
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC305831 representing NM_080753. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCACCCAGACTCTGCTGCCTGCCTGTTCTCTGTGTGCTGCTGCTGCAGGCCAGGGAGGATAC
CGTGACAAGAAGAGGATGCAGAAAACAGCTATCCCCAGAAATCAAAGTCTGCCAGCAGCAGCCTAAA
CTATATCTATGCAAACACTTATGTGAATCTCACCGAGATTGTCAAGCAAATAACATATGCTGTTCTACC
TACTGTGGGAATGTTTGCATGAGCATCCTGTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-Mlul



[View online »](#)

Plasmid Map:


ACCN: NM_080753

Insert Size: 240 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_080753.2](#)

RefSeq Size: 482 bp

RefSeq ORF: 240 bp

Locus ID: 140832

UniProt ID: [Q9H1F0](#)

Cytogenetics: 20q13.12

Protein Families: Secreted Protein, Transmembrane

MW: 8.9 kDa

Gene Summary: This gene encodes a member of the WAP-type four-disulfide core (WFDC) domain family. The WFDC domain, or WAP signature motif, contains eight cysteines forming four disulfide bonds at the core of the protein, and functions as a protease inhibitor. Most WFDC gene members are localized to chromosome 20q12-q13 in two clusters: centromeric and telomeric. This gene belongs to the telomeric cluster. [provided by RefSeq, Jul 2008]