

Product datasheet for **SC305909**

Menin (MEN1) (NM_130801) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Menin (MEN1) (NM_130801) Human Untagged Clone
Tag:	Tag Free
Symbol:	MEN1
Synonyms:	MEAI; SCG2
Vector:	<u>pCMV6 series</u>



[View online »](#)

Fully Sequenced ORF:	<p>>NCBI ORF sequence for NM_130801, the custom clone sequence may differ by one or more nucleotides</p> <pre> ATGGGGCTGAAGGCCGCCAGAACGCTGTTCCCGCTGCGCTCCATCGACGACGTGGTG CGCCTGTTTGCTGCCGAGCTGGGCCGAGAGCCGACCTGGTCTCCTTTCTTGGTG CTGGGCTTCGTGGAGCATTTTCTGGCTGTCAACCGCGTCATCCCTACCAACGTTCCCGAG CTCACCTTCCAGCCAGCCCCGCCCGCCGACCCGCTGGCGGCTCACCTACTTTCCCGTG GCCGACCTGTCTATCATCGCCGCCCTCTATGCCCGCTTACCCGCCAGATCCGAGGGCC GTCGACCTGTCCCTCTATCCTCGAGAAGGGGTGTCTCCAGCCGTGAGCTGGTGAAGAAG GTCTCCGATGTATATGGAACAGCCTCAGCCGCTCCTACTTCAAGGATCGGGCCACATC CAGTCCCTCTCAGCTTCATCACAGTTGGAGCCAGTAGGCACCAAATTGGACAGCTCC GGTGTGGCCTTTGCTGTGGTTGGGGCTGCCAGGCCCTGGGTCTCCGGGATGTCCACCTC GCCCTGTCTGAGGATCATGCCTGGGTAGTGTGGGCCCAATGGGGAGCAGACAGCTGAG GTCACCTGGCACGGCAAGGGCAACGAGGACCGCAGGGCCAGACAGTCAATGCCGTGTG GCTGAGCGGAGCTGGCTGTACCTGAAAGGATCATACATGCGCTGTGACCGCAAGATGGAG GTGGCGTTCATGGTGTGTGCCATCAACCCTTCCATTGACCTGCACACCGACTCGTGGAG CTTCTGCAGCTGCAGCAGAAGCTGCTCTGGCTGCTCTATGACCTGGGACATCTGGAAGG TACCCCATGGCCTTAGGGAACCTGGCAGATCTAGAGGAGCTGGAGCCACCCCTGGCCGG CCAGACCCACTCACCTCTACCAAGGGCATTGCCTCAGCCAAGACCTACTATCGGGAT GAACACATCTACCCTACATGTACCTGGCTGGCTACCACTGTGCGAACCGCAATGTGCGG GAAGCCCTGCAGGCTGGGCGGACACGGCCACTGTATCCAGGACTACAATACTGCGCG GAAGACGAGGAGATCTACAAGGAGTTCTTTGAAGTAGCCAATGATGTATCCCAACCTG CTGAAGGAGGCAGCCAGTTGCTGGAGGCGGGCAGGAGCGGCCGGGGAGCAAAGCCAG GGCACCCAGAGCCAAGGTTCCGCCCTCCAGGACCCTGAGTGCTTCGCCACCTGTGCGA TTCTACGACGGCATCTGCAAATGGGAGGAGGGCAGTCCCACGCCGTGTGCTGCACGTGGGC TGGGCCACCTTTCTGTGCACTCCTAGGCCGTTTTGAGGGACAGGTGCGGCAGAAGGTG CGCATAGTGAGCCGAGAGGCCGAGGCGCCGAGGCCGAGGAGCCGTGGGCGAGGAAGCC CGGAAGGCCGCGGGCGGGGCCACGGCGGGAGTCCAAGCCAGAGGAGCCCCCGCGCCC AAGAAGCCAGCACTGGACAAGGGCTGGGCACCGGCCAGGGTGCAGTGTGAGGACCCCCC CGGAAGCCTCCTGGGACTGTGCTGGCACAGCCGAGGCCCTGAAGGTGGCAGCAGCGCT CAGGTGCCAGCACCCGACGATCACCACCGCCGAGGGTCCAGTGCTCACTTTCCAGAGT GAGAAGATGAAGGCATGAAGGAGCTGCTGGTGGCCACCAAGATCAACTCGAGCGCCATC AAGCTGCAACTACGGCACAGTCGCAAGTGCAGATGAAGAAGCAGAAAGTGTCCACCCCT AGTGACTACACTCTGTCTTCTCAAGCGGCAGCGCAAAGGCCTCTGA </pre>
Restriction Sites:	Please inquire
ACCN:	NM_130801
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_130801.1 , NP_570713.1
RefSeq Size:	2738 bp
RefSeq ORF:	2736 bp
Locus ID:	4221
UniProt ID:	O00255
Cytogenetics:	11q13.1
Protein Families:	Druggable Genome, Transcription Factors
Gene Summary:	<p>This gene encodes menin, a tumor suppressor associated with a syndrome known as multiple endocrine neoplasia type 1. Menin is a scaffold protein that functions in histone modification and epigenetic gene regulation. It is thought to regulate several pathways and processes by altering chromatin structure through the modification of histones. [provided by RefSeq, May 2019]</p> <p>Transcript Variant: This variant (e1C) has a distinct 5' UTR from variant 1, but shares some 5' UTR sequence with variant e1B. Variants 1, e1B, e1C, e1D, e1E, and e1F1 all encode the same isoform (1).</p>