

Product datasheet for **SC305954**

ANKS3 (NM_133450) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ANKS3 (NM_133450) Human Untagged Clone
Tag:	Tag Free
Symbol:	ANKS3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC305954 representing NM_133450.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGTCCGAGCTCAGCGATGAAGCCAGCGAGCCGGAACCTCTGAACCCGAGCTTGCCATGTGGCACGGG
CTCGGGACACAGGTCAGCGGGGAGGAGCTGGATGTCCTCCGATCTTACACAGCTGCTCCATTGGC
CAGTATGAAGTGGTGAAGGAGTGTGTGCAGCGGAGAGATTAGATTTGAATAAGAAGAATGGTGGTGGC
TGGACCCCGCTGATGTATGCCTCTACATTGGCCACGACACAATCGTGCACCTGTGCTTGGAGCGGGG
GTGAGTGTGAATGTGCCGACCCAGAAGGGCAGACTCCACTGATGCTGGCCTCCAGCTGTGGCAACGAG
AGCATCGCTACTTTCTTCTCCAGCAAGGTGCAGAGCTAGAAATGAAAGACATCCAGGGCTGGACAGCC
CTTTCCACTGTACCAGCGCCGGGCACACAGCAGTGGTCAAGTTCCTCTTGGACAGTGGAGCCAATGCC
AACGTGAGGGAGCCGATATGTGGATTTACTCCCTTGATGGAAGCAGCTGTGCTGGCCATGAGATAATC
GTGCAGTATTTTCTGAATCACGGAGTCAAGGTGGACGCGAGAGACCACAGTGGAGCCACAGCCCGGATG
CTGGCCAAGCAGTACGGACACATGAAGATCGTGGCCTTGATGGACTTACTCGCCCTCTCTGCCAAG
AGCCTCTATCGGAGCCCAGAAAAGTACGAAGATCTGAGCTTCTGACGAGTCTGCCCTGCTCCTCAG
AGACAGAGGCCTTGCCGGAAGAAGGGTGTGAGCATCCACGAGGGACCGGAGCCCTGGCCAGGATCACA
GGCATTGGCCTGGGCGGAGAGCCACGGCCTCGCTATGAGCAGGCTCCTCCCGTGGCTATGTCACC
TTCAACAGCAGTGGCGAGAACCCTGGAAGAAGAGGGCTCTGCTGCCGGGATGTACCTCCCCATC
AATGAGCGGGATGTGGAGAGCAGCAGCAGCAGCAGTGGGAGGAACATGCTTTCTGTGCCAACCTG
GGGCCCCGTCAGAGCAGCAGCAGCAGCAGGAGGCTGGCCAGAGCCAGGGGCTCAGCAGCGAAGCTTCT
GTGGAGAGCAACGAGGACTCGGATCATGCCTGTAAAAGCTCAGCTCGAAACAAGCTAAAAGTTACATG
AAGACCAAGAATCCTGACAGCCAGTGGCCTCCCGCGCTGCAACTGACAGGGAAGGCTTTCTCGTGG
TCCAGCCCCCAGACTCAGAGGGCCCCCTACTCAGGACCCAGGACCTTGCCGCACTGCTGGAGCAGATC
GGGTGTCTGAAGTACCTGCAGGTGTTTGGAGAGCAGGACGTGGACCTCCGCATCTTTCTGACCTCACT
GAGAGCGACCTGAAGGAAATGGCATCACGCTGTTTGGGCCAAGAGGAAGATGACGTCCGCCATTGCC
CGCTGGCACAGCAGTGCCCGCCACCCGGGGATGCCCTGGAGTGGCCTACGCCGACCGGCTGGAGGCT
GAGATGCAGGAGCTCGCCATCCAGCTGCACAAGCGCTGCGAGGAGGTAGAGGCCACGCGGGGCCAGGTG
TGTGAGGAGCAGGAGTGCAGCGCTGGTGGAGAGTGCCTGCTGGAGCAGGACCGCCCGGAGGAC
CTCCAGGCCCGGCTGCGGAGAGCTGGGCCCTGGCCCGGATGCTGCCCTCGTCTGGACCAGCTGCGA
GCCTGTCAAGCTGAGCTGTCATCTCGAGTGAAGCAGGACCCCTGGTGCAGCCACTCTGGGCTA
GCCGTCAGGAGCTGACTCCAAGGGCTGGCAAGCGTCCCTGCAGGCCATGAGCCTCCCGAGCTCTCG
GGAGCCCTGGAGGACCGTGTCCGTGAGATGGGGCAAGCACTGTGCTTAGTGACCCAGAGCCTGGAGAAG
CTGCAGGTGCTGAACGGGAAGAAGTGGCGGGAGACCTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-Mlul

Plasmid Map:


ACCN: NM_133450

Insert Size: 1971 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_133450.3](#)

RefSeq Size: 2664 bp

RefSeq ORF: 1971 bp

Locus ID: 124401

UniProt ID: [Q6ZW76](#)

Cytogenetics: 16p13.3

MW: 72 kDa

Gene Summary: May be involved in vasopressin signaling in the kidney.[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (1) encodes the longer isoform (1).