

Product datasheet for **SC305959**

ZNF418 (NM_133460) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF418 (NM_133460) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF418
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC305959 representing NM_133460.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGAATACGACTCACTATAGGGCGGCCGGGAATTCGTGCACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGCAGGGCACTGTGGCATTGGAAGATGTGGCTGTGAACCTTTCCAGGAGGAGTGGAGTCTCCTTAGT
GAGGTTTCAGAGATGCCTTTACCATGACGTGATGCTGGAGAACTGGGTACTTATATCCTCCCTGGGTTGT
TGGTGTGGATCAGAAGATGAGGAGGCACCTTCTAAGAAGAGCATTCTATACAAAGAGTGTCTCAGGTC
AGCACTCTGGGGCAGGTGTGTCTCCCAAGAAGGCCCACTCTTGTGAAATGTGTGGCGCATCTTGGGA
GACATTTTGCACCTTGGCAGATCATCAGGGGACACATCACAAGCAGAACTGCACAGGTGTGAGGCATGG
GGGAATAAATTGTATGATAGTTCAAACCGTCCGACCAGAATCAGTACCTTGGAGAGAAACCTATAGA
AGCAGTGTGAGGAAGCATTGTTTGTGAAGAGGTGAAGTTCCATGTGTGAGGAGTGTCTATCTTC
ATTGAGTGTGAGGAAGCATTGTTTGTGAGGAGTGTGAGTGTGAGTGTGAGTGTGAGTGTGAGTGTGAGT
AAGTCAAACAGCAAACCTGAGTGTGAGTGTGAGTGTGAGTGTGAGTGTGAGTGTGAGTGTGAGTGTGAGT
TGCATGAAACATTCTAGCACCAAACACGTATTTGTTCAACAGCAGAGACTTCCCTCTAGAGAGGAATGT
TATTGCTGGGAATGTGGGAAATCCTTTAGCAAATATGATAGCGTCAGTAATCATCAGAGAGTTCACACT
GGGAAAAGACCTTATGAATGTGGAGAATGTGGGAAATCTTTAGTCATAAGGGCAGCCTTGTTCAGCAT
CAGCGAGTTCACACTGGGAAAAGACCTTATGAATGTGGAGAATGTGGGAAATCTTTAGTCATAAGGGC
AGCCTTGTTCAGCATCAGCGAGTTCATACTGGAGAAAGACCTTATGAGTGTGGAGAATGTGGGAAATCT
TTAGTCAAAATGGTACTCTCATTAAACATCAACGAGTTCACACTGGAGAAAGACCTTATGAGTGTGAA
GAATGTGGGAAATGTTTTACTCAGAAGGGCAATCTCATTCAACATCAACGAGTTCACACTAGTGAAGA
CCTTATGAGTGTGAAGAATGTGGGAAATGTTTTAGTCAAAAGGGCACCTAACTGAACATCATCGAGTT
CACACTAGAGAACGACCTTATGAGTGTGGAGAATGTGGGAAATCTTTAGTCGAAAGGGACACCTTAGG
AACCATCAGCGAGGTTCACACTGGAGAAAGACCTTACGAGTGTGGAGAATGTGGGAAATCTTTAGTCGA
AAGGGCAACCTCATTGAGCATCAGCGAAGCCACACTGGAGAAAGCCTTATGAGTGTAGAGAGTGTAGG
AAATTATTTAGGGCAAGTCCCACCTCATTGAACACCAGAGAGTTCACACTGGAGAAAGGCCATATGAA
TGTAATGAATGTGGGAAATCATTCAAGACAGCTCTGGGTTTCGTGTTTCATCAGAGAGTTCACACTGGA
GAAAAACCGTTTGTAGTGTAGTGAATGTGGGAAGTCATTTCTCAAAGCTGTTCCCTCCTTCGACATCGG
AGAGTTCATACTGGAGAAAGCCTTATGAATGTGGAGAATGTGGGAAGTCATTTTCATCAGAGCTCTTCC
CTCCTTCGACATCAGAAACTCACACTGCAGAAAGACCTTATGAGTGCAGAGAATGTGGGAAATCTTC
TCCAGTCTCCTTGAACACAGGAGAGTTCACACTGGAGAAAGCCTTATGAATGCAGGGAATGTGGAAAA
ACATTTACTCGAAGGTCTGCGCATTTTAAACATCAGAGACTTCATACTCGAGGAAAGCCTTACGAGTGC
AGCGAATGTGGGAAATCCTTTGCTGAAACCTTCAGTCTTACTGAACACAGGAGAGTACACACTGGAGAA
AGGCCTTATGAGTGCAGTGAATGTGGGAAATCATTTCATCGAAGCTCTTCTCCTTCGACATCAGAGA
GTTTCACACAGAAAGAAGTCCTTACAAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_133460

Insert Size: 2031 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_133460.2</u>
RefSeq Size:	3727 bp
RefSeq ORF:	2031 bp
Locus ID:	147686
UniProt ID:	<u>Q8TF45</u>
Cytogenetics:	19q13.43
MW:	77.9 kDa
Gene Summary:	Transcriptional repressor.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (3) lacks an alternate 5' exon and uses an alternate in-frame splice junction compared to variant 1. The resulting isoform (c) has a shorter and distinct N-terminus compared to isoform a. Variants 3 and 4 both encode the same isoform (c).