

## Product datasheet for **SC305989**

### CD99L2 (NM\_134445) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	CD99L2 (NM_134445) Human Untagged Clone
Tag:	Tag Free
Symbol:	CD99L2
Synonyms:	CD99B; MIC2L1
Vector:	<u>pCMV6 series</u>
Fully Sequenced ORF:	>NCBI ORF sequence for NM_134445, the custom clone sequence may differ by one or more nucleotides ATGGTGGCCTGGCGCTCGGCGTTCCTTGCTGCCTCGCTTCTCCTTGGCCACCCTGGTC CAGCGAGGATCTGGGACTTTGATGATTTTAACTGGAGGATGCAGTGAAAGAACTTCC TCAGTAAAGCGAAATGATTTTGACTTGGCTGATGCCCTGGATGATCGAAATGATCGAGAT GATGGCCGAGGAAACCAATTGCTGGAGGAGGAGGTTTTTCAGACAAGGATCTTGAAGAC ATAGTAGGGGTGGAGAATACAAACCTGACAAGGGTAAAGGTGATGGCCGGTACGGCAGC AATGACGACCCTGGATCTGGCATGGTGGCAGAGCCTGGCACCATTGCCGGGTGGCCAGC GCCCTGGCCATGGCCCTCATCGGTGCCGTCTCCAGCTACATCTCCTACCAGCAGAAGAAG TTCTGCTTCAGCATTAGCAGGGTCTCAACGCAGACTACGTGAAGGGAGAGAACCCTGGAA GCCGTGGTATGTGAGGAACCCCAAGTAAATACTCCACGTTGCACACGCAGTCTGCAGAG CCGCCGCCGCCCGCCGAACCAGCCCGGATCTGA
Restriction Sites:	Please inquire
ACCN:	NM_134445
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_134445.2](#), [NP\\_604394.1](#)

**RefSeq Size:** 3419 bp

**RefSeq ORF:** 573 bp

**Locus ID:** 83692

**UniProt ID:** [Q8TCZ2](#)

**Cytogenetics:** Xq28

**Protein Families:** Transmembrane

**Gene Summary:** This gene encodes a cell-surface protein that is similar to CD99. A similar protein in mouse functions as an adhesion molecule during leukocyte extravasation. Alternate splicing results in multiple transcript variants. [provided by RefSeq, May 2010]  
Transcript Variant: This variant (3) lacks three consecutive in-frame exons in one region and an in-frame exon in another region, compared to variant 5. The resulting isoform (3, also referred to as E4) lacks two internal segments, compared to isoform 5.