

## Product datasheet for **SC306044**

### CTIP2 (BCL11B) (NM\_138576) Human Untagged Clone

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | Expression Plasmids  |
| Product Name:             | CTIP2 (BCL11B) (NM_138576) Human Untagged Clone  |
| Tag:                      | Tag Free   |
| Symbol:                   | BCL11B   |
| Synonyms:                 | ATL1; ATL1-alpha; ATL1-beta; ATL1-delta; ATL1-gamma; CTIP-2; CTIP2; hRIT1-alpha; IDDFSTA; IMD49; RIT1; ZNF856B |
| Mammalian Cell Selection: | Neomycin   |
| Vector:                   | pCMV6-Entry (PS100001)   |
| E. coli Selection:        | Kanamycin (25 ug/mL)   |




[View online »](#)

**Fully Sequenced ORF:** >SC306044 representing NM\_138576.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGTCCCGCCGCAAAACAGGGCAACCCGACAGCATTGTCCCAGAGGGAGCTCATACCCCAGAGGCTGAC
CATGTGGAGGCCCATCTCGAAGAAGACGAGGGTCTGGAGATAGAGGAGCCAAGTGGCCTGGGGCTG
ATGGTGGGTGGCCCCGACCCTGACTGCTCACCTGTGGCCAGTGTCAAATGAACTTCCCCTGGGGGAC
ATCCTGGTTTTTATAGAGCACAAAAGGAAGCAGTGTGGCGGCAGCTTGGGTGCTGCTATGACAAGGCC
CTGGACAAGGACAGCCCGCCACCCTCTCACGCTCCGAGCTCAGGAAAGTGTCCGAGCCGGTGGAGATC
GGGATCCAAGTACCCCCGACGAAGATGACCACCTGCTCTACCCACGAAAGGCATCTGTCCAAGCAG
GAGAACATTGCAGGGCCGTGCAGGCCTGCCAGCTGCCAGCGGTGGCCCCATAGCTGCCTCTCCAC
CCTCACTCATCCGTGATCACTTCACTCTGCGTGCCTGGGCGCTCTCCCGCCCTGCCTCCCCTGCCG
TGCTGCAGCGCGCCCGGTCTCGGGTACGGGACTCAGGGTGGGGTGCAGACGGAGGCTCCCTTTGGA
TGCCAGTGTGAGTGTGTCAGGTAAGATGAGCCTTCCAGCTACATTTGCACAACATGCAAGCAGCCCTTC
AACAGCGCGTGGTTCTGCTGCAGCACGGCAGAACACGCACGGCTTCCGCATCTACCTGGAGCCCGGG
CCGGCCAGCAGCTCGCTACGCCCGGGCTCACCATCCCGCCGCCCTCGGGCCGGAGGCCGTGGCGCAG
TCCCCGCTCATGAATTTCTGGGCGACAGCAACCCCTTCAACCTGCTGCGCATGACGGGCCCATCTCTG
CGGGACCACCCGGGCTTCGGCGAGGGCCGCTGCCGGGCACGCCGCTCTTTCAGTCCCCCGCCGCGC
CACCACCTGGACCCGACCCGCTCAGTGCCGAGGAGATGGGGCTCGTCCGCCAGCACCCAGTGCCTTC
GACCGAGTCATGCGCCTGAACCCATGGCCATCGACTCGCCGCCATGGACTTCTCGGGCGGCTCCGC
GAGCTGGCGGGCAACAGCTCCACGCCCGCCCGTGTCCCCGGCCGGCAACCCATGACACCGCTC
CTGAACCCCTTCCAGCCAGCCCAAGTCCCGTTCTGAGCACGCCCGCTGCCGCCATGCCCTT
GGCGGCACGCGCCCGCCGACGCCAGCCAAGAGCAAGTGTGCGAGTTCTGCGGCAAGACCTTCAAG
TTCCAGAGCAATCTCATGTCACCGGCCAGTACACGGGCGAGAAGCCCTACAAGTGCAGCTGTGC
GACCACGCGTGTGCGCAGGCCAGCAAGCTCAAGCGCCACATGAAGACGCACATGCACAAGGCCGGCTCG
CTGGCCGGCCGCTCCGACGACGGGCTCTCGGCCGACGCTCCCCGAGCCGGCACGAGGAGTGGCG
GGCGAGGGCTCAAGGGCCGACGGTACTTCCGCCACCAGAGGCGACCCGTCGCTGGGCCACGAG
CCGGAGGAGGAGGACGAGGAGGAGGAGGAGGAGGAGGAGGAGTGTACTGGAGAACGAGAGCCGGCC
GAGTCGAGCTTACGATGACTCGGAGCTGAGCCGCAACCGCGAGAACGGCGGTGGTGGGTGCCGGG
GTCCCCGGCGCGGGGGCGCGGCCAAGGCGTGGCTGACGAGAAGGCGCTGGTGTGGCAAGGTC
ATGGAGAAGCTGGGCCTAGGCGCACTGCCGCACTACGGCGAGTCTGGCCGACAAGCAGAAGCGCGGC
GCCTTCTGAAGCGTGGCGGGCGGGGGACGCGGGGACGACGACGACGCGGGCGGCTGCGGGGAC
GCGGGCGCGGGCGCGGGTCAACGGGCGCGGGGGCGGCTTCGCGCCAGGCACCGAGCCCTTCCCAGG
CTTTCCCAGCAAGCCCGCGCCGCTGCCAGCCCCGGGCTCAACAGCGCCGCAAGCGCATCAAGGTG
GAGAAGGACCTGGAGTCCGCCCGCCGCGTCACTCCCGTCCGAGAACGTGTACTCGCAGTGGTGGTG
GGTACGCGGCGTCCGCGCACTTCAAGGACCCCTTCTGGGTTTACGGACGACGACAGTCCGCC
TTCGCCACGTGTCGAGCACTCGTCCGAGAACGGCAGCTGCGTTCTCACGCCCGCCGGGACCTG
CTGGACGGCGGCTCTCGGGCCGACGGCAGCCAGCGGAGGACACCCCGACCTGGGCGGCCCG
GGCCCCGGGCGGCCAGCTCCAAGGAGGGCCGCCGACGACAGTGTGAGTACTGCGGCAAGGTGTTT
AAGAACTGCAGCAACTTACGGTGCACCGCGGAGCCACACCGCGAGCGGCTTACAAGTGCAGCTG
TGCAACTACGCGTGCAGCAGAGCAGCAAGCTCACGCGCCACATGAAGACGCACGGGAGATCGGCAAG
GAGGTGTACCGTGCAGCATCTGCCAGATGCCCTTACGCGTCTACAGCACCTGGAGAAACATGAAA
AAGTGGCACGGGAGCACTTGTGACTAACGACGTCAAATCGAGCAGGCCGAGAGGAGCTAA
AGCGGACCGACGCTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGAT
ATCCTGGATTACAAGGATGACGACGATAAGGTTTAA
  
```

**Restriction Sites:** SgfI-RsrII

**Plasmid Map:**


**ACCN:** NM\_138576

**Insert Size:** 2685 bp

**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in *E. coli* are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_138576.3](#)

**RefSeq Size:** 8044 bp

**RefSeq ORF:** 2685 bp

**Locus ID:** 64919

**UniProt ID:** [Q9C0K0](#)

**Cytogenetics:** 14q32.2

**Protein Families:** Druggable Genome

**MW:** 95.5 kDa

**Gene Summary:** This gene encodes a C2H2-type zinc finger protein and is closely related to BCL11A, a gene whose translocation may be associated with B-cell malignancies. Although the specific function of this gene has not been determined, the encoded protein is known to be a transcriptional repressor, and is regulated by the NURD nucleosome remodeling and histone deacetylase complex. Four alternatively spliced transcript variants encoding distinct isoforms have been found for this gene. [provided by RefSeq, Aug 2013]  
Transcript Variant: This variant (1) encodes the longest isoform (1).