

Product datasheet for **SC306133**

TGIF2LY (NM_139214) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TGIF2LY (NM_139214) Human Untagged Clone
Tag:	Tag Free
Symbol:	TGIF2LY
Synonyms:	TGIFLY
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC306133 representing NM_139214. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCC**CGATCGCC**
 ATGGAGGCGCTGCAGACGGCCCGCTGAGACCCAAAGCCCGGTGGAAAAAGACGCCCGCGGAAGACC
 CAAAGCCAGCCCAAGACACCTCAATCATGTCGAGAAATAACGCAGATACAGGCAGAGTTCTTGCCCTTA
 CCAGAGCACAAGAAGAAGCGCAAGGAACTTGCCAGCCGAGTCCGTTAAGATCCTCCGCGACTGGATG
 TATAAGCATCGGTTTAAGGCCTACCCTTCAGAAGAAGAGAAGCAAATGCTGTGAGAGAAGACCAATTTG
 TCTTTGTTGCGGATTTCTAACTGGTTTATCAATGCTCGCAGACGCATTCTCCCGGATATGCTTCAACAG
 CGTAGAAACGACCCCATCATTGGCCACAAAACGGGCAAAGATGCCCATGCCACCCACCTGCAGAGCACC
 GAGGCGTCTGTGCCGGCCAAGTCAGGGCCAGTGGTCCAGACAATGTACAAAGCCTGCCCTGTGGCCCT
 TGCCAAAGGGCCAGATGTCAAGAGAGAAGCAACCAGATCCGGAGTCGGCCCTAGCCAGAAGCTCACCG
 GAA**TAG**
 A**CGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Restriction Sites:	Sgfl-MluI
ACCN:	NM_139214
Insert Size:	558 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).


[View online »](#)

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_139214.2</u>
RefSeq Size:	880 bp
RefSeq ORF:	558 bp
Locus ID:	90655
UniProt ID:	<u>Q8IUJ0</u>
Cytogenetics:	Yp11.2
Protein Families:	Transcription Factors
MW:	20.8 kDa
Gene Summary:	This gene encodes a member of the TALE/TGIF homeobox family of transcription factors. This gene lies within the male specific region of chromosome Y, in a block of sequence that is thought to be the result of a large X-to-Y transposition. The C-terminus of this protein is divergent from that of its chromosome X homolog (TGIF2LX), suggesting that this protein may act as a regulator of TGIF2LX. [provided by RefSeq, Jul 2008]