

Product datasheet for **SC306138**

LGI3 (NM_139278) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LGI3 (NM_139278) Human Untagged Clone
Tag:	Tag Free
Symbol:	LGI3
Synonyms:	LGIL4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC306138 representing NM_139278.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCGGGGCTGCGGGCCAGGGGGGCCGGGGCCGGGGCTGCTGGCGCTCTCCGCGCTCGGCTTCTGC
CTCATGCTGCAAGTCAGCGCTAAGAGGCCCCCAAGACGCCCCCTGCCGCCAGCTGCTCTTGCAACC
AGGGACACCGCCTTCTGCGTGGACTCAAAGGCGGTGCCAGGAACCTGCCCTCGGAGGTCACTCCCTG
ACCCTGGTGAATGCCGCCTTCTCAGAGATCCAGGATGGAGCGTTCTCCACCTCCCGCTGCTGCAGTTC
TTGTTACTCAACTCCAACAAGTTTACACTGATTGGAGACAACGCCTTCACAGGACTGTGCGACCTGCAG
TATCTCTTCATTGAGAACAATGACATCTGGGCACTATCCAAGTTCACCTTCGAGGACTCAAGTCTTG
ACTCACCTCTCGTGGCCAACAATAACCTGCAGACACTGCCAGAGACATCTCCGGCCCTGGACATC
CTGAATGACTTGGACCTGCGGGGCAACTCACTCAACTGTACTGCAAGGTGAAGTGGTGGTGGAGTGG
CTGGCACACACCAACACCACGGTGGCACCACATCTACTGCGCCAGCCGCCCGCTTCCAGGAGCACAAG
GTGCAGGACCTGCCGCTGCGGGAGTTGATTGCATCACCACAGATTTTGTGTTGTACCAGACCCTGGCC
TTCCCAGCAGTGTGCGGTGAGCCCTTCTCTACTCCAGTGACCTCTATTTGGCTCTGGCCAGCCAGGA
GTCAGTGCCTGCACCATCTGAAGTGGGACTATGTTGAGCGGCAGCTTCGAGACTATGATAGAATCCCA
GCCCCCTCTGCACTGCAAGCCGATGGTGGTGGACAGCCAGCTGTACGTGGTGGTGGCCAGCTG
TTTGGCGGCTCTTACATTTACCACTGGGATCCCAACACCACGCGCTTCACCAGGTGCAAGACATTGAC
CCGACGCGCTGCGCAAGCCTAACGACCTAGAAGCCTTCCGCATCGACGGTGACTGGTACTTTGCCGTG
GCTGACAGCTCCAAGGACGGCGCCACCAGCCTCTACCGCTGGCACCAGAATGGCTTCTACTCCCACCAG
GCACTGCACCCCTGGCACCCTGACACCGACCTGGAGTTTGTGGACGGCGAGGGCAAGCCACGGCTGATT
GTGTCCAGCAGCTCCCAGGCACCGCTCATCTATCAGTGGAGTCGACCCAGAAGCAGTTTGTGGCCAG
GGTGAGGTGACCCAGGTGCCTGATGCCCAAGCTGTGAAACACTTTCGTGCGCGGCCGACAGCTACCTG
TGCTCAGCCGCTACATTGGCGACTCCAAGATCTGCGCTGGGAGGGTACCCGCTTCTCGGAGGTGCAG
GCCCTGCCCTCCCGGGGCTCGTGGCCCTGCAGCCCTTCTTGTGGTGGCCGCCCTACCTGGCACTG
GGCAGTGATTTCTCTTACCCAGATCTACCAGTGGGATGAGGGACGACAGAAGTTGTACGGTTCCAG
GAGCTGGCTGTGAGGCTCCTCGGGCTTCTGCTACATGCCTGCTGGGACGCCAGCTACTCTGGCC
CCCAGCTTCAAGGGACAGACGCTGGTGTATAGACACATTGTGGTGGATCTCAGTGCCTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_139278

Insert Size: 1647 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_139278.2</u>
RefSeq Size:	3285 bp
RefSeq ORF:	1647 bp
Locus ID:	203190
UniProt ID:	<u>Q8N145</u>
Cytogenetics:	8p21.3
MW:	61.7 kDa
Gene Summary:	May participate in the regulation of neuronal exocytosis.[UniProtKB/Swiss-Prot Function]