

Product datasheet for **SC306183**

Calpain 12 (CAPN12) (NM_144691) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Calpain 12 (CAPN12) (NM_144691) Human Untagged Clone
Tag:	Tag Free
Symbol:	CAPN12
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC306183 representing NM_144691.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGATCGCC
ATGGCATCCAGCAGTGGGAGGGTACCATCCAGCTCGTGGATGAGGAGGCTGGGGTGGAGCCGGGGCC
CTGCAGCTTTTTCGGGGCCAGAGCTATGAGGCAATTCGGGCAGCCTGCCTGGATTCGGGGATCCTGTTC
CGCGACCTTACTTCCCTGCTGGCCCTGATGCCCTTGGCTATGACCAGCTGGGGCCGGACTCGGAGAAG
GCCAAAGGCGTGAAATGGATGAGGCCCATGAGTTCTGTGCTGAGCCGAAGTTCTCTGTGAAGACATG
AGCCGCACAGACGTGTGCAGGGGAGCCTGGGTAAGTGTGGTTCCTGCAGCTGCCGCCTCCCTTACT
CTGTATCCCCGGCTCCTGCGCCGGTGGTCCCTCCTGGACAGGATTTCCAGCATGGCTACGAGGCGTC
TTCCACTTCCAGCTCTGGCAGTTTGGCCGCTGGATGGACGTCGTGGTGGATGACAGGCTGCCCGTGCCT
GAGGGGAAGCTGATGTTCTGCTGCGCTCGGAACAGCGGAATGAGTTCTGGGCCCCACTCCTGGAGAAGGCC
TACGCCAAGCTCCACGGCTCCTATGAGGTGATGCGGGGCGCCACATGAATGAGGCTTTTGTGGATTC
ACAGGCGCGTGGGCGAGGTGCTCTATCTGAGACAAAACAGCATGGGGCTTTCTGCCCCTGCGCCAT
GCCCTGGCCAAGGAGTCCCTCGTGGGCGCCACTGCCCTGAGTGATCGGGGTGAGTACCGCACAGAAGAG
GGCCTGGTAAAGGGACACGCGTATTCCATCACGGGCACACACAAGGTGTTCTGGGCTTACCAAGGTG
CGGCTGCTGCGGCTGCGGAACCCATGGGGCTGCGTGGAGTGGACGGGGGCTGGAGCGACAGCTGCCCA
CGCTGGGACACACTCCCCACCGAGTGCCCGCATGCCCTGCTGGTAAAAAGGAGGATGGCGAGTTCTGG
ATGGAGCTGCGGGACTTCTCCTCCATTTGACACCGTGCAGATCTGCTCGTGAGCCCGGAGGTGCTG
GGCCCCAGCCCGAGGGGGGGCTGGCACGTCCACACCTTCAAGGCCGCTGGGTGCGTGGCTTCAAC
TCCGGCGGGAGCCAGCCTAATGCTGAAACCTTCTGGACCAATCCTCAGTCCGTTAACGCTGCTGGAG
CCTGATGAGGAGGATGACGAGGATGAGGAAGGGCCCTGGGGGGCTGGGGGGCTGCAGGGGCACGGGGC
CCAGCGCGGGGGGCGCACGCCCAAGTGACGGTCTTCTGTCCCTCATCCAGCGCAACCGCGCGGC
CTGAGAGCCAAGGGCCTCACTTACCTCACCGTTGGCTTCCACGTGTTCCAGATTCCAGAGGAGTGTG
GGCCTCTGGGATCCCCGCGCAGCCATGCGCTCCTGCCCGGCTGCTGCGCGCCGACCGCTCGCCCTC
AGCGCCCCCGCGACGTGACCCGCGCTGCTGCCTGCGTCCAGGCCACTACCTGGTGGTCCGAGCACC
GCCACGCGCGGACGAGGCTGACTTCACTCTGCGTGTCTTCCGAGCGCCGCCACAGGCCGTGGAG
ATCGACGACGTGATCAGCGCAGACCTGCAGTCTCTCCAGGGCCCCTACCTGCCCTGGAGCTGGGGTTG
GAGCAGCTGTTTCAGGAGCTGGCTGGAGAGGAGGAAGAACTCAATGCCTCTCAGTCCAGGCCTTACTA
AGCATTGCCCTGGAGCCTGCCAGGGCCATACCTCCACCCAGAGAGATCGGGCTCAGGACCTGTGAG
CAGCTGCTGCAGTGTTCGGGCATGGGCAAGCCCTGGCCTTACACCACTCCAGCAGCTCTGGGGCTAC
CTCCTGGAGTGGCAGGCCATATTAACAAGTTCGATGAGGACACCTCTGGAACATGAACTCTACGAG
CTGAGGCTGGCACTGAATGCAGCAGGCTTCCACCTGAACAACAGCTGACCCAGACCCCTACCCAGCCG
TACCGGGATAGCCGTCTGCGTGTGGACTTCGAGCGGTTCTGTGCTGTGGCCACCTCACCTGCATC
TTCTGCCACTGCAGCCAGCACCTGGATGGGGGTGAGGGGGTCACTGCTGACCCACAGACAGTGGATG
GAGGTGGCCACCTTCTCCTAG
AGCGGACCGACGCGTACGCGGCCCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGAT
ATCCTGGATTACAAGGATGACGACGATAAGGTTTAA
  
```

Restriction Sites: Sgfl-RsrII

Plasmid Map: □

ACCN: NM_144691

Insert Size: 2160 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_144691.4
RefSeq Size:	3129 bp
RefSeq ORF:	2160 bp
Locus ID:	147968
UniProt ID:	Q6ZSI9
Cytogenetics:	19q13.2
Protein Families:	Druggable Genome
MW:	81 kDa
Gene Summary:	The calpains, calcium-activated neutral proteases, are nonlysosomal, intracellular cysteine proteases. The mammalian calpains include ubiquitous, stomach-specific, and muscle-specific proteins. The ubiquitous enzymes consist of heterodimers with distinct large, catalytic subunits associated with a common small, regulatory subunit. This gene encodes a member of the calpain large subunit family. [provided by RefSeq, Jun 2012]